**AIA-100**

***KIP0091***

For Informational/Research Purposes Only

**Distributed By:**

Immuno-Biological Laboratories, Inc. (IBL-America)  
8201 Central Ave. NE, Suite P  
Minneapolis, Minnesota 55432, USA  
Phone: (888) 523-1246  
Fax.: (763) 780-2988  
Email: [info@ibl-america.com](mailto:info@ibl-america.com)  
Web: [www.ibl-america.com](http://www.ibl-america.com)

For Informational/Research Purposes Only

---

**LOT** : 110217/1

en

Read entire protocol before use.

## AIA-100

### I. INTENDED USE

Radioimmunoassay for the *in vitro* semi-quantitative measurement of human free anti-insulin antibodies (AIA) in serum and plasma.

### II. GENERAL INFORMATION

- A. Proprietary name : DIAsource AIA-100 Kit
- B. Catalog number : KIP0091 : 100 tests
- C. Manufactured by : DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :  
Tel : +32 (0)10 84.99.11      Fax : +32 (0)10 84.99.90

### III. CLINICAL BACKGROUND

#### A. Biological activity

The presence of circulating anti-insulin antibodies (AIA) in diabetics treated with insulin has been recognized as early as 1955. The highly purified insulin preparations, presently available, are less immunogenic than some of the previously used, less pure, preparations. Bovine insulin is more immunogenic than the porcine hormone. Also, it has recently been recognized that AIA may develop, in patients treated with human insulin.

The determination of circulating anti-insulin antibodies is of clinical importance for the following reasons :

- § The presence of free anti-insulin antibodies in plasma interferes with the determination of insulin by radioimmunoassay;
- § At very high titers, the anti-insulin antibodies may induce a state of insulin resistance;
- § Anti-insulin antibodies may influence the quality of the glycemic control, in diabetic patients, by prolonging the half life of insulin.

#### B. Clinical applications

- § Evaluation of the presence of free anti-insulin antibodies prior to the determination of insulin levels, by radioimmunoassay, in patients having received previous insulin therapy ;
- § Evaluation of states of insulin resistance;
- § Adjunct to the diagnosis of surreptitious insulin injections (factitious insulin-induced hypoglycaemia);
- § Monitoring of the evaluation of anti-insulin antibodies, in patients receiving the newly formulated human insulin preparations (new cases; patients in whom a previous treatment with porcine or bovine insulin is replaced by the administration of human insulin).

#### IV. PRINCIPLES OF THE METHOD

The presence of circulating anti-insulin antibodies in insulin treated diabetics is estimated on a semi-quantitative basis, by the determination of the binding of  $^{125}\text{I}$ -Tyr-A14-insulin to the serum fraction precipitated by the polyethylene glycol (PEG) (gamma globulins).

#### V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution
<b>Ag</b> $^{125}\text{I}$	1 vial 1 ml 74 kBq	red	<b>Add</b> 10 ml distilled water
TRACER: $^{125}\text{I}$ odine labelled INSULINE in phosphate buffer with bovine serum albumin , thymol and sodium azide (<0.1%)			
<b>CONTROL</b> <b>N</b> Controls - N = 1 to 3 (see exact values on vial labels) in human serum with thymol	3 vials lyophilised	yellow	<b>Add</b> 1ml distilled water
<b>PEG</b>	1 vial 105 ml	green	<b>Ready</b> for use
PEG: Polyethylene glycol (16%) in phosphate buffer with bovine serum albumin, Tween 20 and sodium azide (0.5%)			

**Negative Control:** The first control contains no anti-insulin antibodies. It allows the determination of the non-specific tracer binding.

**Positive Controls:** The two other controls respectively contain low and high levels of bovine free anti-insulin antibodies.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 100  $\mu\text{l}$  and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
3. Disposable polystyrene tubes (12 x 75 mm)
4. Plastic or aluminium foil
5. Incubator at 37°C
6. Vortex mixer
7. Magnetic stirrer
8. Centrifuge operating at 1500 g
9. Aspiration system (optional)
10. Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- A. **Tracers:** Reconstitute the tracer with 10 ml distilled water.
- B. **Controls:** Reconstitute the controls with 1 ml distilled water.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, tracer and controls are stable for 8 days at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid subsequent freeze-thaw cycles.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- After thawing, the samples should be mixed and centrifuged.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling.

Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure

1. Label polystyrene tubes in duplicate for each control, sample and total counts.
2. Briefly vortex controls and samples and dispense 100  $\mu\text{l}$  of each into the respective tubes.
3. Dispense 100  $\mu\text{l}$  of  $^{125}\text{I}$ odine labelled Insuline into each tube, including the tubes for total counts.
4. Shake the tube rack gently by hand to liberate any trapped air bubbles.
5. Cover the tubes (with plastic or aluminium foil) and incubate for 2 hours at 37°C.
6. Add 1 ml of the PEG solution (at room temperature) into each tube, except the total counts. Briefly vortex the tubes.
7. Incubate for 15 minutes at room temperature.
8. Centrifuge for 15 minutes at 1500 g. The use of a refrigerated centrifuge is not necessary, provided that the temperature does not rise up to 25°C.
9. Immediately aspirate (or decant) the supernatants carefully from each tube (except total counts). Be careful not to disturb the precipitate.
10. Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the total counts according to the following formula :

$$\text{B/T (\%)} = \frac{\text{Counts (Control or sample)}}{\text{Total Counts}} \times 100$$

#### XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time values.

AIA	cpm	B/T (%)
Total count	33780	
Controls:		
Negative Control	1265	3.7
Low Positive Control	6877	20.4
High Positive Control	19435	57.5
Sample 1	4054	12.0
Sample 2	7904	23.4
Sample 3	12363	36.3

#### XIII. PERFORMANCE AND LIMITATIONS

##### A. Precision

###### INTRA-ASSAY PRECISION

Serum	N	$\text{<}X\text{>} \pm \text{SD}$ (B/T x 100)	CV (%)	Serum	N	$\text{<}X\text{>} \pm \text{SD}$ (B/T x 100)	CV (%)
A	15	$3.2 \pm 0.4$	12.5	A	38	$3.5 \pm 0.7$	20
B	15	$20.4 \pm 0.6$	2.9	B	38	$20.8 \pm 1.0$	4.8
C	15	$59.6 \pm 1.0$	1.6		38	$60.0 \pm 2.0$	3.3

SD: Standard Deviation; CV: Coefficient of variation

#### XIV. INTERNAL QUALITY CONTROL

- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

#### XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

The assay was performed on 80 sera of patients who never received insulin therapy. The observed binding percentage of  $^{125}\text{I}$ -Tyr-A14-human-insulin was the following:  $5.5\% \pm 0.9\%$  (mean value  $\pm 1$  standard deviation). Consequently, one can consider that a binding percentage of  $^{125}\text{I}$ -Tyr-A14-human-insulin higher than the mean value  $+ 3$  standard deviations (8.2 %) corresponds to the presence of circulating anti-insulin antibodies.

#### XVI. PRECAUTIONS AND WARNINGS

##### Safety

For *in vitro* diagnostic use only.

This kit contains  $^{125}\text{I}$  (half-life: 60 days), emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations. This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals. All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

#### XVII. BIBLIOGRAPHY

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978)  
**The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.**  
Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976)  
**Equilibrium binding studies of insulin antibodies in diabetic subjects.**  
Clin. Chem. 22/17, 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970)  
**The antigenicity of pig insulin.**  
Diabetologia 6,15-20.
4. DESBUQUOIS B. and AURBACH G.D. (1971)  
**Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.**  
J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974)  
**Measurement of antibodies to insulin in serum.**  
Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975)  
**Insulin antibodies and the control of diabetes.**  
Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985)  
**Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.**  
Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973)  
**Radioimmuno-electrophoresis in the determination of insulin binding to IgG. Methodological studies.**  
Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983)  
**Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.**  
Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983)  
**Effects of species of origin, purification levels and formulation on insulin immunogenicity.**  
Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983)  
**Immunogenicity of recombinant DNA human insulin.**  
Diabetologica 25: 465-469.

#### XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS $\mu\text{l}$	SAMPLE(S) CONTROLS $\mu\text{l}$
Samples, controls	-	100
Tracer	100	100
Incubation		
PEG	-	1000
Incubation Centrifugation Separation		
Counting	Count tubes for 60 seconds	

DIAsource Catalogue Nr : KIP0091	P.I. Number : 1700439/en	Revision nr : 110217/1
-------------------------------------	-----------------------------	---------------------------

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#### XVII. BIBLIOGRAPHIE

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978)  
**The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.**  
Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976)  
**Equilibrium binding studies of insulin antibodies in diabetic subjects.**  
Clin. Chem. 2217, 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970)  
**The antigenicity of pig insulin.**  
Diabetologia 6,15-20.

4. DESBUQUOIS B. and AURBACH G.D. (1971)  
**Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.**  
J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974)  
**Measurement of antibodies to insulin in serum.**  
Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975)  
**Insulin antibodies and the control of diabetes.**  
Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985)  
**Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.**  
Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973)  
**Radioimmunolectrophoresis in the determination of insulin binding to IgG. Methodological studies.**  
Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983)  
**Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.**  
Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983)  
**Effects of species of origin, purification levels and formulation on insulin immunogenicity.**  
Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983)  
**Immunogenicity of recombinant DNA human insulin.**  
Diabetologica 25: 465-469.

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#### XVII. LITERATUR

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978) **The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.** Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976) **Equilibrium binding studies of insulin antibodies in diabetic subjects.** Clin. Chem. 22/17, 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970) **The antigenicity of pig insulin.** Diabetologia 6,15-20.

4. DESBUQUOIS B. and AURBACH G.D. (1971) **Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.** J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974) **Measurement of antibodies to insulin in serum.** Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975) **Insulin antibodies and the control of diabetes.** Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985) **Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.** Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973) **Radioimmunolectrophoresis in the determination of insulin binding to IgG. Methodological studies.** Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983) **Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.** Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983) **Effects of species of origin, purification levels and formulation on insulin immunogenicity.** Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFORD S. (1983) **Immunogenicity of recombinant DNA human insulin.** Diabetologia 25: 465-469.

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#### XVII. BIBLIOGRAFIE

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978)  
**The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.**  
Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976)  
**Equilibrium binding studies of insulin antibodies in diabetic subjects.**  
Clin. Chem. 22/1, 1089-1094.

3. DECKERT T. and GRUNDALL E. (1970)  
**The antigenicity of pig insulin.**  
Diabetologia 6, 15-20.
4. DESBUQUOIS B. and AURBACH G.D. (1971)  
**Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.**  
J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974)  
**Measurement of antibodies to insulin in serum.**  
Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975)  
**Insulin antibodies and the control of diabetes.**  
Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985)  
**Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.**  
Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973)  
**Radioimmunolectrophoresis in the determination of insulin binding to IgG. Methodological studies.**  
Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983)  
**Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.**  
Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983)  
**Effects of species of origin, purification levels and formulation on insulin immunogenicity.**  
Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983)  
**Immunogenicity of recombinant DNA human insulin.**  
Diabetologica 25: 465-469.

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#### XVII. BIBLIOGRAFIA

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978)  
**The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.**  
Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976)  
**Equilibrium binding studies of insulin antibodies in diabetic subjects.**  
Clin. Chem. 22(17), 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970)  
**The antigenicity of pig insulin.**  
Diabetologia 6, 15-20.

4. DESBUQUOIS B. and AURBACH G.D. (1971)  
**Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.**  
J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974)  
**Measurement of antibodies to insulin in serum.**  
Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975)  
**Insulin antibodies and the control of diabetes.**  
Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985)  
**Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.**  
Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973)  
**Radioimmuno-electrophoresis in the determination of insulin binding to IgG. Methodological studies.**  
Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983)  
**Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.**  
Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983)  
**Effects of species of origin, purification levels and formulation on insulin immunogenicity.**  
Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983)  
**Immunogenicity of recombinant DNA human insulin.**  
Diabetologica 25: 465-469.

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1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978)  
**The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.**  
Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976)  
**Equilibrium binding studies of insulin antibodies in diabetic subjects.**  
Clin. Chem. 22(7), 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970)  
**The antigenicity of pig insulin.**  
Diabetologia 6, 15-20.

4. DESBUQUOIS B. and AURBACH G.D. (1971)  
**Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.**  
J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974)  
**Measurement of antibodies to insulin in serum.**  
Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975)  
**Insulin antibodies and the control of diabetes.**  
Quarterly Journal of Medicine, New Series, XLIV, n° 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985)  
**Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.**  
Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973)  
**Radioimmuno-electrophoresis in the determination of insulin binding to IgG. Methodological studies.**  
Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983)  
**Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.**  
Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983)  
**Effects of species of origin, purification levels and formulation on insulin immunogenicity.**  
Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983)  
**Immunogenicity of recombinant DNA human insulin.**  
Diabetologica 25: 465-469.

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1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978) *The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.* Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976) *Equilibrium binding studies of insulin antibodies in diabetic subjects.* Clin. Chem. 22(17), 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970) *The antigenicity of pig insulin.* Diabetologia 6, 15-20.
4. DESBUQUOIS B. and AURBACH G.D. (1971) *Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.* J. Clin. Endocrinol. Metab. 33: 732-38.
5. DIXON K. (1974) *Measurement of antibodies to insulin in serum.* Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975) *Insulin antibodies and the control of diabetes.* Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-53.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985) *Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.* Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973) *Radioimmuno-electrophoresis in the determination of insulin binding to IgG. Methodological studies.* Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983) *Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.* Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983) *Effects of species of origin, purification levels and formulation on insulin immunogenicity.* Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983) *Immunogenicity of recombinant DNA human insulin.* Diabetologia 25: 465-469.

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### XVII. BIBLIOGRAFIA

1. ASLIN C.M., HARTOG M. and GOLDIE D.J. (1978) *The relationship between circulating free and bound insulin, insulin antibodies, insulin dosage and diabetic control in insulin treated diabetics.* Acta Endocrinologica 87, 330-338.
2. BAXTER R.C., YVE D.K. and TURTLE T.R. (1976) *Equilibrium binding studies of insulin antibodies in diabetic subjects.* Clin. Chem. 22(17), 1089-1094.
3. DECKERT T. and GRUNDALL E. (1970) *The antigenicity of pig insulin.* Diabetologia 6, 15-20.
4. DESBUQUOIS B. and AURBACH G.D. (1971) *Use of polyethylene glycol to separate free and antibody-bound peptides hormones in radioimmunoassays.* J. Clin. Endocrinol. Metab. 33: 732-738.
5. DIXON K. (1974) *Measurement of antibodies to insulin in serum.* Clin. Chem. 20/10, 1275-1281.
6. DIXON K., EXON P.D. and MALINS H. (1975) *Insulin antibodies and the control of diabetes.* Quarterly Journal of Medicine, New Series, XLIV, n<sup>o</sup> 176, 543-553.
7. WALDHAUSL, W.K., BRATUSCH-MARRAIN P., KRUSE V., JENSEN Y., NOWOTNY Y. and VIERHAPPER H. (1985) *Effect of insulin antibodies on insulin pharmacokinetics and glucose utilization in insulin-dependent diabetic patients.* Diabetes 34: 166-173.
8. CHRISTIANSEN A.H. (1973) *Radioimmuno-electrophoresis in the determination of insulin binding to IgG. Methodological studies.* Horm. Metab. Res. 5: 147-154.
9. FALHOLT K., HOSKAM J.A.M., KARAMANOS B.G., SUSSTRUNK H., VISWANATHAN M. and HEDING L.G. (1983) *Insulin-specific IgE in serum of 67 diabetic patients against human insulin (NOVO), porcine insulin, and bovine insulin. Four cases reports.* Diabetes care 6: 61-65.
10. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S. and GOLDMAN J. (1983) *Effects of species of origin, purification levels and formulation on insulin immunogenicity.* Diabetes 32:592-599.
11. FINEBERG S.E., GALLOWAY J.A., FINEBERG N.S., RATHBUN M.J. and HUFFERD S. (1983) *Immunogenicity of recombinant DNA human insulin.* Diabetologica 25: 465-469.

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	<u>Used symbols</u>
	Consult instructions for use
	Storage temperature
	Use by
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>CONTROL</b>	Control
<b>IVD</b>	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
	Wash solution concentrated
	Zero calibrator
	Calibrator #
	Control #
	Tracer
	Tracer
	Tracer concentrated
	Tracer concentrated
	Tubes
	Incubation buffer
	Acetonitrile
	Serum
	Specimen diluent
	Dilution buffer
	Antiserum
	Immunoadsorbent
	Calibrator diluent
	Reconstitution solution
	Polyethylene glycol
	Extraction solution
	Elution solution
	Bond Elut Silica cartridges
	Pre-treatment solution
	Neutralization solution
	Tracer buffer
	Microtiterplate
	HRP Conjugate
	HRP Conjugate
	HRP Conjugate concentrate
	HRP Conjugate concentrate
	Conjugate buffer
	Chromogenic TMB concentrate
	Chromogenic TMB solution
	Substrate buffer
	Stop solution
	Incubation serum
	Buffer
	AP Conjugate
	Substrate PNPP
	Biotin conjugate concentrate
	Avidine HRP concentrate
	Assay buffer
	Biotin conjugate
	Specific Antibody
	Streptavidin HRP concentrate
	Non-specific binding
	2nd Antibody
	Acidification Buffer
	Distributor

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