

Code No. 27766

## Bisphenol A (BPA) Assay Kit - IBL

### INTRODUCTION

Bisphenol A (BPA), principal raw material for epoxy resins and polycarbonate is suspected as one of the endocrine disrupting chemicals. BPA is a chemical building block for making polycarbonate plastic used for food containers and feeding bottles, which touches mouths directly, and is also used in production of epoxy resins for food and beverage can linings. In Japan, more than 0.48 million tons of BPA are produced each year and the water survey report from ministry of the environment says that there is 0.11 µg/L of BPA in rivers. So, we are anxious about the adverse impact for aquatic organisms and ecosystems. Vom Saal reported that in pregnant mice dosed at 2 µg per kg of body weight with a level of 1/25 of the threshold limit value per day, the prostate of the male offspring was enlarged.

This product is made for research use only by measuring BPA in serum of plasma; it is based on an ELISA using anti-rabbit IgG antibody coated solid-phase method.

### PRINCIPLE

This kit is based on a competitive ELISA using anti-rabbit IgG antibody solid-phase method. BPA standard or sample is taken on a anti-rabbit IgG antibody coated microplate and enzyme-labeled BPA and anti-BPA serum are added to cause a reaction. After the reaction is completed, the plate is washed and the enzyme reaction is allowed to proceed with addition of a substrate. Tetra Methyl Benzidine (TMB) is used as a coloring agent (Chromogen). The absorbance is then measured by using a plate reader. A standard curve is drawn, and the BPA concentration in the sample is obtained.

### MEASUREMENT RANGE

0.3 ~ 100 ng/mL

### INTENDED USE

This IBL's assay kit is capable for the quantitative determination BPA in human serum and plasma.

### KIT COMPONENT

|   |                                                                                         |             |
|---|-----------------------------------------------------------------------------------------|-------------|
| 1 | Precoated plate : Anti-Rabbit IgG Antibody                                              | 96 Well x 1 |
| 2 | Labeled BPA : HRP conjugated BPA                                                        | 6 mL x 1    |
| 3 | Standard : BPA<br>(0 ng/mL, 0.3 ng/mL, 1 ng/mL, 3 ng/mL, 10 ng/mL, 30 ng/mL, 100 ng/mL) | 0.3 mL x 7  |
| 4 | Anti-BPA Antiserum                                                                      | 6 mL x 1    |
| 5 | Unused number                                                                           |             |
| 6 | Chromogen : TMB solution                                                                | 15 mL x 1   |
| 7 | Stop solution : 1N H <sub>2</sub> SO <sub>4</sub>                                       | 12 mL x 1   |
| 8 | Wash buffer Conc. : (40X) 0.05% Tween20 in phosphate buffer                             | 50 mL x 1   |

### OPERATION MANUAL

#### 1. Materials needed but not supplied

- Plate reader (450nm)
- Graduated cylinder and beaker
- Plate mixer
- Paper towel
- Washing bottle for precoated plate
- Disposable test tube for "6, Chromogen"
- Micropipette and tip
- Deionized water
- Graph paper (semi-logarithmic)

#### 2. Preparation

Preparation of wash buffer  
 "8, Wash buffer Conc." is a concentrated (40X) buffer. Adjust the temperature of "8, Washing buffer Conc." to room temperature and then, mix it gently and completely before use. Dilute 50 mL of "8, Wash buffer Conc." with 1,950 mL of deionized water and mix it. This is the wash buffer for use. This prepared wash buffer shall be stored in refrigerator and used within 2 weeks after dilution.

The other reagents are ready for use.

#### 3. Measurement procedure

All reagents shall be brought to room temperature before use. Then mix it gently and completely before use. Make sure of no change in quality of the reagents.

Standard curve shall be prepared simultaneously with the measurement of test samples.

Usually, duplicate measurement shall be conducted.

Hereinafter, refer to Fig. 1 and Table 1.

- 1) Prepare necessary amounts of "1, Precoated plate" strips for measurement. 16 wells (2 strips) for preparation of the standard curve including blanks and 2 wells for measuring one sample are needed.
- 2) Add 20 µL of "3, Standard" (0 ng/mL) into well A-1 and A-2, and add "3, Standard" (0.3 ng/mL) into wells B-1 and B-2. In the same manner, add each concentration "3, Standard" into wells C-1, C-2, D-1, D-2, E-1, E-2, F-1, F-2, G-1 and G-2. (Leave wells H-1 and H-2 open for Blanks.) Add 20 µL of test sample into well of 3rd - 12th column.
- 3) Add 50 µL of "2, Labeled BPA" to each of the wells into which the standard or sample was put. Also, add 50 µL of "4, Anti-BPA Antiserum" to each of the wells into which the standard or sample was put.
- 4) Seal the precoated plate and stir it for 5 minutes with the plate mixer, and then incubate the precoated plate for 55 minutes at room temperature (20-30 °C).
- 5) Wash each well of the precoated plate vigorously with wash buffer using the washing bottle. Then, fill each well with wash buffer and leave the precoated plate laid for 15-30 seconds. Remove wash buffer completely from the precoated plate by snapping. This procedure must be repeated more than 7

times. Then, remove the remaining liquid from all wells completely by snapping the precoated plate onto paper towel.

In case of using a plate washer, after 4 times washing with plate washer, washing with above washing bottle must be repeated 3 times.

- 6) Take the required quantity of "6, Chromogen" into a disposable test tube. Then, pipette 100 µL from the test tube into the wells. Please do not return the rest of the test tube to "6, Chromogen" bottle to avoid contamination.
- 7) Incubate the precoated plate for 30 minutes at room temperature in the dark.
- 8) Add 100 µL of "7, Stop solution" to all wells. Mix the liquid by tapping the side of precoated plate. The liquid will turn yellow from blue by addition of "7, Stop solution".
- 9) Remove any dirt or drop of water on the bottom of the precoated plate and confirm there is no bubble on the surface of the liquid. Then, run the plate reader and conduct measurement at 450 nm against one of Blanks (H-1 or H-2) within 30 minutes after addition of "7, Stop solution".

Fig. 1




Table 1

| Sample                                                                                                                        | Standard | Test samples | Blank |
|-------------------------------------------------------------------------------------------------------------------------------|----------|--------------|-------|
|                                                                                                                               | Standard | Test sample  | —     |
| Labeled BPA                                                                                                                   | 20 µL    | 20 µL        | —     |
| Anti-BPA Antiserum                                                                                                            | 50 µL    | 50 µL        | —     |
| Mix for 5 minutes by plate mixer after sealing plate.                                                                         |          |              |       |
| Incubate for 55 minutes at room temperature (20-30°C).                                                                        |          |              |       |
| Washing 7 times                                                                                                               |          |              |       |
| Chromogen                                                                                                                     | 100 µL   |              |       |
| Incubate for 30 minutes at room temperature (shielded).                                                                       |          |              |       |
| Stop solution                                                                                                                 | 100 µL   |              |       |
| Tap the plate for mixing, and then read the plate at 450nm against a blank within 30 minutes after addition of Stop solution. |          |              |       |

### SPECIAL ATTENTION

- 1) Use serum or plasma as a sample. When samples are to be stored, keep them in a cool place or at -20°C.
- 2) Duplicate measurement of test samples and standard is recommended.
- 3) Use test samples in neutral pH range. The contaminations of organic solvent may affect the measurement.
- 4) Use only wash buffer contained in this kit for washing the precoated plate. Insufficient washing may lead to the failure in measurement.
- 5) Remove the wash buffer completely by tapping the precoated plate on paper towel. Do not wipe wells with paper towel.
- 6) "6, Chromogen" should be stored in the dark due to its sensitivity against light. "6, Chromogen" should be avoided contact with metals.
- 7) Measurement should be done within 30 minutes after addition of "7, Stop solution".

### CALCULATION OF TEST RESULT


- 1) Calculate the average of the absorbance values of "0 ng/mL BPA standard" wells, A-1 and A-2, (B<sub>0</sub>).
- 2) For the other wells, calculate the binding rates to B<sub>0</sub>, B / B<sub>0</sub> (%), as following formula.

$$B / B_0 (\%) = \frac{\text{absorbance of each well}}{B_0} \times 100$$

- 3) Plot the value (B/B<sub>0</sub>) of each standard on semi logarithmic graph paper, with the concentrations of standard on the logarithmic scale (abscissa) and the binding rates (B/B<sub>0</sub>) on the integer scale (ordinate), and draw a standard curve.
- 4) Obtain the binding rates (B / B<sub>0</sub>) of other wells which were used for the samples, in the same way. And read the BPA concentration of samples from the standard curve.

Example of measurement values and standard curve

| Conc. (ng/mL) | Absorbance (450nm)      | B / B <sub>0</sub> (%) |
|---------------|-------------------------|------------------------|
| 0             | 1.828 (B <sub>0</sub> ) | 100.0                  |
| 0.3           | 1.554                   | 85.0                   |
| 1             | 1.349                   | 73.8                   |
| 3             | 1.058                   | 57.9                   |
| 10            | 0.685                   | 37.5                   |
| 30            | 0.343                   | 18.8                   |
| 100           | 0.163                   | 8.9                    |



\* The typical standard curve is shown above. This curve can not be used to derive test results. Please run a standard curve for each assay.

### PERFORMANCE CHARACTERISTICS

#### 1. Sensitivity

- 1) The absorbance of BPA at 0 ng/mL ( $B_0$ ) is 1.5 or more.
- 2) The absorbance ratio of standard 100 ng/mL ( $B_{100}$ ) to standard 0 ng/mL ( $B_0$ ), ( $B_{100}/B_0 \times 100$ ) is less than 20 %.

#### 2. Intra Assay

| Measurement Value (ng/mL) | SD value | CV value (%) | n |
|---------------------------|----------|--------------|---|
| 1.32                      | 0.19     | 14           | 8 |
| 4.49                      | 0.38     | 8.5          | 8 |
| 14.73                     | 1.03     | 7            | 8 |
| 31.15                     | 1.72     | 5.5          | 8 |

#### 3. Specificity

| Compound                 | Cross-reactivity |
|--------------------------|------------------|
| BPA                      | 100 %            |
| BPA-Glucronide           | 85.0 %           |
| BPA-Na-Sulfate           | 68.0 %           |
| Bisphenol B              | 8.3 %            |
| Bisphenol F              | 0.2 %            |
| Diethylstilbestrol       | ≤ 0.02 %         |
| Hexesterol               |                  |
| 17β-Estradiol            |                  |
| 4-Heptylphenol           |                  |
| 4-n-Nonylphenol          |                  |
| 4-Propylphenol           |                  |
| 4-Hexyloxyphenol         |                  |
| 4-Pentylphenol           |                  |
| 4-Hexylphenol            |                  |
| 4-Butylphenol            |                  |
| 2-ter-Butylphenol        |                  |
| 4-Dodecylphenol          |                  |
| Di-n-Butyl-Phthalate     |                  |
| Benzyl-n-Butyl Phthalate |                  |
| Daidzein                 |                  |
| Genistein                |                  |
| Bis-GMA                  |                  |

### PRECAUTION FOR INTENDED USE AND/OR HANDLING

1. Bring each reagent and samples to room temperature and mix it to homogeneity without foaming.
2. "7, Stop solution" is a strong acid substance. Therefore, be careful not to have your skin and clothes contact "7, Stop solution" and pay attention to the disposal of "7, Stop solution".
3. Since the measurement results are affected by the time and temperature of incubation, the samples and standards must be incubated under the same conditions.
4. "1, Precoated plate" and "3, Standard" contain sodium azide. Therefore, dispose these materials after diluting them with large quantity of water to avoid production of explosive metallic azide.
5. After testing, disposal must be made in accordance with national and local regulation separating the infectious waste.
6. Wash hands after handling reagents.
7. Do not mix the reagents with the reagents from a different lot or kit.
8. All reagents should be stored at 2 - 8°C and do not use expired reagents. Put unused precoated plate strips into a sealed bag with a desiccant, and use by expiration date.
9. This kit is for research purpose only. Do not use for clinical diagnosis.
10. Handle samples as if there is capable of infection (HBvirus, HIV). The implements used for the assay should be treated in 0.1 % sodium hypochlorite solution or autoclave.

### STORAGE AND THE TERM OF VALIDITY

Storage Condition : 2 - 8°C  
 The term of validity : 12 months  
 (The expiry date is specified on outer box.)

### REFERENCE

1. Yan H, Takamoto M, Sugane K. Exposure to Bisphenol A Prenatally or in Adulthood Promotes T(H)2 Cytokine Production Associated with Reduction of CD4CD25 Regulatory T Cells. *Environ Health Perspect.* 2008 Apr;116(4):514-9.
2. Yamanaka H, Moriyoshi K, Ohmoto T, Ohe T, Sakai K. Efficient microbial degradation of bisphenol A in the presence of activated carbon. *J Biosci Bioeng.* 2008 Feb;105(2):157-60.
3. Mahalingaiah S, Meeker JD, Pearson KR, Calafat AM, Ye X, Petrozza J, Hauser R. Temporal variability and predictors of urinary bisphenol A concentrations in men and women. *Environ Health Perspect.* 2008 Feb;116(2):173-8.
4. Bonefeld-Jørgensen EC, Long M, Hofmeister MV, Vinggaard AM. Endocrine-disrupting potential of bisphenol A, bisphenol A dimethacrylate, 4-n-nonylphenol, and 4-n-octylphenol in vitro: new data and a brief review. *Environ Health Perspect.* 2007 Dec;115 Suppl 1:69-76.
5. Le HH, Carlson EM, Chua JP, Belcher SM. Bisphenol A is released from polycarbonate drinking bottles and mimics the neurotoxic actions of estrogen in developing cerebellar neurons. *Toxicol Lett.* 2008 Jan 30;176(2):149-56.
6. Vom Saal FS, Hughes C. An extensive new literature concerning low-dose effects of bisphenol A shows the need for a new risk assessment. *Environ Health Perspect.* 2005 Aug;113(8):926-33.
7. Inoue H, Tsuruta A, Kudo S, Ishii T, Fukushima Y, Iwano H, Yokota H, Kato S. Bisphenol a glucuronidation and excretion in liver of pregnant and nonpregnant female rats. *Drug Metab Dispos.* 2005 Jan;33(1):55-9.
8. Hiroi H, Tsutsumi O, Takeuchi T, Momoeda M, Ikezuki Y, Okamura A, Yokota H, Taketani Y. Differences in serum bisphenol a concentrations in premenopausal normal women and women with endometrial hyperplasia. *Endocr J.* 2004 Dec;51(6):595-600.
9. Inoue H, Yuki G, Yokota H, Kato S. Bisphenol A glucuronidation and absorption in rat intestine. *Drug Metab Dispos.* 2003 Jan;31(1):140-4.
10. Shibata N, Matsumoto J, Nakada K, Yuasa A, Yokota H. Male-specific suppression of hepatic microsomal UDP-glucuronosyl transferase activities toward sex hormones in the adult male rat administered bisphenol A. *Biochem J.* 2002 Dec 15;368(Pt 3):783-8.
11. Matsumoto J, Yokota H, Yuasa A. Developmental increases in rat hepatic microsomal UDP-glucuronosyltransferase activities toward xenoestrogens and decreases during pregnancy. *Environ Health Perspect.* 2002 Feb;110(2):193-6.
12. Inoue H, Yokota H, Makino T, Yuasa A, Kato S. Bisphenol a glucuronide, a major metabolite in rat bile after liver perfusion. *Drug Metab Dispos.* 2001 Aug;29(8):1084-7.
13. Yokota H, Iwano H, Endo M, Kobayashi T, Inoue H, Ikushiro S, Yuasa A. Glucuronidation of the environmental oestrogen bisphenol A by an isoform of UDP-glucuronosyltransferase, UGT2B1, in the rat liver. *Biochem J.* 1999 Jun 1;340 ( Pt 2):405-9.
14. Kamiura T, *et al.* Determination of Bisphenol A in Air. *J. Environ. Chem.* 1997 7:275-9.
15. Howdeshell KL, Hotchkiss AK, Thayer KA, Vandenberg JG, vom Saal FS. Exposure to bisphenol A advances puberty. *Nature.* 1999 Oct 21;401(6755):763-4.
16. Krishnan AV, Stathis P, Permuth SF, Tokes L, Feldman D. Bisphenol-A: an estrogenic substance is released from polycarbonate flasks during autoclaving. *Endocrinology.* 1993 Jun;132(6):2279-86.
17. Vom Saal FS. Sexual differentiation in litter-bearing mammals: influence of sex of adjacent fetuses in utero. *J Anim Sci.* 1989 Jul;67(7):1824-40.
18. Ohkuma H, Abe K, Ito M, Kokado A, Kambegawa A, Maeda M. Development of a highly sensitive enzyme-linked immunosorbent assay for bisphenol A in serum. *Analyst.* 2002 Jan;127(1):93-7.
19. Y. Usuki: *Bio. Clinica.* 2000 15: 55-8. (in Japanese)
20. Y. Takai: Proceedings of the 51st Annual meeting of the Japan Society for Obstetrics and Gynecology, Tokyo, April, p.S300, 1999 (in Japanese)
21. A. Kokado: Proceedings of the 120th Annual meeting of the Japan Society for Pharmaceutical, Gifu, March, p.179 (4), 2000 (in Japanese)
22. H. Ohkuma: Book of abstract of the 5th Annual meeting of the Immunochemical Society of Japan at Kobe, June, p. 21, 2000 (in Japanese)
23. H. Kakishima: Book of abstract of the 4th Annual meeting of the Japan Society of Endocrine Disrupter Research at the Japanese Society City of Tsukuba, December, p.421, 2001 (in Japanese)
24. N. Kawagoe: Book of abstract of the 4th Annual meeting of the Japan Society of Endocrine Disrupter Research at the Japanese Society City of Tsukuba, December, p.147, 2001 (in Japanese)

Version 1.