



User's Manual

Salivary Progesterone ELISA



IB79304



96 Wells

02/08

Contents / Inhaltsverzeichnis

1	INTRODUCTION	3
2	PRINCIPLE.....	3
3	WARNINGS AND PRECAUTIONS.....	3
4	REAGENTS.....	4
5	SPECIMEN COLLECTION AND PREPARATION	5
6	ASSAY PROCEDURE	6
7	EXPECTED NORMAL VALUES.....	7
8	QUALITY CONTROL.....	7
9	PERFORMANCE CHARACTERISTICS	8
10	LIMITATIONS OF PROCEDURE	10
11	REFERENCES	10
12	SYMBOLS USED WITH IBL-AMERICA ELISA 'S.....	12

1 INTRODUCTION

1.1 Intended Use

Enzyme immunoassay for the *in vitro diagnostic* quantitative measurement of active free progesterone (a female hormone) in saliva.

Measurements obtained by this device may be used in the diagnosis and treatment of disorders of the ovaries or placenta and can be used as an aid for prediction of ovulation.

1.2 Summary and Explanation

Progesterone (4-pregnene-3, 20-dione) is a C21 steroid hormone containing a keto-group (at C-3) and a double bond between C-4 and C-5. Like other steroids, it is synthesized from cholesterol via a series of enzyme-mediated steps (1)

The steroid hormone Progesterone is a female sex hormone which, in conjunction with estrogens, regulates the accessory organs during the menstrual cycle and it is particularly important in preparing the endometrium for the implantation of the blastocyte and in maintaining pregnancy (2)

In non-pregnant women progesterone is mainly secreted by the corpus luteum whereas in pregnancy the placenta becomes the major source (3,4). Minor sources for progesterone are the adrenal cortex for both sexes and the testes for males.

The Progesterone level in saliva represents the concentration of the active free Progesterone.

2 PRINCIPLE

The IBL-AMERICA Salivary Progesterone ELISA kit is based on the competition principle and the microplate separation.

An unknown amount of Progesterone present in the sample and a fixed amount of Progesterone conjugated with horseradish peroxidase compete for the binding sites of a rabbit polyclonal Progesterone -antiserum coated onto the wells. After one hour incubation the microplate is washed to stop the competition reaction. After adding the substrate solution, the concentration of Progesterone is inversely proportional to the optical density measured.

3 WARNINGS AND PRECAUTIONS

1. For in-vitro diagnostic use only. For professional use only.
2. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
3. Do not mix reagents of different lots. Do not use expired reagents.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 – 8°C in the sealed foil pouch and used in the frame provided.
5. Avoid contact with Stop Solution, 0.5M H₂SO₄. It may cause skin irritation and burns.
6. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
7. Use separate pipette tips for each sample, control and reagent to avoid cross contamination.
8. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
9. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
10. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
11. Kit calibrators have been checked for HIV-1/2 and HCV antibodies and HBsAg and found to be negative, but the calibrators and patient samples should be handled as potentially infectious.
12. Some reagents contain Proclin, BND and MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
13. Allow the reagents to reach room temperature (21-26°C) before starting the test. Temperature will affect the absorbance readings of the assay. However, Values for the patient samples will not be affected.
14. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12x8 (break apart) strips, 96 wells;
Wells coated with a anti-Progesterone antibody (polyclonal).
2. **Standard (Standard 0-6)**, 7 vials, 1 mL each, ready to use;
Concentrations: 0; 10; 50; 100; 500; 1000; 5000 pg/mL
Conversion: 1000 pg/mL x 3.18 = nmol/L.
* contain 0.03% Proclin 300, 0.015% BND and 0.010% MIT as preservatives.
3. **Control**, 2 vials, 1.0 mL each, ready to use;
Control values and ranges please refer to vial label or QC-Datasheet.
* contain 0.03% Proclin 300, 0.015% BND and 0.010% MIT as preservatives.
4. **Enzyme Conjugate**, 1 vial, 26 mL, ready to use;
Progesterone conjugated to horseradish peroxidase;
* contain 0.03% Proclin 300, 0.015% BND and 0.010% MIT as preservatives.
5. **Substrate Solution**, 1 vial, 25 mL, ready to use;
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use;
contains 1 N acidic solution.
Avoid contact with the stop solution. It may cause skin irritations and burns.
7. **Wash Solution**, 1 vial, 30 mL (40X concentrated);
see „Preparation of Reagents“.

- * BND = 5-bromo-5-nitro-1,3-dioxane
MIT = 2-methyl-2H-isothiazol-3-one

Note: Additional *Standard 0* for sample dilution is available upon request.

4.2 Materials required but not provided

1. Calibrated EIA reader adjusted to read at 450 nm
2. Precision pipettes (100 and 200 µl)
3. Distilled or Deionized water
4. Timer (60 min. range)
5. Reservoirs (disposable)
6. Test tube or microtube rack in a microplate configuration
7. semi logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for two months if stored as described above.

4.4 Reagent Preparation

Bring all reagents to room temperature before use.

Wash Solution

Add deionized water to the 40X concentrated *Wash Solution*.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL deionized water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 2 weeks at room temperature.

5 SPECIMEN COLLECTION AND PREPARATION

Samples containing sodium azide should not be used in the assay. The saliva samples should be completely colorless. Even the slightest red color shows blood contamination. Such blood contamination will give falsely elevated concentration values. In case of visible blood contamination the patient should discard the sample, rinse the sampling device with tap water, also rinse the mouth with (preferably) cold water, wait for 10 minutes and take a new sample. Do not chew anything during the sampling period. Any pressure on the teeth may result in falsely elevated measurements due to an elevated content of gingival liquid in the saliva sample.

5.1 Specimen Collection

It is recommended to collect saliva samples with commercially available equipment (e.g. SALIVA-SET 5, cat.-no. DE7269-5). Do not use any PE devices or Salivettes for sampling; this in most cases will result in significant interferences. Glass tubes can be used as well, but in this case special attention is necessary for excluding any interference caused by the stopper. As food might contain significant amounts of steroid hormones samples preferably should be taken while fasting. If fasting should be a problem at least any food of animal origin (meat or dairy products) should be avoided prior to finalizing the collection. In the morning breakfast should be done only after finalizing the collection procedure. During the day the collection period should be timed just before an anticipated meal. As the steroid hormone secretion in saliva as well in serum shows an obvious dynamic secretion pattern throughout the day it is important to always collect 5 samples during a 2 hour period; this means every 30 minutes one sample. If possible the volume of each single sample should be a minimum of 0.5 ml (better 1 ml). Saliva flow may be stimulated by drinking water. This is allowed and even recommended before and during the collection period. Drinking of water is not allowed during the last 5 minutes before taking the single samples. The typical timing for a morning collection period would be as follows. Wake-up at 6:00 AM, drinking water and brushing teeth, 1st sample at 6:15 AM, followed by samples at 6:45 AM, 7:15 AM, 7:45 AM, and 8:15 AM, followed by breakfast at 8:25 AM. The typical timing for an afternoon collection period would be like: 1st sample at 5:00 PM, followed by samples at 5:30 PM, 6:00 PM, 6:30 PM, 7:00 PM, followed by dinner at 7:10 PM. Modest variation in the collection timing will not be critical, and the collection time-frame can be extended up to 3 hours.

5.2 Specimen Storage and Preparation

Saliva samples in general are stable at ambient temperature for several days. Therefore mailing of such samples by ordinary mail without cooling will not create a problem. Storage at 4°C can be done for a period of up to one week. Whenever possible samples preferable should be kept at a temperature of -20°C. Even repeated thawing and freezing is no problem. Each sample has to be frozen, thawed, and centrifuged at least once anyhow in order to separate the mucins by centrifugation. Upon arrival of the samples in the lab the samples have to stay in the deep freeze at least overnight. Next morning the frozen samples are warmed up to room temperature and mixed carefully. Then the samples have to be centrifuged for 5 to 10 minutes. Now the clear colorless supernatant is easy to pipette. If the sample should show even a slight reddish tinge it should be discarded. Otherwise the concentration value most probably will be falsely elevated. Due to the episodic variations of the steroid secretion we highly recommend the strategy of multiple sampling. If such a set of multiple samples has to be tested the lab (after at least one freezing, thawing, and centrifugation cycle) has to mix the aliquots of the 5 single samples in a separate sampling device and perform the testing from this mixture.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Standard 0* solution and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) Dilution 1:10: 10 µl saliva + 90 µl *Standard 0* (mix thoroughly)
- b) Dilution 1:100: 10 µl of dilution a + 90 µl *Standard 0* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Assay Procedure

Each run must include a standard curve.

1. Secure the desired number of coated strips in the frame holder.
2. Dispense **100 µl** of each Progesterone Standards and Controls into appropriate wells.
3. Dispense **100 µl** of each sample into selected wells.
4. Dispense **200 µl** of Enzyme Conjugate into each sample and standard well and mix the plate for thoroughly for 10 seconds.
5. Incubate for **60 minutes** at room temperature.
6. Briskly shake out the contents of the wells and rinse the wells 3 times with diluted Wash Solution (400 µl per well). Strike the inverted wells sharply on absorbent paper towel to remove residual droplets.

Important note:

The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!

7. Add **200 µl** of Substrate Solution to each well.
8. Incubate for **15 minutes** at room temperature.
9. Stop the reaction by adding **100 µl** of Stop Solution to each well.
10. Determine the absorbance of each well at 450 ± 10 nm.
It is recommended to read the wells within 10 minutes.

6.3 Calculation of Results

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical(Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred calculation method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted. For the calculation of the concentrations this dilution factor has to be taken into account.

Conversion:

Progesterone Conversion: $1000 \text{ pg/mL} \times 3.18 = \text{nmol/L}$.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

Standard	Absorbance Units (450 nm)
Standard 0 (0 pg/mL)	1.89
Standard 1 (10 pg/mL)	1.71
Standard 2 (50 pg/mL)	1.54
Standard 3 (100 pg/mL)	1.39
Standard 4 (500 pg/mL)	0.95
Standard 5 (1000 pg/mL)	0.71
Standard 6 (5000 pg/mL)	0.44

7 EXPECTED NORMAL VALUES

In order to determine the normal range of SLV Progesterone, saliva samples from 80 adult male and 120 female apparently healthy subjects, ages 21 to 75 years, were collected in the morning and analyzed using the IBL-AMERICA SLV Progesterone ELISA kit.

The following ranges were calculated from this study.

	Age group	Salivary progesterone pg/mL
Women	21 - 50 yrs. Follicular phase n = 40	19.6 – 86.5 pg/mL
	21 - 50 yrs. Luteal phase n = 40	99.1 – 332.6 pg/mL
	51 - 75 yrs. Postmenopausal n = 40	6.0 – 56.4 pg/mL
Men	21 - 50 yrs. n = 40	12.7 – 57.4 pg/mL
	51 - 75 yrs. n = 40	15.2 – 65.1 pg/mL

Therapy should not be decided based on results alone. The results should be correlated to other clinical observations and diagnostic tests.

Furthermore, we recommend that each laboratory establish its own range for the population tested, because the values differ between age, new born, children, adolescents and adults.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or IBL-AMERICA directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Sensitivity

The lowest detectable level of progesterone that can be distinguished from the Zero Standard is 3.8 pg/mL at the 95 % confidence limit.

9.2 Specificity

The following materials have been evaluated for cross reactivity. The percentage indicates cross reactivity at 50% displacement compared to Progesterone.

Steroid	% cross reactivity
Progesterone	100.0
Desoxycorticosterone	1.1
Pregnenolone	0.35
17 α -Hydroxyprogesterone	0.3
Corticosterone	0.2
11-Desoxycortisol	0.1
Estriol	<0.1
Estradiol 17 β	<0.1
Testosterone	<0.1
Cortisone	<0.1
DHEA-S	<0.02
Cortisol	<0.02

9.3 Reproducibility

Intra-Assay

The intra-assay variation was determined by 20 replicate measurements of 5 saliva samples within one run. The within-assay variability is shown below:

Mean (pg/mL)	1328.7	650.1	293.8	186.9	23.3
SD (pg/mL)	75.9	30.4	15.8	10.6	1.8
CV (%)	5.7	4.7	5.4	5.7	7.6
n =	20	20	20	20	20

Inter-Assay

The inter-assay (between-day) variation was determined by duplicate measurements of 5 saliva samples over 10 days.

Mean (pg/mL)	55.5	565.4	1338.0	567.9	630.6
SD (pg/mL)	4.3	41.5	70.9	40.5	45.1
CV (%)	7.7	7.3	5.3	7.1	7.2
n =	20	20	20	20	20

Inter-Lot

The Inter-Lot (between-lot) variation was determined by triplicate measurements of five saliva samples in three different kit lots. The between lot variability is shown below:

Mean (pg/mL)	60.4	630.3	1574.8	360.3	635.1
SD (pg/mL)	4.3	37.6	64.1	28.7	33.4
CV%	7.1	6.0	4.1	8.0	5.3
n=	9	9	9	9	9

9.4 Recovery

Recovery of the IBL-AMERICA Progesterone ELISA was determined by adding increasing amounts of the analyte to five different saliva samples containing different amounts of endogenous analyte. Each sample (native and spiked) was assayed and analyte concentrations of the samples were calculated from the standard curve. The percentage recoveries were determined by comparing expected and measured values of the samples

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5
Concentration pg/mL	13.9	87.9	147.8	283.0	386.8
Average % recovery	96.4	103.2	100.3	97.3	102.2
Range of from % recovery to	91.1 101.6	101.1 104.8	94.9 103.8	92.7 106.6	97.8 105.7

9.5 Linearity

In total six saliva samples containing different amounts of analyte were serially diluted with Standard 0 and assayed with the IBL-AMERICA ELISA.

Three of these samples were serially diluted directly, and the other 3 samples at first were spiked with progesterone and then serially diluted up to 1:128.


The percentage recovery was calculated by comparing the expected and measured values for progesterone. A linearity of 3.8 – 4600 pg/mL has been identified as the usable range for this assay. Samples above this range must be diluted and re-run.

	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6
Concentration (pg/mL)	58.3	98.7	1073	2802	6000	5000
Average % Recovery	104.1	97.5	95.79	98.96	102.6	103.5
Range of from Recovery % to	99.9 106.7	89.9 100.8	91.2 105.5	93.6 107.3	99.1 106.5	94.4 107.7

9.6 Comparison Studies

A study was performed that evaluated 306 saliva samples collected from adult men and women ages 20 – 75. The samples were run on the IBL-AMERICA test and a commercially available LIA method to determine the concentration of free progesterone in the saliva samples. A correlation of 0.9373 and regression formula of $y = 0.912x + 6.066$ was obtained versus this method.

An additional study was performed using 101 saliva samples from adult men and women ages 20 - 75. These samples were compared to the LIA method, and yielded a correlation of $r = 0.9923$ with the following regression formula.



10 LIMITATIONS OF PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Blood contamination of more than 0.16% in saliva samples will affect results, and usually can be seen by eye.











Concentrations of Sodium Azide $\geq 0.02\%$ interferes in this assay and may lead to false results.











11 REFERENCES

1. Authors: Charles D. West, Damodar K. Mahajan, Virginia J. Chavré et.al.
Title: Simultaneous Measurement of Multiple Plasma Steroids by Radioimmunoassay Demonstrating Episodic Secretion.
Published in: Journal of Clinical Endocrinology and Metabolism **1973** (36 No.6) pages 1230 – 1236.
2. Authors: Ross F Vining, Robynne A McGingley, Richard G. Symons
Title: Hormones in Saliva: Mode of Entry and Consequent Implications for Clinical Interpretation
Published in: Clinical Chemistry **1983**, Vol 20 (10), pages 1752 - 1756
3. Authors: Jung K. Choe, Firyal S. Khan-Dawood, M. Yusoff Dawood
Title: Progesterone and Estradiol in the saliva and plasma during the menstrual cycle.
Published in: Amer. J. of Obstetrics & Gynecology **1983**, Vol 147 (5), pages 557 - 562
4. Authors: L. Cedard, A. Guichard, Y. Janssens, G. Tanguy, P. Boyer, JR Zorn
Title: Progesterone and Estradiol in saliva after in vitro fertilization and embryo transfer.
Published in: Fertility and Sterility **1987**, Vol 47 (2), pages 278 - 283
5. Authors: Ross F. Vining, Robynne A. McGinley
Title: The Measurement of Hormones in Saliva: Possibilities and Pitfalls
Published in: J. Steroid Biochem. **1987** (27 no.1-3) pages 81-94
6. Authors: Marguerite J. McNeely, Michael R. Soules
Title: The diagnosis of luteal phase deficiency: a critical review
Published in: Fertility and Sterility, **1988** (50 No.1) pages 1 – 15
7. Authors: Michael R. Soules, Donald K. Clifton, Robert A. Steiner et.al.
Title: The Corpus Luteum: Determinants of Progesterone Secretion in the Normal Menstrual Cycle.
Published in: Obstetrics and Gynaecology **1988** (71 no.5) pages 659 – 666
8. Authors: Taru Vuorento, Ilpo Huhtaniemi
Title: Daily levels of salivary Progesterone during menstrual cycle in adolescent girls.
Published in: Fertility and Sterility **1992** (58:4) pages 685 - 690.
9. Authors: M.M. Finn, J.P. Gosling, D.F. Tallon, S. Baynes, F.P. Meehan et.al.
Title: The frequency of salivary Progesterone sampling and the diagnosis of luteal phase insufficiency.
Published in: Gynecol. Endocrinol. **1992** 6 pages 127 - 134
10. Authors: T. Blom, A. Ojanotko-Harri, M. Laine, I. Huhtaniemi
Title: Metabolism of Progesterone and Testosterone in human parotid and submandibular salivary glands *in vitro*.
Published in: J Steroid Biochem Molec Biol **1993**, Vol. 44 No.1, pages 69-76
11. Authors: Torsten M. Delfs, Olaf G. Naether, Susann Klein, Freimut A. Leidenberger, Patrick Fottrell, Ralf C. Zimmermann
Title: 24-Hour profiles of salivary Progesterone
Published in: Fertility and Sterility **1994** (62 no.5) pages 960 – 966

12. Authors: Peter O'Leary, Peter Feddema, Katharine Chan, Mario Taranto, Sharon Evans
Title: Salivary, but not serum or urinary levels of Progesterone are elevated after topical application of Progesterone cream to pre- and postmenopausal women.
Published in: *Clinical Endocrinology* **2000**, Vol 53, pages 615 - 620.
13. Authors: M. Gröschl, M. Rauh, P. Schmid, H-G Dörr
Title: Relationship between salivary Progesterone, 17OH-Progesterone, and Cortisol levels throughout the normal menstrual cycle of healthy postmenarche girls.
Published in: *Fertility & Sterility* **2001**, Vol 76 (3), pages 615 - 617
14. Authors: M. Gröschl, M. Rauh, P. Schmid, HG Dörr
Title: Relationship between salivary Progesterone, 17OH-Progetserone, and Cortisol levels throughout the normal menstrual cycle of healthy postmenarcheal girls.
Published in: *Fertility and Sterility* **2001**, Vol 76 (3), pages 615 - 617
15. Authors: T. Vimpeli, H. Tinkanen, H. Huhtala, L. Rönberg, E. Kujansuu
Title: Salivary and serum Progesterone concentrations during two luteal support regimens used in in-vitro fertilization treatments
Published in: *Fertility & Sterility* **2001**, Vol.4, pages 847-848
16. Authors: Nike Stikkelbroeck, C Sweep, Didi Braat, Ad Hermus, Barto Otten
Title: Monitoring of menstrual cycles, ovulation, and adrenal suppression by saliva sampling in female patients with 21-hydroxylase deficiency.
Published in: *Fertility and Sterility* **2003**, Vol 80 (4), pages 1030 - 1036
17. Authors: M Gröschl, M Rauh, HG Dörr
Title: Circadian rhythm of salivary Cortisol, 17OHP, and Progesterone in healthy children.
Published in: *Clin Chem.* **2003**, Vol. 49 (10), pages 1688-1691
18. Authors: RT Chatterton, ET Mateo, N Hou, AW Rademaker, S Acharya, VC Jordan
Title: Characteristics of salivary profiles of Estradiol and Progesterone in postmenopausal women.
Published in: *Journal of Endocrinology* **2005**, Vol. 186, pages 77-84
19. Author: John G. Lewis
Title: Steroid Analysis in Saliva: An overview (Review article)
Published in: *Clin Biochem Rev.*, August **2006**, Vol.27, pages 139-146
20. Authors: L.J. Bond, E.T. Vella, Y. Kiparissis, K.E. Wynne-Edwards
Title: Anthropometry and body composition do not predict bioavailable Androgen or Progesterone concentration in adolescent girls.
Published in: *Am J Hum Biol* **2006**, Vol.18, pages 639-653
21. Authors: V.J. Vitzthum, H. Spielvogel, J. Thornburg, B. West
Title: A prospective study of early pregnancy loss in humans
Published in: *Fertil Steril* **2006**, Vol.86, pages 373-379
22. Authors: MM Wirth, EA Meier, BL Fredrickson, OC Schultheiss
Title: Relationship between salivary Cortisol and Progesterone levels in humans.
Published in: *Biological Psychology* **2007**, Vol.74, pages 104-107
23. Authors: A. Nunez-de-la-Mora, R.T. Chatterton, O.A. Choudhury, G.R. Bentley
Title: Childhood Conditions Influence Adult Progesterone Levels
Published in: *www.PLOSmedicine.org* **2007**, Vol.4 (5) pages 1-8
24. Authors: B.K. Gandara, L. Leresche, L. Mancl
Title: Patterns of salivary Estradiol and Progesterone across the menstrual cycle.
Published in: *Ann. NY Acad. Sci.* **2007**, Vol. 1098, pages 446-450

12 SYMBOLS USED WITH IBL-AMERICA ELISA'S

Symbol	English	Deutsch	Français	Espanol	Italiano
	European Conformity	CE-Konformitätskennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di Cat.
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributor	Vertreiber	Distributeur	Distribuidor	Distributore
<i>Content</i>	Content	Inhalt	Contenu	Contenido	Contenuto
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Numéro	Volumen/Número	Volume/Quantità
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Plaques de microtitration	Placas multipocillo	Micropozzetti
<i>Antiserum</i>	Antiserum	Antiserum	Antisérum	Antisero	Antisiero
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Conjugué enzymatique	Conjugado enzimático	Tracciante enzimatico
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complexe enzymatique	Complex enzimático	Complesso enzimatico
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Solution substrat	Solución de sustrato	Soluzione di substrato
<i>Stop Solution</i>	Stop Solution	Stopplösung	Solution d'arrêt	Solución de parada	Soluzione d'arresto
<i>Zero Standard</i>	Zero Standard	Nullstandard	Zero Standard	Estándar cero	Standard zero
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Contrôle	Control	Controllo
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampon d'essai	Tampón de ensayo	Tampone del test
<i>Wash Solution</i>	Wash Solution	Waschlösung	Solution de lavage	Solución de lavado	Soluzione di lavaggio
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)
<i>1 N HCl</i>	1 N HCl	1 N HCl	1N HCl	1 N HCl	
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Solution pour dilution de l'échantillon	Solución para dilución de la muestra	Diluyente dei campioni
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungsmedium	Solution pour dilution du conjugué	Solución para dilución del conjugado	Diluyente del tracciante

Symbol	Portugues	Dansk	Svenska	Ελληνικά
	Conformidade com as normas europeias	Europaeisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση
	Consulte as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Εγχειρίδιο χρήστη
	Diagnóstico in vitro	In vitro diagnostik	Diagnostik in vitro	in vitro διαγνωστικό
	Catálogo n.º	Katalognummer	Katalog nummer	Αριθμός καταλόγου
				
	No do lote	Lot nummer	Batch-nummer	Αριθμός Παρτίδος
		Indeholder tilstrækkeligt til "n" test	Innehåller tillräckligt till "n" tester	Περιεχόμενο επαρκές για «n» εξετάσεις
	Temperatura de conservação	Opbevaringstemperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης
	Prazo de validade	Udløbsdato	Bäst före datum	Ημερομηνία λήξης
	Fabricante	Producent	Tillverkare	Κατασκευαστής
<i>Distributed by</i>				
<i>Content</i>	Conteúdo	Indhold	Innehåll	Περιεχόμενο
<i>Volume/No.</i>	Volume/Número	Volumen/antal	Volym/antal	Όγκος/αριθ..
<i>Microtiterwells</i>	Alvéolos de microtitulação	Mikrotiterbrønde	Brunnar i Mikrotiterplatta	Πηγαδάκια Μικροπιλοδοτήσεως
<i>Antiserum</i>	Anti-soro	Antiserum	Antiserum	Αντιπρός
<i>Enzyme Conjugate</i>	Conjugado enzimático	Enzymkonjugat	Enzymkonjugat	Συζευγμένο ενζυμο
<i>Enzyme Complex</i>	Complexo enzimático	Enzymkompleks	Enzymkomplex	Σύμπλοκο ενζύμου
<i>Substrate Solution</i>	Solução de substrato	Substratopløsning	Substratlösning	Διάλυμα υποστρώματος
<i>Stop Solution</i>	Solução de paragem	Stopopløsning	Stopp lösning	Διάλυμα τερματισμού
<i>Zero Standard</i>	Padrão zero	Standard 0	Standard 0	Πρότυπο Μηδέν
<i>Standard</i>	Calibrador	Standard	Standard	Πρότυπα
<i>Control</i>	Controlo	Kontrol	Kontroll	Έλεγχος
<i>Assay Buffer</i>	Tampão de teste	Assay buffer	Assay Buffer	Ρυθμιστικό Διάλυμα Εξέτασης
<i>Wash Solution</i>	Solução de lavagem	Vaskebuffer	Tvätt lösning	Διάλυμα πλύσεως
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl	1 N HCl	1 N HCl
<i>Sample Diluent</i>				
<i>Conjugate Diluent</i>				