


# MedFrontier® FGF23 Intact

## Intact Fibroblast Growth Factor 23

The better choice to measure the biologically active form of FGF23. The MedFrontier® FGF23 Intact assay kit is a reliable, reproducible, and accurate serum blood test delivering optimum results quickly.



### ASSAY KEY FEATURES

- ➔ Measures intact FGF23, the full-length bioactive form of FGF23
- ➔ Superior dynamic range of 10-3000 pg/mL
- ➔ Excellent reproducibility, CVs below 5%
- ➔ Microtiter plate assay that every lab can perform
- ➔ CLEIA technology
- ➔ Only 20 µL serum needed



### AREAS OF FGF23 RESEARCH\*



#### Osteomalacia

Tumor-induced osteomalacia (TIO)



#### Kidney Disease

Chronic kidney disease (CKD), Acute kidney injury (AKI)



#### Rickets

X-linked hypophosphatemic rickets (XLH), Autosomal dominant hypophosphatemic rickets (ADHR)

#### Cardiovascular Disease

Left ventricular hypertrophy (LVH)



#### Vitamin D Metabolism

Control of 1,25(OH)<sub>2</sub> vitamin D homeostasis



#### Mineral Bone Disease

Hyperphosphatemia, Iron deficiency



## CHEMILUMINESCENCE ENZYME IMMUNOASSAY (CLEIA)

- ➡ CLEIA known for better sensitivity, specificity, reliability and reproducibility
- ➡ Serum intact FGF23 concentration determined using CLEIA
- ➡ Antibody-coated microtiter plate configuration for easy use
- ➡ Use of two mouse monoclonal antibodies measures only intact FGF23

## ASSAY PERFORMANCE / ACCURACY & REPRODUCIBILITY

Intact FGF23 Concentration in Serum (pg/mL)	Expected Value (pg/mL)	Measured Value (pg/mL)	Accuracy (%)	Reproducibility (CV : N=10)
Low	22.8	22.7	99.6	3.1%
Medium	249.9	271.6	108.7	1.4%
High	1828.8	1946.3	106.4	2.4%

## INTACT FGF23 AND VITAMIN D REGULATE PHOSPHATE HOMEOSTASIS

Intact FGF23 is the biologically active protein encoded by the FGF23 gene responsible for phosphate and Vitamin D metabolism.

