

Code No. 17166

# Rat GRO/CINC-2 $\beta$ Assay Kit - IBL

#### INTRODUCTION

Cytokine-induced neutrophil chemoattractant-1 (CINC-1) was originally purified from media conditioned by IL-1 $\beta$  stimulated rat kidney epithelioid cells (NRK-52E). Watanabe's group at Toyama Medical and Pharmaceutical University identified amino acid sequence that encodes for rat CINC-1 in 1989. CINC-1 is a member of the alpha (CXC) subfamily of chemokines. Three additional rat CXC chemokines (CINC-2  $\alpha$ , CINC-2  $\beta$ , CINC-3/MIP-2) have been identified. The protein sequence of CINC-1 is 63 - 67% identical to that of CINC-2 $\alpha$ , CINC-2 $\beta$ , CINC-3/MIP-2. In addition, GRO  $\alpha$ , GRO  $\beta$  and GRO $\gamma$  is sharing 68%, 71% and 69%, identity with CINC-1. This has been suggested that CINCs are the rat counterpart of human GROs.

# PRINCIPLE

This kit is a solid phase sandwich ELISA using 2 kinds of high specific antibodies. Tetra Methyl Benzidine (TMB) is used as coloring agent (Chromogen). The strength of coloring is in proportion to the quantities of Rat GRO/CINC-2  $\beta$ .

### MEASUREMENT RANGE

9.38 ~ 600 pg/mL

### INTENDED USE

- The IBL's GRO/CINC-2 β EIA Kit is a complete kit for the quantitative determination of GRO/CINC-2 β in serum and supernatant of cell culture media.
- Both recombinant and native forms of GRO/CINC-2 β can be detected with the kit.

#### **KIT COMPONENT**

| 1 | Precoated plate        | : Anti- Rat GRO/CINC-2 $\beta$ Rabbit IgG Affinity Purify                 | 96Well x 1 |
|---|------------------------|---|------------|
| 2 | Labeled antibody       | : HRP conjugated Anti-Rat GRO/CINC-2     Rabbit IgG Fab' Affinity Purifiy | 10.5mL x 1 |
| 3 | EIA buffer             | : 1% BSA, 0.05% Tween 20 in PBS   | 30mL x 1   |
| 4 | Standard               | : Recombinant Rat GRO/CINC-2 β  | 0.5mL x 1  |
| 5 | Substrate Buffer       | : ① White cap /Solution for (7) Chromogen TMB                             | 5mL x 1    |
|   |                        | : ② Pink cap /0.01% H <sub>2</sub> O <sub>2</sub> solution                | 5.5mL x 1  |
| 6 | Stop solution          | : 1N H <sub>2</sub> SO <sub>4</sub>                                       | 11mL x 1   |
| 7 | Chromogen              | : Tetra Methyl Benzidine (TMB)  | 1mg x 2    |
| 8 | Wash buffer Conc.      | : 0.05% Tween20 in phosphate buffer (X40)                                 | 50mL x 1   |
| 9 | Solution for Labeled a | antibody : 1% BSA, 0.05%Tween20 in PBS                                    | 10.5mL x 1 |
|   |                        |   |            |

#### **OPERATION MANUAL**

#### 1. Materials needed but not supplied

- 2. Preparation
  - 1) Preparation of wash buffer

"8, Wash buffer Conc." is a concentrated (X40) buffer. The temperature of "8, Wash buffer Conc." shall be adjusted to room temperature and then, mix it gently and completely before use. Dilute 50mL of "8, Wash buffer Conc." with 1,950mL of distilled water and mix it. This is the wash buffer for use. This prepared wash buffer shall be stored in refrigerator and used within 2 weeks after dilution.

2) Solution of labeled antibody

Put 10.5mL of "9, Solution for Labeled antibody" into the vial of "2, Labeled antibody" and place it at room temperature for 5 minutes. Then mix it gently and completely.

This operation should be done just before the application of labeled antibody.

3) Preparation of TMB buffer

Put two tablets of TMB ("7, Chromogen") into the vial of "5-①, Substrate buffer" and mix it gently and completely. This prepared mixture is a white colored suspension. Then put "5-②, Substrate buffer" into this mixture and mix it gently and completely. This is TMB buffer for use. This operation should be done just before the application of TMB buffer.

37.5 pg/mL 18.75 pg/mL 9.38 pg/mL 0 pg/mL (Test Sample Blank)

Put 230  $\mu$  L of Standard solution into tube–1 and mix it gently. Then, put 230  $\mu$  L of tube-1mixture into tube-2. Dilute two times standard solution in series to set up 7 points of diluted standard between 600 pg/mL and 9.38 pg/mL. Tube-8 is the test sample blank as 0 pg/mL. See following picture.



6) Dilution of test sample

Tube-5

Tube-6

Tube-7

Tube-8

Test sample may be diluted with "3, EIA buffer" if the need arises. If the concentration of Rat GRO/CINC-2 $\beta$  in samples may not be estimated in advance, the pre-assay with several different dilutions will be recommended to determine the proper dilution of samples.

## 3. Measurement procedure

All reagents shall be brought to room temperature approximately 30 minutes before use. Then mix it gently and completely before use. Confirm no change in quality of the reagents. Standard curve shall be prepared simultaneously with the measurement of test samples.

|  | Test Sample            | Standard                                     | Test Sample<br>Blank              | Reagent<br>Blank      |  |  |
|--|------------------------|--|-----------------------------------|-----------------------|--|--|
| Reagents   | Test sample<br>100 μ L | Diluted<br>standard<br>(Tube 1~7)<br>100 µ L | EIA buffer<br>(Tube-8)<br>100 μ L | EIA buffer<br>100 μ L |  |  |
| Incubation for 1 hour at 37°C with plate lid             |                        |  |                                   |                       |  |  |
|  | Washing 7 times        |  |                                   |                       |  |  |
| Labeled<br>Antibody                                      | 100 µ L                | 100 μ L                                      | 100 μ L                           | -                     |  |  |
| Incubation for 30minutes at 37°C with plate lid          |                        |  |                                   |                       |  |  |
| Washing 9 times  |                        |  |                                   |                       |  |  |
| TMB buffer   | 100 μ L                | 100 μ L                                      | 100 μ L                           | 100 μ L               |  |  |
| Incubation for 30 minutes at room temperature (shielded) |                        |  |                                   |                       |  |  |
| Stop<br>solution   | 100 µ L                | 100 µ L                                      | 100 µ L                           | 100 µ L               |  |  |
| Read the pla   | te at 450nm with       | nin 30 minutes af                            | ter application of                | Stop solution.        |  |  |

- 1) Determine wells for reagent blank. Put  $100 \,\mu$  L each of "3, EIA buffer" into the wells.
- 2) Determine wells for test sample blank, test sample and diluted standard. Then, put  $100 \,\mu$  L each of test sample blank (tube-8), test sample and dilutions of standard (tube-1~7) into the appropriate wells.
- Incubate the precoated plate for 1 hour at 37℃ after covering it with plate lid.
- 4) Wash each well of the precoated plate vigorously with wash buffer using washing bottle. Then, fill each well with wash buffer and place the precoated plate for 15~30 seconds. Remove wash buffer completely from the precoated plate by snapping.

This procedure must be repeated more than 7 times.

Then, remove the remaining liquid from all wells completely by snapping the precoated plate onto paper towel.

In case of using plate washer, after 4 times washing with plate washer, washing with above washing bottle must be repeated 3 times.

- 5) Pipette  $100 \,\mu$  L of labeled antibody solution into the wells of test samples, diluted standard and test sample blank.
- 6) Incubate the precoated plate for 30 minutes at 37°C after covering it with

#### 4) Preparation of Standard

Put just 0.5mL of distilled water into the vial of "4, Standard" and mix it gently and completely. This solution is 1,200 pg/mL Rat GRO/CINC-2  $\beta$  standard.

### 5) Dilution of Standard

Prepare 8 tubes for dilution of "4, Standard". Put  $230 \,\mu$ L each of "3, EIA buffer" into the tube.

Specify the following concentration of each tube.

| Tube-1 | 600 pg/mL |
|--------|-----------|
| Tube-2 | 300 pg/mL |
| Tube-3 | 150 pg/mL |
| Tube-4 | 75 pg/mL  |

#### plate lid.

- 7) Wash the precoated plate 9 times in the same manner above 4).
- 8) Pipette  $100 \,\mu$  L TMB buffer into the wells.
- 9) Incubate the precoated plate for 30 minutes at room temperature in the dark. The liquid will turn blue by the addition of TMB buffer.
- 10) Pipette  $100 \,\mu$  L of "6, Stop solution" into the wells. Mix the liquid by tapping the side of precoated plate. The liquid will turn yellow by the addition of "6, Stop solution".
- Remove any dirt or drop of water on the bottom of the precoated plate and confirm there is no bubble on the surface of the liquid. Then, run the plate reader and conduct measurement at 450nm.
   The measurement shall be done within 30minutes after the addition of "6.

The measurement shall be done within 30minutes after the addition of "6, Stop solution".

IBL Co., Ltd. 5-1 Aramachi, Takasaki-shi, Gunma, 370-0831, JAPAN URL: <u>http://www.ibl-japan.co.jp</u> E-mail: <u>do-ibl@ibl-japan.co.jp</u> TEL: +81 (0)27-310-8040 FAX : +81 (0)27-310-8045



### SPECIAL ATTENTION

- Test samples should be measured soon after the collection. In case of the storage of test samples, they should be stored under frozen conditions and do not repeat freeze/thaw cycles. Thaw the test samples at low temperature and mix them completely before measurement.
- 2. Test samples should be diluted with "3, EIA buffer", if the need arises.
- 3. The measurement of test samples and standard in duplicate is recommended.
- 4. Use test samples in neutral pH range. The contaminations of organic solvent may affect the measurement.
- 5. Use only wash buffer contained in this kit for washing the precoated plate. Insufficient washing may lead to the failure in measurement.
- Remove the wash buffer completely by tapping the precoated plate on paper towel.
   Do not wipe wells with paper towel.
- 7. TMB buffer should be prepared just before use. Do not use the colored TMB buffer. TMB tablets must be put to only "5-①, Substrate buffer". TMB tablets can not be dissolved in "5-②, Substrate buffer". TMB buffer should be stored in the dark due to its sensitivity against light. TMB buffer should be avoided contact with metals.
- 8. Measurement should be done within 30 minutes after addition of "6, Stop solution".
- 9. Storage of HRP conjugated antibody is not recommended. However, if the HRP conjugates do not use at one time, please store it at below -20℃.

### CALCULATION OF TEST RESULT

Subtract the absorbance of test sample blank from all data, including standards and unknown samples before plotting. Plot the subtracted absorbance of the standards against the standard concentration on log-log graph paper. Draw the best smooth curve through these points to construct the standard curve. Read the concentration for unknown samples from the standard curve.

Example of standard curve



\* The typical standard curve is shown above. This curve can not be used to derive test results. Please run a standard curve for each assay.

#### PERFORMANCE CHARACTERISTICS

#### 1. Titer Assay

| Specimen  | Titer<br>(X) | Measurement<br>Value (pg/mL) | Theoretical<br>Value (pg/mL) | %     |
|-----------|--------------|------------------------------|------------------------------|-------|
|           | 4            | 276.1                        | 300                          | 92.0  |
|           | 8            | 143.3                        | 150                          | 95.5  |
| RPMI-1640 | 16           | 68.1                         | 75                           | 90.8  |
|           | 32           | 33.8                         | 37.5                         | 90.1  |
|           | 64           | 19.2                         | 18.8                         | 102.1 |
|           | 4            | 315.4                        | 300                          | 105.1 |
|           | 8            | 165.6                        | 150                          | 110.4 |
| Rat Serum | 16           | 84.9                         | 75                           | 113.2 |
|           | 32           | 38.1                         | 37.5                         | 101.6 |
|           | 64           | 18.9                         | 18.8                         | 100.5 |

3. Inter - Assay

| Measurement<br>Value (pg/mL) | SD value | CV value<br>(%) | n |
|------------------------------|----------|-----------------|---|
| 285.3                        | 7.8      | 2.7             | 7 |
| 51.6                         | 3.4      | 6.6             | 7 |
| 17.8                         | 0.3      | 1.7             | 7 |

4. Intra - Assay

| Measurement<br>Value (pg/mL) | SD value | CV value<br>(%) | n |
|------------------------------|----------|-----------------|---|
| 253.0                        | 7.4      | 2.9             | 3 |
| 37.2                         | 2.8      | 7.5             | 3 |

<sup>5.</sup> Specificity

| Compound        | Cross Reactivity |
|-----------------|------------------|
| Rat GRO/CINC-2β | 100.0%           |
| Rat GRO/CINC-2α | ≦0.1%            |
| Rat GRO/CINC-1  | ≦0.1%            |
| Rat GRO/CINC-3  | ≦0.1%            |
| Rat MCP-1       | ≦0.1%            |
| Rat Rantes      | ≦0.1%            |
| Rat MIP-1 α     | ≦0.1%            |
| Rat IL-1 β      | ≦0.1%            |

#### PRECAUTION FOR INTENDED USE AND/OR HANDLING

- 1. All reagents should be stored at 2~8°C. All reagents shall be brought to room temperature approximately 30 minutes before use.
- 2. HRP conjugated antibody and standard are lyophilized products. Be careful to open these vials.
- 3. "6, Stop solution" is a strong acid substance. Therefore, be careful not to contact your skin and clothes with "6, Stop solution" and pay attention to the disposal of "6, Stop solution".
- 4. "1, Precoated plate" and "4, Standard" contain sodium azide. Therefore, dispose these materials after diluting them with large quantity of water to avoid the production of explosive metallic azide.
- 5. Cetylpyridinium chloride is used as preservative for "2, Labeled antibody", "3, EIA buffer" and "9, Solution for Labeled antibody".
- 6. Wash hands after handling reagents.
- 7. Do not mix the reagents with the reagents from different lot or different kit.
- 8. Do not use the reagents expired.
- 9. This kit is for research purpose only. Do not use for clinical diagnosis.

#### STORAGE AND THE TERM OF VALIDITY

| Storage Condition    | : 2 ~ 8℃                                     |
|----------------------|--|
| The term of validity | : 6 months                                   |
|                      | (The expiry date is specified in outer box.) |

#### REFERENCES

 Makita, H. et al. Effect of anti-macrophage migration inhibitory factor antibody on lipopoly saccharide-induced pulmonary neutrophil accumulation. AMERICAN JOURNAL OF RESPIRATORY AND CRITICAL CARE MEDICINE., 1998: 158 (2), 573-579

| Version |         |                                   |
|---------|---------|-----------------------------------|
| 990901  | Revised |                                   |
| 000401  | Revised | (Wash buffer conc)                |
| 000601  | Revised | (Up-dated layout)                 |
| 001024  | Revised | (Up-dated layout)                 |
| 010411  | Revised | (Additional performance of a kit) |
| 030415  | Revised |                                   |
| 040201  | Revised | (Up-dated layout)                 |
| 040501  | Revised |                                   |



Distributed By: **IBL-America, Inc.** 8201 Central Ave NE, Suite P Minneapolis, MN 55432, USA <u>info@ibl-america.com</u> (888) 523 1246

# 2. Added Recovery Assay

| Specimen      | Theoretical<br>Value(pg/mL) | Measurement<br>Value (pg/mL) | %     |
|---------------|-----------------------------|------------------------------|-------|
|               | 300                         | 278.1                        | 92.7  |
|               | 150                         | 149.8                        | 99.9  |
| 10% FCS added | 75                          | 74.6                         | 99.5  |
| RPMI-1640     | 37.5                        | 37.1                         | 98.9  |
|               | 18.8                        | 19.2                         | 102.1 |
|               | 9.4                         | 9.4                          | 100.0 |

 IBL Co., Ltd.
 5-1 Aramachi, Takasaki-shi, Gunma, 370-0831, JAPAN
 TEL: +81 (0)27-310-8040
 FAX : +81 (0)27-310-8045

 URL:
 <u>http://www.ibl-japan.co.jp</u>
 E-mail: <u>do-ibl@ibl-japan.co.jp</u>
 FAX : +81 (0)27-310-8040
 FAX : +81 (0)27-310-8045