

Code No. 27770

## Rat LRG Assay Kit - IBL

### INTRODUCTION

Leucine rich alpha2-glycoprotein (LRG) is one of leucine-rich repeat family proteins. LRG has been shown to be involved in protein-protein interaction, signal transduction, and cell adhesion and development. Recently, it was reported that serum levels of LRG were significantly elevated in Rheumatoid arthritis (RA) patients compared to healthy controls and decreased after anti-TNF therapy (ref. 1). And another report indicate that LRG levels in cerebrospinal fluid in iNPH (idiopathic Normal Pressure Hydrocephalus) patients are remarkably elevated in compared with controls (ref. 2) and it may be possible to differentiate iNPH patients from others showing similar symptoms like Alzheimer's disease (ref. 3).

This product is an ELISA kit for measuring of Rat LRG.

### PRINCIPLE

This kit is a solid phase sandwich ELISA using 2 kinds of highly specific antibodies. Tetra Methyl Benzidine (TMB) is used as a coloring agent (Chromogen). The strength of coloring is proportional to the quantities of Rat LRG.

### MEASUREMENT RANGE

6.25 - 400 ng/mL

### INTENDED USE

For research use only, not for use in diagnostic procedures.

- This assay kit is capable for the quantitative determination of rat LRG in serum, EDTA-plasma, urine and cell culture supernatant.
- The guide line of dilution for serum or plasma samples is 4 or 8-fold with "4, EIA buffer" in this kit. LRG in blood samples cannot be measured reasonably by doubling dilution.
- The guide line of dilution for urine samples is about 8-fold with "4, EIA buffer."

### KIT COMPONENT

1	Precoated plate	:		
	Anti-Rat LRG (191) Rabbit IgG Affinity Purify		96Well x 1	
2	Labeled antibody Conc.	:		
	(30X) HRP conjugated Anti-Rat LRG (312) Rabbit IgG Fab' Affinity Purify		0.4mL x 1	
3	Standard	:	Recombinant Rat LRG	0.5mL x 2
4	EIA buffer*			30mL x 1
5	Solution for Labeled antibody*			12mL x 1
6	Chromogen	:	TMB solution	15mL x 1
7	Stop solution*			12mL x 1
8	Wash buffer Conc.*			50mL x 1

### OPERATION MANUAL

#### 1. Materials needed but not supplied

- Plate reader (450nm)
- Graduated cylinder and beaker
- Refrigerator (as 4°C)
- Paper towel
- Incubator (37°C ± 1°C)
- Washing bottle for precoated plate
- Disposable test tube for "2, Labeled antibody Conc." and "6, Chromogen"
- Micropipette and tip
- Deionized water
- Graph paper (log/log)
- Tube for dilution of Standard

#### 2. Preparation


- 1) Preparation of wash buffer  
"8, Wash buffer Conc." is a concentrated (40X) buffer. Adjust the temperature of "8, Washing buffer Conc." to room temperature and then, mix it gently and completely before use. Dilute 50 mL of "8, Wash buffer Conc." with 1,950 mL of deionized water and mix it. This is the wash buffer for use. This prepared wash buffer shall be stored in refrigerator and used within 2 weeks after dilution.
- 2) Preparation of Labeled antibody  
"2, Labeled antibody Conc." is a concentrated (30X). Dilute "2, Labeled antibody Conc." with "5, Solution for Labeled antibody" in 30 times according to required quantity into a disposable test tube. Use this resulting solution as Labeled antibody.  
Example)  
In case you use one strip (8 well), the required quantity of Labeled antibody is 800 µL. (Dilute 30 µL of "2, Labeled antibody Conc." with 870 µL of "5, Solution for Labeled antibody" and mix it. And use the resulting solution by 100 µL in each well.)  
This operation should be done just before applying labeled antibody.  
The remaining "2, Labeled antibody Conc." should be stored at 4°C in firmly sealed vial.
- 3) Preparation of Standard  
Put just 0.5 mL of deionized water into the vial of "3, Standard" and mix it gently and completely. This solution is 800 ng/mL Rat LRG standard.
- 4) Dilution of Standard  
Prepare 8 tubes for dilution of "3, Standard". Put 230 µL each of "4, EIA buffer" into the tube.  
Specify the following concentration of each tube."

Tube-1	400 ng/mL
Tube-2	200 ng/mL
Tube-3	100 ng/mL
Tube-4	50 ng/mL
Tube-5	25 ng/mL
Tube-6	12.5 ng/mL
Tube-7	6.25 ng/mL
Tube-8	0 ng/mL (Test Sample Blank)

Put 230 µL of Standard solution into tube-1 and mix it gently. Then, put 230 µL of tube-1 mixture into tube-2. Dilute two times standard solution in series to set up 7

points of diluted standard between 400 ng/mL and 6.25 ng/mL. Tube-8 is the test sample blank as 0 ng/mL.

See following picture.



#### 5) Dilution of test sample

Test samples should be diluted with "4, EIA buffer" as necessary.

In the case of blood (serum or EDTA-plasma) samples, the specimens have to be diluted 4-fold or more because "4, EIA buffer" is an essential factor for the assay.

If the concentration of Rat LRG in samples may not be estimated in advance, the pre-assay with several different dilutions will be recommended to determine the proper dilution of samples.

### 3. Measurement procedure

All reagents shall be brought to room temperature approximately 30 minutes before use. Then mix it gently and completely before use. Make sure of no change in quality of the reagents. Standard curve shall be prepared simultaneously with the measurement of test samples.

Reagents	Test Sample	Standard	Test Sample Blank	Reagent Blank
	Test sample 100 µL	Diluted standard (Tube 1-7) 100 µL	EIA buffer (Tube-8) 100 µL	EIA buffer 100 µL
Incubation for 90 minutes at 37°C with plate lid				
4 times (wash buffer more than 350 µL)*				
Labeled Antibody	100 µL	100 µL	100 µL	-
Incubation for 30 minutes at 4°C with plate lid				
5 times (wash buffer more than 350 µL)*				
Chromogen	100 µL	100 µL	100 µL	100 µL
Incubation for 30 minutes at room temperature (shielded)				
Stop solution	100 µL	100 µL	100 µL	100 µL
Read the plate at 450nm against a Reagent Blank within 30 minutes after addition of Stop solution.				

- 1) Determine wells for reagent blank. Put 100 µL each of "4, EIA buffer" into the wells.
- 2) Determine wells for test sample blank, test sample and diluted standard. Then, put 100 µL each of test sample blank (tube-8), test sample and dilutions of standard (tube-1-7) into the appropriate wells.
- 3) Incubate the precoated plate for 90 minutes at 37°C after covering it with plate lid.
- 4) Wash the plate with the prepared wash buffer and remove all liquid.\*
- 5) Pipette 100 µL of labeled antibody solution into the wells of test samples, diluted standard and test sample blank.
- 6) Incubate the precoated plate for 30 minutes at 4°C after covering it with plate lid.
- 7) Wash the plate with the prepared wash buffer and remove all liquid.\*
- 8) Take the required quantity of "6, Chromogen" into a disposable test tube. Then, pipette 100 µL from the test tube into every well. Please do not return the rest of used chromogen in the test tube into "6, Chromogen" bottle in order to avoid contamination.
- 9) Incubate the precoated plate for 30 minutes at room temperature in the dark. The solution of Chromogen will turn blue.
- 10) Add 100 µL of "7, Stop solution" to all wells. Mix the solution by tapping the side of precoated plate. The solution will turn yellow by addition of "7, Stop solution".
- 11) Remove any dirt or drop of water on the bottom of the precoated plate and confirm there is no bubble on the surface of the solution. Then, run the plate reader and conduct measurement at 450 nm against a reagent blank. The measurement shall be done within 30 minutes after addition of "7, Stop solution".

### SPECIAL ATTENTION


- 1) Test samples should be measured soon after collection. For the storage of test samples, store them frozen and do not repeat freeze/thaw cycles. Thaw the test samples at a low temperature and mix them completely before measurement.
- 2) Test samples should be diluted with "4, EIA buffer", as the need arises.
- 3) Duplicate measurement of test samples and standard is recommended.
- 4) Use test samples in neutral pH range. The contaminations of organic solvent may affect the measurement.
- 5) Use only wash buffer in this kit for washing the precoated plate. Insufficient washing may lead to the failure in measurement.
- 6) Remove the wash buffer completely by tapping the precoated plate on paper towel. Do not wipe wells with paper towel.
- 7) "6, Chromogen" should be stored in the dark due to its sensitivity against light. Avoid contact of Chromogen with metals.
- 8) Measurement should be done within 30 minutes after addition of "7, Stop solution".

## CALCULATION OF TEST RESULT

Subtract the absorbance of test sample blank from all data, including standards and unknown samples before plotting. Plot the subtracted absorbance of the standards against the standard concentration on log-log graph paper. Draw the best smooth curve through these points to construct the standard curve. Read the concentration for unknown samples from the standard curve.

Example of standard curve





Conc. (ng/mL)	Absorbance (450nm)
400	2.933
200	2.022
100	1.353
50	0.887
25	0.550
12.5	0.327
6.25	0.225
0 (Test Sample Blank)	0.102



\* The typical standard curve is shown above. This curve can not be used to derive test results. Please run a standard curve for each assay.

## PERFORMANCE CHARACTERISTICS

### 1. Dilution linearity



### 2. Added Recovery Assay

Specimen	Additive Amount (ng/mL)	Theoretical Value (ng/mL)	Measured Value (ng/mL)	%
Rat Serum (SD) (x4)	152.90	200.71	166.95	83.2
	84.88	132.70	131.15	98.8
Rat Plasma (EDTA) (SD) (x4)	152.90	228.29	191.24	83.8
	84.88	160.27	170.57	106.4
Rat Urine (SD) (x8)	60.00	83.49	92.10	110.3
	30.00	53.49	57.68	107.8
	15.00	38.49	38.59	100.3
Medium with 10% FBS (x4)	152.90	152.90	128.80	84.2
	84.88	84.88	71.90	84.7
	43.48	43.48	39.21	90.2

### 3. Intra - Assay

Mean Value (ng/mL)	SD (ng/mL)	CV (%)	n
84.20	4.70	5.6	24
44.52	2.96	6.6	24
7.87	0.55	7.0	24

### 4. Inter - Assay

Mean Value (ng/mL)	SD (ng/mL)	CV (%)	n
81.16	4.14	5.1	10
42.96	2.94	6.8	10
7.71	0.90	11.7	10

### 5. Sensitivity

1.68 ng/mL

The sensitivity for this kit was determined using the guidelines under the National Committee for Clinical Laboratory Standards (NCCLS) Evaluation Protocols. (National Committee for Clinical Laboratory Standards Evaluation Protocols, SC1, (1989) Villanova, PA: NCCLS.)

### PRECAUTION FOR INTENDED USE AND/OR HANDLING

- All reagents should be stored at 2 - 8°C. All reagents shall be brought to room temperature approximately 30 minutes before use.
- "3, Standard" is lyophilized products. Be careful to open this vial.
- "7, Stop solution" is a strong acid substance. Therefore, be careful not to have your skin and clothes contact "7, Stop solution" and pay attention to the disposal of "7, Stop solution".
- Dispose used materials after rinsing them with large quantity of water.
- Precipitation may occur in "2, Labeled antibody Conc.", "4, EIA buffer" or "8, Wash buffer Conc.", however, there is no problem in the performance.
- Wash hands after handling reagents.
- Do not mix the reagents with the reagents from a different lot or kit.
- Do not use expired reagents.
- This kit is for research purpose only. Do not use for clinical diagnosis.

### STORAGE AND THE TERM OF VALIDITY

Storage Condition : 2 - 8°C

The expiry date is specified on outer box.

### REFERENCE

- Serada S, Fujimoto M, Ogata A, Terabe F, Hirano T, Iijima H, Shinzaki S, Nishikawa T, Ohkawara T, Iwahori K, Ohguro N, Kishimoto T, Naka T. iTRAQ-based proteomic identification of leucine-rich alpha-2 glycoprotein as a novel inflammatory biomarker in autoimmune diseases. *Ann Rheum Dis.* 2010 Apr;69(4):770-4.
- Li X, Miyajima M, Mineki R, Taka H, Murayama K, Arai H. Analysis of potential diagnostic biomarkers in cerebrospinal fluid of idiopathic normal pressure hydrocephalus by proteomics. *Acta Neurochir (Wien).* 2006 Aug;148(8): 859-64.
- Nakajima M, et. al/Leucine-rich α-2-glycoprotein is a marker for idiopathic normal pressure hydrocephalus. *Acta Neurochir (Wien).* 2011 Jun;153(6): 1339-46.

Version 2.

February 2017 \*

Made in Japan.

### IBL LRG Product Lines:

Code No.	Name	Volume
27770	Rat LRG Assay Kit -IBL	96 Well
27769	Human LRG Assay Kit -IBL	96 Well
27785	Mouse LRG Assay Kit - IBL	96 Well



Distributed By:

**IBL-America, Inc.**

8201 Central Ave NE, Suite P

Minneapolis, MN 55432, USA

[info@ibl-america.com](mailto:info@ibl-america.com)

(888) 523 1246