


IGF-I ELISA

Enzyme Immunoassay for Quantitative Determination of
human Insulin-like Growth Factor I (IGF-I)
(IGFBP-blocked)
English

**For Research Use Only.
Not for use in diagnostic procedures.**



h **E20**



Distributed By:
IBL-America, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432, USA
info@ibl-america.com
(888) 523 1246













M: Aspenhastr. 25 • D-72770 Reutlingen / Germany
Phone: + 49 - (0) 7121 51484-0 • Fax: + 49 - (0) 7121 51484-10
E-Mail: contact@mediagnost.de • <http://www.mediagnost.de>

EN/ DE/ FR/ IT/ ES/ PT/ NL/ DK/ SE/ PL/ HU/ SK/ CZ/ BG/ EE/ GR/ RO/ SI/ FI

Symbols/ Symbole/ Symboles/ Simboli/ Símbolos/ Símbolos/ Symbolen/ Symboler/ Symboler/ Symbole/ Szimbólumok/ Symboly/ Symboly/ Символи/ Sümbolid/ Σύμβολα/ Simboluri/ Simboli/ Symbolit

DIN EN ISO 15223-1

	Expiry date/ Verfallsdatum/ Date de péremption/ Data di scadenza/ Fecha de caducidad/ Data de validade/ Uiterste gebruiksdatum/ Udløbsdato/ Bäst före-datum/ Termin ważności/ Lejárati idő/ Čas expirácie/ Doba expirace/ Срок на годност/ Aegumiskuupäev/ Ημερομηνία λήξης/ Data de expirare/ Rok uporabe/ Viimeinen käyttöpäivä
	Consider instructions for use/ Bitte Gebrauchsanweisung beachten/ Consultez la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consulte las instrucciones de uso/ Respeitar as instruções de utilização/ A.u.b de gebruiksaanwijzing volgen/ Se brugsanvisningen/ Läs anvisningarna före användning/ Proszę przeczytać instrukcję obsługi/ Vegye figyelembe a használati utasításban foglaltakat/ Postupujte podľa pokynov na použitie/ Dodržujte návod k použití/ Моля, спазвайте инструкцията за употреба/ Palun järgige kasutusjuhendit/ Λάβετε υπόψη σας τις οδηγίες χρήσης/ Vá rugám szá respectați instrucțiunile de utilizare/ Upoštečajte navodila za uporabo/ Lue käyttöohje huolellisesti!
	Lot-Batch Number/ Charge-Chargennummer/ Lot-Code du lot/ Lotto-Numero di lotto/ Lote-Código de lote/ Lote-Código do Lote/Lot-Partijnummer/ Lot-Batchkode/ Partisatskod/ Numer serii/ Tétel-sarzs szám/ Číslo šarže/ Číslo šarže/ Партиден номер/Partii number/ Παρτίδα-αριθμός παρτίδας/ Lot-număr lot/ Številka serije/ Erä
	Manufactured by/ Hergestellt von/ Fabriqué par/ Prodotto da/ Fabricado por/ Fabricado por/ Vervaardigd door/ Fabrikation af/ Tillverkad av/ Wyprodukowane przez/ Gyártotta/ Vyrobene/ Vyrobeno v/ Производител/ Τοοτja/ Κατασκευάζεται από/ Proodus de/ Proizvajalec/ Valmistaja
	Catalogue Number/ Bestellnummer/ Numéro de référence/ Numero di riferimento/ Número de referencia/ Número de Referência/ Referentienummer/ Referencennummer/ Bestellningsnummer/ Numer katalogowy/ Rendelési szám/ Katalógovné číslo/ Objednací číslo/ Каталоген номер/ Tellimisnumber/ Αρ. παραγγελίας/Număr comandă/ Številka naročila/ Viite tai tilausnumero
	Store at between/ Lagerung bei zwischen/ Conserver à entre/ Conservare a tra/ Conservar a temp. entre/ Armazena entre/ Bewaar bij tussen/ Opbevares mellem/ Förvaras vid/ Przechowywać w/ Tárolási tartomány/ Skladujte v rozsahu / Skladujte v rozmezí/ Температурно ограничение/ Säilätada temperatuuridel/ Φύλαξη σε θερμοκρασία/ Depozitare între/ Skladiščenje med/ Säilytys x-y Celsiusasteen lämpötilassa
	Contains sufficient for x tests/ Inhalt ausreichend für x Tests/ Contient suffisant pour x tests/ Contenuto sufficiente per x test/ Contenido suficiente para x pruebas/ Conteúdo suficiente para x testes/ Bevat voldoende voor x bepalingen/ Indeholder tilstrækkeligt til x prøver/ Innehållet räcker till x analyser/ Zawartość na x testów/ Tartalma x teszt elvégzésére elegendő/ Obsahuje materiál pre x testov/ Obsah dostatočný pro x testů/ Съдържание достатъчно за x тестове/ Sisust jätkub x katse jaoks/ Το περιεχόμενο επαρκεί για x δοκιμές/ Conținut suficient pentru x teste/ Vsebina zadostuje za x preizkusov/ Sisältö riittää x testille
	Incubation time/ Inkubationszeit/ Temps d'incubation/ Tempo d'incubazione/ Tiempo de incubación/ Tempo de incubação/ incubatietijd/ Inkubationstid/ inkubationstid/ Czas inkubacji/ Inkubációs idő/ Inkubačná lehota/ Inkubační doba/ Инкубационен период/ Inkubatsiooniaeg/ Χρόνος επώασης/ Timp de incubare/ Inkubacijska doba/ inkubaatioaika
	Incubate at/ Inkubation bei/ Incuber à/ Incubare a/ incubar a/ Incubar a/ incubatietemperatuur/ Inkubation ved/ inkubation vid/ Inkubacja przy/ Inkubáció hőmérséklete/ Inkubácia pri/ Inkubace při/ Инкубира се при/ Inkubatsioon temperatuuril/ Επώαση στους/ Incubare la/ Inkubacija pri/ Inkubaatiolämpötila
	Shaking/ Schütteln/ Mélanger/ Agitare/ Agitar/ Agitação/ Schudden/ Ryster/ Skaka/ Wstrząsanie/ Rázás/ Pretrepāt/ Profrepāt/ Разклащане/ Raputada/ Ανακινήστε/ Vibrare/ Stresite/ Sekoita
MTP	Mikrotiterplate/ Mikrotiterplatte/ plaque de microtitrage/ Piastra di microtitolazione/ Placa de microtitulación/ Placa de Microtitulação/ Microtiterplaat/ Mikrotiterplade/ mikrotiterplatta/ microtiterplaat/ Płytká microtiter/ Mikrotiter lap/ Mikrotitračná podložka/ Mikrotitrační podložka/ Mikrotitryná plaka/ Mikrotiterplaat/ Τρυβλίο μικροτιτλοδότησης/ Microplacă/ Mikrotitrská plošča/ Mikrotitruslevy
Rec in	Reconstitute in/ Rekonstituieren in/ Reconstituer dans/ Ricostituire nel/ Reconstituier en/ Reconstituier em/ Reconstituieren in/ Rekonstituér i/ Rekonstituera/ Rekonstytuować w/ Helyreállítás/ Znovu pripravit za/ Znovu pripraviti za/ Разтваряне в/ Moodustada uuesti/ Ανασυστήστε σε/ Reconstituire în/ Predelava v/ Rekonstituo
SPE	Sample/ Probe/ Echantillon/ Campione/ Muestra/ Amostra/ Monster/ Prøve/ prov/ Próbka/ Minta/ Vzorka/ Vzorek/ Проба/ Proov/ Δείγμα/ Probă/ Vzorec/ Näyte
DET	Antibody Conjugate/ Antikörperkonjugat/ Anticorps conjuguée/ Coniugato di anticorpo/ Conjugado de anticuerpos/ Conjugado anticorpo/ Antilichaamconjugaat/ Antistoffer-konjugat/ Antikropps-konjugat/ Koniugat antycial/ Antitest páros/ Protílátkový konjugát/ Protílátkový konjugát/ Антигяло конюгат/ Antikehad konjugaat/ Σύμπλοκο αντισώματος/ Compuși din anticorpi/ Antitelesa konjugat/ Vasta-aine konjugaatti
EC	Enzyme Conjugate/ Enzymkonjugat/ Conjugué enzymatique/ Coniugato di enzima/ Conjugado de enzimas/ Conjugado Enzima/ Enzymconjugaat/ Enzym-konjugat/ Enzymkonjugat/ Koniugat enzymów/ Enzim páros/ Enzymatický konjugát/ Enzymatický konjugát/ ензим конюгат/ Ensüümi konjugaat/ Σύμπλοκο –ενζύμου/ Compuși din enzime/ Encima konjugat/ Entsyymi konjugaatti
SB	Sample Buffer/ Probenpuffer/ Tampon d'échantillon/ Buffer campione/ Tampón de muestra/ Tampão de amostra/ Monsterbuffer/ Prøvebuffer/ Provbuffer/ Bufor próbki/ Mintapuffer/ Pufir na vzorku/ Vzorkovací pufir/ Примерен бyфep/ Proovipuhver/ Ρυθμιστικό διάλυμα δείγματος/ Tampon de probă/ Vzorční puffer/ Näytepuskuri
X:X	Dilute/ Verdünnen / Diluer / Diluire / Diluir / Diluir / Verdunnen / Fortyndes / Späd / Rozcieńczenie / Hígítás / Riedit / Ředit / Разреждане / Lahjendada / Αραιώστε / Diluați / Razredčiti / Laimennetaan

CAL	Calibrator X/ Kalibrator X/ calibrateur X/ calibratore X/ calibrador X/ calibrador X/ kalibrator X/ kalibrator X/ kalibrator X/ kalibrator X/ kalibrátor X/ kalibrátor X/ kalibrátor X/ kalibrátor X/ калибратор X/ kalibraator X/ Βαθμονομητής X/ calibrator X/ kalibrator X/ kalibraattori X
CTR1 / CTR2	Control X/ Kontrolle X/ Contôle X/ controllo X/ control X/ Contrôle X/ controle X/ Kontrol X/ Kontroll X/ kontrolne X/ Ellenőrző X/ Kontrolné X/ Kontrolní X/ Контролен X/ Kontroll X/ ελέγχου X/ control X/ Kontrolni X/ Kontrolli X
WB	Washing Buffer Concentrate/ Waschpufferkonzentrat/ Tampon de lavage conc./ Tampone di lavaggio concentrato/ Tampón de lavado concentrado/ Tampão de Lavagem Concentrado/ wasbuffer, geconcentreerd/ Vaskebufferkonzentrat/ Vaskebufferkonzentrat/ tvättbuffertkoncentrat/ Bufor płukania koncentrat/ Mosópuffer koncentrátum/ Koncentrát vymývacieho pufra/ Концентрат на промивен бυфер/ Pesurpuhvri kontsentraat/ Συμπύκνωμα ρυθμιστικού διαλύματος έκπλυσης/ Concentrat pentru tamponul de spălare/ Koncentrat izpiralne pufra/ Pesuliuositiiviste
WB 1:20	Washing Buffer/ Waschpuffer/ Tampon de lavage/ Tampone di lavaggio/ Tampón de lavado/Tampão de Lavagem/ wasbuffer/ Vaskebuffer/ tvättbuffert/ Bufor płukania/ Mosópuffer/ Vymývaci puffer/ Vymývaci pufri/ Промивен бυфер/ Pesurpuhver/ Ρυθμιστικό διάλυμα έκπλυσης/ Tampon pentru spălare /Izpiralni puffer/ Pesuliuos
S	Substrate/ Substrat/ Substrat/ Substrato/ Substrato/ Substrato/ substraat/ Substrat/ Substrat/ Substrat/ Substrat/ Szubsztrátum/ Substrát/ Substrát/ Субстрат Substraat/ Υπόστρωμα/ Substrat/ Substrat/ Substraattiliuos
STP	Stop Solution/ Stopplösung/ Stop Solution/ Soluzione di stop/ Stop Solución/ Solução Stop/ stopoplossing/ Stopopløsning/ Stopplösning/ Stop roztwór/ Megállító oldat/ Roztok na ukončenie/ Roztok pro ukončeni/ Стопираци разтвор/ Stopp-lahus/ Διάλυμα διακοπής/ Soluție de oprire/ Stop roztopina/ Pysäytysliuos
TAPE	Cover Plate with sealing tape/ Platte abkleben/ Recouvrir la microplaque avec bande adhésive/ Coprire la piastra con nastro adesivo/ Cubrir la placa con una cinta adhesiva/ Cobrir a Placa com fita adesiva/ plaatje met tape afdekken/ Afdækningsplade med tape/ maskera platta/ Odkleic płytke/ Tányér leragasztása/ Oblepit' podložku lepiacou páskou/ Olepit podložku lepící páskou/ Плака с лента за запечатване/ Katta plaat isoleerklleeplindiga/ Κολληστε το πλακίδιο με κολλητική ταινία/ Аоперити placa cu o bandă adezivă/ Prelepiti ploščo/ Peitä mikrotitrauslevy oheisella teipillä
MEASURE	Measure plate within 30 min at 450 nm (Referencefilter ≥ 590 nm) / Ausmessung innerhalb von 30 min bei 450 nm (Referenzfilter ≥ 590 nm) / Measure lábsorbance en léspace de 30 min à 450 nm avec ≥ 590 nm longueur d'onde pour référence/ Misurazione entro 30 min. a 450 nm (filtro di riferimento ≥ 590 nm) / Medición de la placa dentro de los siguientes 30 min a 450 nm (filtro de referencia ≥ 590 nm) / Medir a placa dentro de 30 min a 450 nm (Filtro de referência ≥ 590 nm) / Binnen 30 minuten bij 450 nm meten (referentiefilter ≥ 590 nm). / Mål plade i løbet af 30 min ved nm (referencefilter ≥ 590 nm) / Mät inom 30 min vid 450 nm (referensfilter ≥ 590 nm). / Pomiar w ciągu 30 min przy 450 nm (filtr odniesienia ≥ 590 nm) / Ki mérés 30 percen belül 450 nm-nél (referenciaszűrő ≥ 590 nm) / Merat' 30 minut pri 450 nm (Referenčných filtrov ≥ 590 nm) / Měřit 30 minut při 450 nm (Referenční filtr ≥ 590 nm) / Отчитане в рамките на 30 min при 450 nm (референтен филтър ≥ 590 nm) / Mõõtmise 30 min jooksul 450 nm korral (võrdlusfilter ≥ 590 nm) / Μέτρηση εντός 30 min στα 450 nm (φίλτρο αναφοράς ≥ 590 nm) / Măsurare în decurs de 30 min la 450 nm (filtru de referință ≥ 590 nm) / Izmerite ploščico v 30 min pri 450 nm (referenčni filter ≥ 590 nm) / Mittaa 30 minuutin aikana 450 nm:ssä (referenssi suodatin ≥ 590 nm)
Literature	Literature/ Literatur/ Bibliographie/ Letterario/ Bibliografía/ Literatura documentação/ literatuur/ Litteratur/ litteratur/ Literatura/ Irodalom/ Literatura/ Literatura/ Литература// Kirjandus/ Βιβλιογραφία/ Bibliografie/ literatura/ Lähdeluettelo
International Test description	International test description/ internationale Testanleitung/ description internationale de test/ Istruzioni per il test internazionali/ Descripción de ensayo internacional/ Descrição internacional do teste/ internationale testbeschrijving/ internationell testbeskrivning/ Opis testu międzynarodowego/ nemzetközi teszt-útmutató/ Medzinárodný návod k testu/ Mezinárodní návod k testu/ rahvusvaheline katse kirjeldus/ Διαθνεις οδηγίες για εργαστηριακές δοκιμές/ instructiuni internaționale pentru testare/ mednarodna navodila za preizkus/ Kansainvälinen käyttöohje
End	in all required wells/ in allen benötigten Vertiefungen/ dans tous les godets requis/ in tutti i pozzetti richiesti/ en todos los pozos requeridos/ em todos os tubos necessários/ in alle nodige putjes/ i alle nødvendige brønde/ i alla nödvändiga brunnar/ we wszystkich potrzebnych wgłębieniach/ minden szükséges forrásban/ vo všetkých potrebných miestach/ ve všech potřebných místech/ във всички необходими ямки/ kõigis vajalikes süvendites/ σε όλες τις απαραίτητες κοιλότητες/ în toate cavitățile necesare/ v vseh zahtevanih vdolbinah/ kaikkii tarvittaviin mikrotitrauslevyn syvennyksiin

TABLE OF CONTENTS

1	INTENDED USE	5
2	INTRODUCTION	5
3	PRINCIPLE.....	6
4	WARNINGS AND PRECAUTIONS	7
5	SAMPLES	8
6	MATERIALS.....	9
7	TECHNICAL NOTES.....	10
8	ASSAY PROCEDURE.....	11
9	EVALUATION OF RESULTS.....	12
10	LIMITATIONS OF PROCEDURE	13
11	EXEMPLARY VALUES	13
12	PERFORMANCE CHARACTERISTICS.....	15
13	ASSAY COMPARISON.....	17
14	REFERENCES	18
15	INTERNATIONAL ASSAY DESCRIPTION.....	19
16	ASSAY PROCEDURE.....	20

For Research Use Only.

Not for use in diagnostic procedures.

CAUTION: Not for human or animal therapeutic or diagnostic use.

For in vitro use only.

For professional use only.

Instructions for Use

IGF-I ELISA E20	96 Determinations
Principle of the test	Enzyme-linked Immunoassay
Duration (incubation period)	1.75 h
Antibodies	specific, monoclonal antibody and high-affinity polyclonal antiserum
Cross reactivity with IGF-II, Insulin, C-Peptide	< 0.1%
Buffers	Ready for use and 20fold concentrate
Calibrators	5 single Calibrators: 2 -50 ng/mL, recombinant human IGF-I
Reference material	International Standard WHO/NIBSC 02/254
Assay Range	0.09 – 1050 ng/mL
Controls	2 Controls, freeze-dried
Samples	human serum / plasma
Required sample volume	10 µL
Sample dilution	1:21
Analytical sensitivity	∅ 0.09 µg/L
Intra- / Interassay Variance	∅ < 10 %
Reference values	Blum W.F., Schweizer R. Insulin-Like Growth Factors and Their Binding Proteins. In: Ranke MB (ed): Diagnostics of Endocrine Function in Children and Adolescents. Basel, Karger, 2003, pp.166-199.

1 INTENDED USE

This enzyme immunoassay kit is for **research use only** and measures IGF-I in human serum or plasma or in various body fluids and in cell culture media.

2 INTRODUCTION

Insulin-like growth factors (IGF) I and II play a pivotal role in regulating the proliferation, differentiation and specific functions of many cell types (1-3). IGF-I is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 Dalton (5). Its major regulators are growth hormone (GH) and nutrition (6), although its production in specific tissues is affected by a multitude of tropic hormones and other peptide growth factors. IGFs are avidly bound to specific binding proteins (IGFBP). The seven classes of IGFBPs which are known at present (7,8,22). A major problem of IGF-I measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-I/IGFBP-3 complexes during the assay incubation.

IGF-I is present in low concentrations in various body fluids and in conditioned cell culture media of many cell lines. However, the determination of IGF-I in these specimens is particularly difficult due to the presence of IGFBPs usually in excessive amounts. The IGFBP-blocked IGF-I ELISA avoids these problems and allows the simple, correct and sensitive IGF-I determination in numerous samples with minimal expenditure of time.

Therefore, various techniques were applied to physically separate IGF-I from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2,12,13). These

techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries. The most widely used method is the acid-ethanol extraction (13,14) with a recovery of only 70-80 % of IGFBP-bound IGF-I as a result of co-precipitation. The absolute results of such an extraction are therefore false low (15). The extraction removes the IGF-BPs only insufficiently and leads to reduction in sensitivity of the assay due to pre-dilution of the samples by the extraction procedure. Furthermore, the remaining IGFBP may still interfere in the assay. In addition, the acid-ethanol extraction is ineffective in specimens other than serum or plasma (e.g. cell culture media), in which determination of IGF-I is already difficult enough due to the fact that IGF-BPs are frequently present at large excess. To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation is not required before measurement.

3 PRINCIPLE




Figure 2 Principle of the IGFBP-blocked IGF-I ELISA

In order to dissociate IGF-I from the IGF-BPs, the samples must be diluted in an acidic buffer (**Sample Buffer SB**) (Figure 2). The diluted samples are then pipetted into the assay wells. The IGF-I antiserum is dissolved in a buffer, which is able to neutralize the acidic samples. After the IGF-I antibody solution has neutralized the samples, the present excess IGF-II occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of the resulting free IGF-I. With this method, the IGF-BPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the extremely low cross-reactivity of the IGF-I antibody with IGF-II, the excess of IGF-II does not disturb the interaction of the first antibody with IGF-I.

The Mediagnost ELISA for IGF-I E20 is a so-called Sandwich-Assay using two specific and high affinity antibodies. The IGF-I in the samples binds to the first antibody coated on the microtiterplate, the second specific anti-IGF-I-antibody binds in turn to the immobilised IGF-I. The second antibody is biotinylated, the subsequently incubated Streptavidin-Peroxidase-Enzyme Conjugate will

bind to it, and thus in the final substrate incubation step colour development will be catalysed quantitatively depending on the IGF-I-level of the samples. The Calibrators of the ELISA E20 are prepared from recombinant IGF-I in concentrations of 2, 5, 15, 30 and 50 ng/mL.

4 WARNINGS AND PRECAUTIONS

For in-vitro use only. For Professional use only.

For Research Use Only. For Professional use only.

The Mediagnost kit is suitable only for in vitro and research use and not for internal use in humans or animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Caution: This kit contains material of human and/or animal origin. Therefore, all components and specimens should be treated as potentially infectious.

Do not use obvious damaged or microbial contaminated or spilled material.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Controls CTR1 and CTR2**

Source human serum for the Controls provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and specimens should be treated as potentially infectious.

Reagents CAL A-E, DET, EC, SB, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution STP

The Stop Solution contains 0.2 M acid sulphur acid (H₂SO₄)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

5.1 Sample type

Serum and Plasma

Serum and Heparin/EDTA Plasma yield comparable values. The IGF-I levels are reduced in citrate plasma samples, because of the relatively high amount of anticoagulant.

Furthermore, suitable samples are: cerebrospinal and urine samples as well as cell culture medium. IGF-I levels can vary considerably, the optimal dilution must be found out by the user.

5.2 Specimen collection

Use standard venipuncture for the blood sampling. Hemolytic reactions are to be avoided.

5.3 Required sample volume: 10 µL

5.4 Sample stability

In firmly closed sample vials:

- Storage at 20-25°C max. 24 hours
- Storage at -20°C : min. 2 years
- Freeze-thaw cycles max. 3

The storage of samples over a period of 2 years at -20°C showed no influence on the reading. Freezing and thawing of samples should be minimized. 3 freeze-thaw cycles showed no effect on samples.

5.5 Interference

Triglyceride, bilirubin and hemoglobin in the sample do not interfere to a concentration of 100 mg/mL and 200 µg/mL or 1 mg/mL. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.


5.6 Sample dilution

- Dilution: **1:21** with **Sample Buffer SB**
- Pipette **200 µL Sample Buffer SB** in PE-/PP-Tube (application of a multi-stepper is recommended in larger series); add **10 µL sample** (dilution 1:21). After mixing use 2 x 20 µL of this dilution in the assay.
- Attention: serum and plasma samples must be diluted at least 1:10 in **Sample Buffer SB** in order to achieve sufficient acidification of the samples.
- Depending on the expected IGF-I values the samples can be diluted higher in **Sample Buffer SB**.
- Sample stability after dilution of the sample: maximum 2 hours at 20-25°C.

6 MATERIALS

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the Calibration Curve.

MTP	Microtiter plate , ready for use, coated with mouse-anti-hIGF-I-antibody. Wells are separately breakable.	(8x12) wells
CAL A-E	Calibrators , lyophilized, (recombinant human hIGF-I), concentrations are given on vial labels and on quality certificate.	5 x 500 µL
CTR1	Control 1 , lyophilised, (human serum), concentration is given on quality certificate.	1 x 500 µL
CTR2	Control 2 , lyophilised, (human serum), concentration is given on quality certificate.	1 x 500 µL
DET	Antibody Conjugate , ready for use, contains goat biotinylated anti-hIGF-I antibody.	1 x 9 mL
EC	Enzyme Conjugate , ready for use, contains HRP (Horseradish-Peroxidase)-labelled Streptavidin.	1 x 12 mL
SB	Sample Buffer , ready for use.	1 x 25 mL
WB	Washing Buffer , 20-fold concentrated solution	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbencidine.	1 x 12 mL
STP	Stop Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x
	Instructions for use	1 x
-	Quality Certificate	1 x

6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ³ 590 nm

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C until its expiry date. The lyophilized reagents should be stored at –20°C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** Calibrators **A-E** and Controls **CTR1** and **CTR2** must be stored at –20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C.

Preparation of reagents

Bring all reagents to room temperature 20-25°C before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Do not use expired reagents. Temperature WILL affect the absorbance readings of the assay. However, correctness of the results will not be affected.

Use separate pipette tips for each sample, control and reagent to avoid cross contamination.

Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a Substrate that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.

Mix the contents of the microtiterplate wells thoroughly to ensure good test results. Do not reuse microtiterplate wells.

Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

Reconstitution

The Calibrators **A–E** and Control **CTR1** and **CTR2** are reconstituted with the Sample Buffer **SB**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Control **CTR1** and **CTR2** with the Sample Buffer **SB** in the same ratio (1:21) as the sample.

The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20-fold concentrate with Aqua dest..

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive—store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a microtiter plate shaker. We recommend 350 rpm. Depending on the design of the shaker, the shaking frequency should be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8 ASSAY PROCEDURE

NOTES: When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody Conjugate **DET**, Enzyme Conjugate **EC** and Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**. All determinations (Blank, Calibrators **A-E**, Control **CTR1** and **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

- 1) Add **80 µL Antibody Conjugate DET** in all wells used.
- 2) Pipette in positions A1/2 **20 µL Sample Buffer SB**.
- 3) Pipette in positions B1/2 **20 µL of the Calibrator A (2 ng/ml)**
Pipette in positions C1/2 **20 µL of the Calibrator B (5 ng/ml)**,
Pipette in positions D1/2 **20 µL of the Calibrator C (15 ng/ml)**,
Pipette in positions E1/2 **20 µL of the Calibrator D (30 ng/ml)**,
Pipette in positions F1/2 **20 µL of the Calibrator E (50 ng/ml)**.
To control the correct accomplishment of the assay **20 µL** of the 1:21 (or in respective dilution ratio of the samples) in **Sample Buffer SB** diluted **Controls CTR1&CTR2** can be pipetted in positions G1/2 and H1/2.
Pipette **20 µL** each of the diluted sample (e.g. dilute 1:21 with Sample Buffer **SB**) in the rest of wells, according to your requirements.
- 4) Cover the wells with sealing tape and incubate the plate for **1 hour at room temperature 20-25°C** (shake at 350 rpm).
- 5) After incubation aspirate the contents of the wells and wash the wells 5 times **300 µL Washing Buffer WB** / well.
- 6) Following the last washing step pipette **100 µL** of the **Enzyme Conjugate EC** in each well.
- 7) Cover the wells with sealing tape and incubate the plate for **30 minutes at room temperature 20-25°C** (shake at 350 rpm).
- 8) After incubation wash the wells 5 times with **Washing Buffer WB** as described in step 5.
- 9) Pipette **100 µL** of the **Substrate S**.
- 10) Incubate the plate for **15 minutes in the dark at room temperature 20-25°C**
- 11) Stop the reaction by adding **100 µL Stop Solution STP** to all wells.
- 12) Measure the absorbance within **30 minutes at 450 nm** (**Reference filter ³ 590 nm**).

9 EVALUATION OF RESULTS

9.1 Establishing of the Calibration Curve

The International Standard for hIGF-I, WHO NIBSC Code 02/254 was used as calibrator material and the following IGF-I concentrations are used. For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of Calibrator **E** should be above 1.00.

Calibrator	A	B	C	D	E
ng/mL	2	5	15	30	50
nmol/L	0.26	0.66	1.96	3.92	6.54

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples, Controls and Calibrators.
- 3) Plot the Calibrator concentrations on the x-axis versus the mean value of the absorbance of the Calibrators on the y-axis.
- 4) Recommendation: Calculation of the Calibration Curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The IGF-I concentration in ng/mL of the samples and controls can be calculated by **multiplication** with the respective **dilution factor**.

9.2 Example of a typical Calibration Curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

	Blank	A	B	C	D	E
ng/mL	0.0	2	5	15	30	50
OD(450-620 nm)	0.00	0.088	0.299	0.985	1.727	2.543

The exemplary shown Calibration Curve in **Figure 3** cannot be used for calculation of your test results. You have to establish a Calibration Curve for each test you conduct!




Figure 3 Exemplary Calibration Curve

9.3 Exemplary calculation of IGF-I concentrations

Sample dilution: 1:21

Measured extinction of the blank	0.0165
Measured extinction of your sample	0.2695

Your measurement program will calculate the IGF-I concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 3rd degree).

In this exemplary case the following equation is solved by the program to calculate the IGF-I concentration in the sample:

$$0.253 = -0.0004x^2 + 0.0727x - 0.0471$$
$$4.57 = x$$

If the dilution factor (1:21) is taken into account the **IGF-I** concentration of the undiluted sample is

$$4.57 \text{ ng/mL} \cdot 21 = 96 \text{ ng/mL}$$

9.4 Interpretation of results

The test is for research use only. Not for use in diagnostic procedures. Therefore, interpretation of test results depends on experimental settings and scientific question.

10 LIMITATIONS OF PROCEDURE

The Mediagnost IGF-I ELISA, E20 is based on a combination of monoclonal capture and polyclonal detection antibodies. Generally, immunological assays are sensible to heterophilic antibodies and rheumatoid factors in the sample. Their influence is reduced by the assay design, but cannot be excluded completely.

11 EXEMPLARY VALUES

IGF-I levels are highly age-dependent in children, less so in adults until the age of about 60. The exemplary ranges in various age groups, which are log-normally distributed, are given in Table 2 by percentiles.

Between 8 and 19 years of age, values are given for boys and girls separately, because the pubertal peak usually occurs approximately 2 years earlier in girls.

Table 1 Exemplary range of serum IGF-I levels given in ng/mL at different pubertal stages according to Tanner. Because no significant difference between boys and girls is observed, both sexes are combined. Only children and adolescents between 7 and 17 years of age are included.

Pubertal Stage	Percentile			
	0.1th	5th	50th	95th
1	61	105	186	330
2	85	156	298	568
3	113	196	352	631
4	171	268	431	693
5	165	263	431	706

Table 2 Serum levels of IGF-I in healthy subjects at various ages. Individuals between 8 and 19 years of age were classified according to gender, as the pubertal peak occurs almost 2 years earlier in girls than in boys.

Age	Percentile													
	0.1	1	5	10	20	30	40	50	60	70	80	90	95	99
0-2 y.	13	20	28	34	43	50	58	66	75	87	102	128	156	220
2-4 y.	20	29	40	48	59	68	77	87	98	111	129	159	189	260
4-6 y.	26	36	50	59	73	85	96	108	122	138	160	196	233	320
6-7 y.	34	46	62	72	87	99	111	124	138	155	176	212	248	332
7-8 y.	45	60	78	90	107	121	134	148	163	181	205	243	281	364
8-9 y. boys	54	71	90	102	119	133	146	160	175	192	214	250	284	362
girls	55	75	99	115	137	156	174	193	214	239	271	324	376	496
9-10 y. boys	63	82	102	115	133	148	162	176	191	209	232	269	304	379
girls	68	89	114	130	152	170	187	205	224	247	276	323	369	469
10-11 y. boys	77	96	117	130	148	162	176	189	203	220	241	274	305	370
girls	81	106	134	153	178	199	219	239	261	287	321	374	426	539
11-12 y. boys	85	106	129	144	163	179	194	209	225	244	267	304	339	413
girls	91	123	160	185	220	248	276	305	337	374	424	503	581	758
12-13 y. boys	88	112	141	159	184	204	223	243	264	289	321	371	419	525
girls	116	155	201	231	274	309	342	377	415	460	519	614	707	914
13-14 y. boys	111	143	179	203	235	261	286	311	339	371	412	477	540	677
girls	163	207	256	287	329	364	395	428	463	504	556	637	716	884
14-15 y. boys	140	182	229	260	303	337	370	404	441	484	539	625	691	896
girls	193	236	284	314	353	385	414	443	474	510	556	628	713	832
15-16 y. boys	176	221	269	299	340	372	402	433	466	504	552	626	697	849
girls	187	231	279	309	350	382	412	442	474	512	559	632	700	845
16-17 y. boys	178	221	267	296	335	366	395	424	455	491	537	607	673	814
girls	183	225	270	298	336	366	394	422	452	486	530	597	660	792
17-18 y. boys	173	207	243	265	294	317	337	358	380	405	436	484	527	618
girls	176	210	246	268	297	320	341	362	384	409	441	488	533	624
18-19 y. boys	167	201	235	256	285	307	327	347	368	393	423	469	512	600
girls	167	199	233	254	281	302	322	341	362	385	414	458	499	583
19-20 y.	158	189	220	240	265	285	304	322	341	363	391	433	471	550
20-30 y.	72	92	115	130	150	167	182	198	215	235	261	302	340	425
30-40 y.	68	87	109	123	142	158	173	188	204	223	248	287	324	404
40-50 y.	64	82	103	116	135	150	164	178	194	212	235	272	310	385
50-60 y.	60	77	97	110	127	142	155	169	184	201	224	260	292	369
60-70 y.	55	72	91	103	120	134	147	161	176	193	215	251	282	362
70-80 y.	25	35	47	55	67	78	88	98	110	124	142	173	207	276
>80 y.	21	30	40	47	58	67	76	85	95	108	125	153	184	245

Serum concentrations are given in ng/ml.

Exemplary values have been evaluated by Prof Blum by a radioimmunoassay identical to Mediagnost IGF-R20. Thus, the age and sex specific exemplary values published in Diagnostics of Endocrine Function in Children and Adolescents. (Edited by Prof Ranke. ISBN-3-335-00496-5) can be applied to all Mediagnost IGF-I assays.

12 PERFORMANCE CHARACTERISTICS

12.1 Sensitivity

Sensitivity was assessed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. The analytical sensitivity of the Mediagnost E20 is 0.091 ng/mL on average, in 19 independent determinations values from 0.03 ng/mL to 0.2 ng/mL were found.

12.2 Specificity

The measurements of E20 cross reactivity with IGF-II, Insulin and C-Peptide. These IGF-related proteins were added to assay buffer in the indicated concentration and the solution was applied as sample without any further dilution. The concentration measured within the blank without any protein was 0.78 µg/L. Thus, neither IGF-II nor Insulin or C-Peptide are measured by the Mediagnost E20 ELISA (see table 3).

Table 3 Specificity. Cross reactivity of the test system with different IGF-I related proteins.

added C-Peptide [µg/L]	measured IGF-I [µg/L]	added Insulin [µg/L]	measured IGF-I [µg/L]	added IGF-II [µg/L]	measured IGF-I [µg/L]
500	0.73	100	0.78	1250	0.77
100	0.78	10	0.77	750	0.73
10	0.77	1	0.76	250	0.77
0	0.78	0	0.78	0	0.78

12.3 Cross Reactions with Animal Samples

Several commercially available animal sera have been used as samples in this assay and therewith it is proven, that the test can be used as heterologous assay for IGF-I measurement in serum samples of primates, cattle, pig, sheep, horse, donkey, goat, dog, cat, rabbit and guinea pig. For the determination of IGF-I in chicken, rat and mouse sera is this kit not usable. Species specific calibration has to be done by the user.

To determine exact species-specific IGF-I concentrations, the test system must be calibrated respectively: for example, using bovine IGF-I in a known concentration.

12.4 Precision

Intra-Assay Variance

Three samples have been measured severalfold in the same assay. The results are shown in table 4. The measured coefficient of variation (CV) is 5.81% on average.

Table 4 Intra-Assay variability

	Number of determinations	Mean value (ng/mL)	Standard deviation (ng/mL)	VC (%)
Sample 1	18	144.8	9.63	6.65
Sample 2	18	140.79	7.15	5.08
Sample 3	18	138.02	7.86	5.69

Inter-Assay and Lot-to-Lot Variance

Serum samples were measured in independent assays. Exemplary results are shown in Table 5. Further, five samples were also tested repeatedly throughout four years in eight different lots. The variability on average is 8.57% (6.8 – 10.5%).

Table 5 Inter-Assay variability

	Number of determinations	Mean value [ng/mL]	Standard deviation [ng/mL]	VC [%]
Sample 1	8	81	5.34	6.56
Sample 2	16	192	12.38	6.43
Sample 3	17	498	27.52	5.53

12.5 Linearity

Linearity was tested by dilution of native sera with different IGF-I contents (Sample 1-3). The amount of measured IGF-I was recalculated and is shown in Figure 4.




Figure 4 Linearity, recalculated IGF-I concentrations of differentially diluted samples. The minimal dilution is 1:10, the recommended dilution is 1:21.

12.6 Recovery and Accuracy

Recombinant IGF-I was added in different amounts to human serum. The IGF-I content of the so enriched samples was measured and recovery in comparison to enriched buffer calculated. Results are shown in Table 6.

Table 6 Recovery of recombinant IGF-I in human serum

IGF-I [µg/L]	Sample 1	Sample 2	Sample 3	Sample 4
Sample	138	172	133	180
+ IGF-I 200	287	372	-	-
+ IGF-I 400	-	-	539	591
% Recovery	85	100	101	102

12.7 Interference

Interference of bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing IGF-I. For comparison the same amount of buffer without any substance was also added to the serum. Table 7 demonstrates that neither bilirubin nor triglycerides exert any influence on the measurement of IGF-I in human serum.

Table 7 Interference of physiologic substances on IGF-I measurement. Human serum samples were enriched with different amounts of triglycerides, bilirubin or hemoglobin and the recovery of IGF-I was measured. Here the relative recovery in [%] of not enriched samples is shown.

	Triglyceride 100 mg/mL	Bilirubin 200 µg/mL	Hemoglobin 10 mg/ml
Serum 1	93	90	97
Serum 2	100	101	110
Serum 3	120	120	104

Influence of binding proteins on IGF-I measurement was exemplarily elucidated with IGFBP-3.

Different amounts of IGF-I and 3 or 6 mg/L of IGFBP-3 were added to sample buffer (pH 2) and phosphate based saline buffer (pH 7.4). After a short incubation of 15 minutes at room temperature these samples were diluted and applied to the Mediagnost E20 as described in

the package insert. In case of sample buffer IGFBP-3 up to 6 mg/L did not interfere with IGF-I measurement. But without acidification of the sample a strong interference of IGFBP-3 with IGF-I measurement was detected (Table 8).

Table 8 Interference of IGFBP-3 with IGF-I measurement

IGFBP-3	Sample Buffer		
	50 µg/L IGF-I	100 µg/L IGF-I	300 µg/L IGF-I
-	46.38	116.14	358.1
3 mg/L	47.33	115.83	384.15
6 mg/L	52.32	123.38	355.41
IGFBP-3	Phosphate buffered Saline		
	50 µg/L IGF-I	100 µg/L IGF-I	300 µg/L IGF-I
-	34.2	90.23	349.04
3 mg/L	7.4	12.16	152.14
6 mg/L	7.2	10.12	48.15

12.8 Assay Calibration

Recombinant human IGF-I produced by E. coli and of >98% purity (SDS-PAGE, Silverstain) is used as Calibrator within the assay. This recombinant hIGF-I devoid of methIGF-I or IGF-I variants with mismatched disulfide bonds is identical to the major authentic IGF-I form in blood. The traceability of this recombinant Calibrator material to the international reference material of the WHO 02/254 has been proven. Results are published by Burns C et al. in Growth Horm IGF Res. 2009 Oct;19(5):457-62. Epub 2009 Mar 20. Mediagnost E20 ELISA is coded by 14c.

The reference material includes **8.5 µg/ampoule** IGF-I measured by amino acid analysis and HPLC. Mediagnost E20 immunoassay (assay No. 14c) measures **11.55 µg/ampoule**. The mean of all tested immunoassays is 11.61 µg/ampoule.

Thus, Mediagnost results are comparable to other immunological tests for measurement of IGF-I and can easily be transformed to WHO 02/254 (25, 26) by a factor of **0.735**.

13 ASSAY COMPARISON

Mediagnost E20 IGF-I was compared with the Mediagnost R20 IGF-I. 196 serum samples were measured in both assays and an excellent coefficient of correlation was shown with $r = 0.95$. Additionally, the Mediagnost IGF-I ELISA E20 was compared with an Enzyme-Immunoassay of other commercially available IGF-I test and a correlation of $R^2 > 0.9$ was shown.




Figure 5 Assay Comparison Mediagnost RIA R20 and Mediagnost ELISA E20

14 REFERENCES

Literature

- 1) Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem.*25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors.* New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology* 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding pro-teins. *J. Clin Endocrinol Metab.* 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
- 15) Ranke MB (ed.): *Diagnostics of Endocrine Function in Children and Adolescents,* Basel, Karger, 2003, pp. 166-199
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med.* 301:1138-1142
- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest.* 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence.* Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab.* 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh).*118:374-380.
- 22) Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. *J Clin Endocrinol Metab* Vol 82, 4:1301-1303
- 23) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Horm Res* 54:60-68
- 24) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Horm Res* 55:155-124
- 25) Address NIBSC: Blance Lane, South Mimms, Potters Bar, Hertford EN 6 3 QG, Great Britain
- 26) Burns C, Rigsby P, Moore M, Rafferty B. (2009) The first International Standard for Insulin-like Growth Factor-1 (IGF-1) for immunoassay: Preparation and calibration in an international collaborative study. *Growth Hormone & IGF Research* (2009)
- 27) Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGFBPs. *Growth Regulation* 4 (Suppl. 1):11-19
- 28) Ranke MB, Osterziel KJ, Schweitzer R, Schuett B., Weber K., Röbbel P, Vornwald A, Blumenstock G, Elmlinger MB (2003) Reference Levels of Insulin-Like-Growth Factor I in the Serum of Healthy Adults: Comparison of Four Immunoassays. *Clin Chem Lab Med* Vol 41(10):1329-1334
- 29) European Medicines Agency, (2009) Guideline on bioanalytical method validation. EMEA/CHMP/192217/2009 Rev.1 Corr.
- 30) U.S. Department of Health and Human Services, Food and Drug Administration (2001) Guidance for Industry Bioanalytical Method Validation.

15 INTERNATIONAL ASSAY DESCRIPTION

International Test Description

CAL A-E	Rec in 500 µL SB	-
CTR1	Rec in 500 µL SB	1:21 SB
CTR2	Rec in 500 µL SB	1:21 SB
WB 20x	-	1:20 A. dest. è WB 1:20
SPE		1:21 SB
-	°C 20-25°C, A ↔ max. 2 h	
80 µL	DET	A1 - End
20 µL	SB	A1/A2
20 µL	CAL A (2 ng/ml)	B1/B2
20 µL	CAL B (5 ng/ml)	C1/C2
20 µL	CAL C (15 ng/ml)	D1/D2
20 µL	CAL D (30 ng/ml)	E1/E2
20 µL	CAL E (50 ng/ml)	F1/F2
20 µL	CTR1 1:21 SB	G1/G2
20 µL	CTR2 1:21 SB	H1/H2
20 µL	SPE 1:21 SB	
TAPE		
A 1 h °C 20-25°C ↔ 350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	EC	
TAPE		
A 0.5 h °C 20-25°C ↔ 350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	S	
A 15 min °C 20-25°C 		
STP		
MEASURE		

16 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 500 µL Sample Buffer SB	-
CTR1	Control 1	in 500 µL Sample Buffer SB	1:21 with Sample Buffer SB
CTR2	Control 2	in 500 µL Sample Buffer SB	1:21 with Sample Buffer SB
WB	Washing Buffer conc.		1:20 with Agua dest. è WB 1:20
Samples + Controls CTR1 and CTR2: dilute 1:21 in Sample Buffer SB , mix immediately, incubate max. 2h. Use 20 µL for each well in the assay.			
Before assay procedure bring all reagents to room temperature (20-25°C).			
Assay procedure in double determination			
Pipette	Reagents		Position
80 µL	Antibody Conjugate DET		in all wells used
20 µL	Sample Buffer SB (Blank)		A1/A2
20 µL	Calibrator A (2 ng/mL)		B1/B2
20 µL	Calibrator B (5 ng/mL)		C1/C2
20 µL	Calibrator C (15 ng/mL)		D1/D2
20 µL	Calibrator D (30 ng/mL)		E1/E2
20 µL	Calibrator E (50 ng/mL)		F1/F2
20 µL	Control CTR1 (1:21 diluted)		G1/G2
20 µL	Control CTR2 (1:21 diluted)		H1/H2
20 µL	Sample SPE (1:21 diluted)		in the rest of the wells according the requirements
Cover the wells with the sealing tape.			
Sample Incubation: 1 h at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 /well		In each well
100 µL	Enzyme Conjugate EC		In each well
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at 20-25°C, 350 rpm			
5x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well		In each well
100 µL	Substrate S		In each well
Incubation: 15 Minutes in the Dark at 20-25°C			
100 µL	Stop Solution STP		In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			