

ASSAY PROCEDURE

Specimen Pretreatment: **None.**

All reagents must reach room temperature before use. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.

1. Prepare working solutions of the free testosterone-HRP conjugate and wash buffer.
2. Remove the required number of well strips. Reseal the bag and return any unused strips to the refrigerator.
3. Pipette 25 µL of each calibrator, control and specimen sample into correspondingly labelled wells in duplicate.
4. Pipette 100 µL of the conjugate working solution into each well (We recommend using a multichannel pipette).
5. Gently shake the plate for 10 seconds.
6. Incubate the plate at 37°C for 1 hour.
7. Wash the wells 3 times with 350 µL of diluted wash buffer per well and tap the plate firmly against absorbent paper to ensure that it is dry. (The use of a washer is recommended.)
8. Pipette 150 µL of TMB substrate into each well at timed intervals.
9. Incubate the plate at 37°C for 10–15 minutes (or until calibrator A attains dark blue colour for desired OD).
10. Pipette 50 µL of stopping solution into each well at the same timed intervals as in step 8.
11. Read the plate on a microplate reader at 450 nm within 20 minutes after addition of the stopping solution.

* If the OD exceeds the upper limit of detection or if a 450 nm filter is unavailable, a 405 or 415 nm filter may be substituted. The optical densities will be lower, however, this will not affect the results of patient/control samples.

CALCULATIONS


1. Calculate the mean optical density of each calibrator duplicate.
2. Draw a calibrator curve on semi-log paper with the mean optical densities on the Y-axis and the calibrator concentrations on the X-axis. If immunoassay software is being used, a 4-parameter or 5-parameter curve is recommended.
3. Calculate the mean optical density of each unknown duplicate.
4. Read the values of the unknowns directly off the calibrator curve.

TYPICAL TABULATED DATA

Sample data only. **Do not** use to calculate results.

Calibrator	Mean OD (450 nm)	Value (pg/mL)
A	2.292	0
B	1.680	0.1
C	1.181	1
D	0.780	5
E	0.426	20
F	0.214	60
Unknown	1.066	1.59
Unknown	0.441	19.6

TYPICAL CALIBRATOR CURVE



PERFORMANCE CHARACTERISTICS

SENSITIVITY

The limit of detection (LoD) was determined from the analysis of 64 replicates of a low value sample and from the LoB. $LoD = LoB + 1.645\sigma_S$, where σ_S is the standard deviation of the low value sample. σ_S was determined to be 0.0093 based on 64 measurements of a low value sample. $LoD = 0.0025 + (1.645 \times 0.0093) = 0.018$ pg/mL.

SPECIFICITY (CROSS-REACTIVITY)

The following compounds were tested for cross-reactivity with the Direct Free Testosterone ELISA kit with testosterone cross-reacting at 100%.

Steroid	% Cross Reactivity
Testosterone	100
5 α -DHT	3.5
Androstenedione	0.17
Progesterone	0.007
Androsterone	0.075
Aldosterone	< 0.008
Cholesterol	< 0.0001
Cortisone	0.0025
DHEA	0.071
DHEAS	0.0014
17 β -Estradiol	0.15
Estriol	< 0.008
Pregnenolone	0.028

INTRA-ASSAY PRECISION

Five samples were assayed 24 times each on the same calibrator curve. The results (in pg/mL) are tabulated below:

Sample	Mean	CV %
1	2.24	6.7
2	3.81	6.4
3	13.6	6.0
4	13.7	5.9
5	23.7	4.8

INTER-ASSAY PRECISION

Three samples were assayed twenty times in duplicate over a period of greater than ten days. The results (in pg/mL) are tabulated below:

Sample	Mean	CV %
1	3.53	8.1
2	13.8	11.5
3	23.3	6.9

COMPARATIVE STUDIES

The Direct Free Testosterone ELISA Kit (y) was compared with a competitor's Free Testosterone Coated Tube RIA Kit (x). The comparison of 60 serum samples yielded the following linear regression results:

$$y \text{ (IBL)} = 0.9362x \text{ (competitor)} + 3.8794, r = 0.97$$

EFFECT OF SEX HORMONE BINDING GLOBULIN (SHBG)

The purpose of this study was to investigate a possible interference caused by the binding of SHBG to the free testosterone-HRP conjugate. A charcoal-stripped human serum pool was spiked precisely with SHBG at concentrations ranging from 6.25–200 µg/mL and was assayed with the Free Testosterone ELISA Kit. Results tabulated below (in pg/mL):

SHBG Added	OD 450 nm	Percent B/B ₀ (%)
0	2.37	100.0
6.25	2.37	99.9
12.5	2.34	98.7
50	2.36	99.5
200	2.27	95.6

The results showed % binding values between 95–100% (B₀ = unspiked serum) even at higher than normal SHBG levels. In conclusion, the results showed that there was no significant binding by SHBG on the free testosterone-HRP conjugate.

EXPECTED NORMAL VALUES

As for all clinical assays each laboratory should collect data and establish their own range of expected normal values. The results of an expected range study with apparently normal healthy subjects yielded the following results (all values are reported in pg/mL):

Cohort Group; Gender/Age	N	95% Confidence Range	Absolute Range
Males / < 13	44	–	ND–1.6
Males / 13–19	37	–	ND–22.3
Males / 20–39	120	9.1–32.2	–
Males / 40–59	120	5.7–30.7	–
Males / ≥ 60	120	5.9–27.0	–
Females / < 13	63	–	ND–1.3
Females / 13–19	17	–	0.2–2.0
Females / 20–39	120	0.1–6.3	–
Females 40–59	120	0.2–4.1	–
Females / ≥ 60	60	0.5–3.9	–

REFERENCES

1. Winter SJ, et al. The Analog Free Testosterone Assay: Are the Results in Men Clinically Useful. *Clin Chem.* 1998; 44(10):2178–82.
2. Ooi DS, et al. Establishing Reference Interval for DPC's Free Testosterone Radioimmunoassay. *Clin Biochem.* 1998; 31(1):15–21.
3. Marcus GJ, Dumford R. A Simple Enzyme-Linked Immunosorbent Assay for Testosterone. *Steroids.* 1985; 46:975–86.

4. Joshi UM, et al. A Sensitive Specific Enzyme Immunoassay for Serum Testosterone. *Steroids.* 1979; 34(1):35–46.
5. Swinkels LM, et al. Salivary and Plasma Free Testosterone and Androstenedione Levels in Women. *Am Clin Biochem.* 1988; 25(Pt 4):354–9.
6. Swinkels LM, et al. A Symetric Dialysis Method for the Determination of Free Testosterone in Human Plasma. *Clin Chem Acta.* 1987; 165(2–3):341–9.
7. Ekins R. Hirsutism: Free and Bound Testosterone. *Ann Clin Biochem.* 1990; 27(Pt 1):91–4.
8. Manni A, et al. Bioavailability of Albumin-Bound Testosterone. *J Clin Endo Metab.* 1985; 61(4):705–10.
9. Ekins R. The Science of Free Testosterone Measurement. *Proc UK NEQAS Meeting.* 1998; 3:35–9.
10. Longcope C, et al. Free Estradiol, Free Testosterone, and Sex Hormone-Binding Globulin in Perimenopausal Women. *J Clin Endo Metab.* 1987; 64(3):513–8.
11. Vermeulen A, et al. The Apparent Free Testosterone Concentration: An Index of Androgenicity. *J Clin Endo Metab.* 1971; 33(5):759–67.
12. Paulson JD, et al. Free Testosterone Concentration in Serum. Elevation is the Hallmark of Hirsutism. *Am J Obs Gyneco.* 1977; 128(8):851–7.
13. Cumming DC, Wall SR. Non-Sex Hormone-Binding Globulin-Bound Testosterone as a Marker for Hyperandrogenism. *J Clin Endo Metab.* 1985; 61(5):873–6.
14. Baxendale PM, et al. Salivary Testosterone: Relationship to Unbound Plasma Testosterone in Normal and Hyperandrogenic Women. *Clin Endocrinol (Oxf).* 1982; 16(6):595–603.
15. Biffignandi P, et al. Female Hirsutism: Pathophysiological Considerations and Therapeutics Implications. *Endocrinol Rev.* 1984; 5(4):498–513.
16. Wu CH. Plasma Free and Protein-Bound Testosterone in Hirsutism. *Obstet Gynecol.* 1982; 60(2):188–94.
17. Bamman BL, et al. Total and Free Testosterone During Pregnancy. *Am J Obstet Gynecol.* 1980; 137(3):293–8.
18. Halpern EP, et al. Production of Anti-Testosterone Antisera in Rabbits. *Clin Chem.* 1980; 26:68.
19. Wheeler MJ. The Determination of Bio-Available Testosterone. *Ann Clin Biochem.* 1995; 32(Pt 4):345–57.
20. Check JH, et al. Falsely Elevated Steroidal Assay Levels Related to Heterophile Antibodies Against Various Animal Species. *Gynecol Obstet Invest.* 1995; 40(2):139–40.

SYMBOLS

European Conformity	In vitro diagnostic device	Consult instructions for use
Contains sufficient for <n> tests	Storage Temperature	Legal Manufacturer
Use by	Catalogue Number	Authorized representative
Lot number	Dilute 1: # Before use	

CONTACT INFORMATION

Manufactured for:

Immuno-Biological Laboratories, Inc.
8201 Central Ave. NE, Suite P
Minneapolis, MN 55432, USA
Phone: +1 (763) - 780-2955
Email: info@ibl-america.com
Web: www.ibl-america.com