

Manufactured for Immuno-Biological Laboratories Inc.
(IBL-America)
8201 Central Avenue, NE, Suite P
Minneapolis, MN 55432
Tel: 763-780-2955
Toll Free: 1-888-523-1246



Instructions for use

Histamine(Cell-cult/Urine)ELISA

Enzyme Immunoassay for the determination of Histamine
in urine and cell culture samples and for the
determination of Total Histamine in whole blood.

For research use only, not for use in diagnostic procedures.

REF

IB89143



Histamine ELISA ^{Fast Track} (Cell culture, Urine and Whole blood)

1. Introduction

1.1 Intended use and principle of the test

Enzyme Immunoassay for the determination of Histamine in urine and cell culture samples and for the determination of Total Histamine in whole blood. For research use only, not for use in diagnostic procedures.

First, Histamine is acylated. The subsequent competitive ELISA kit uses the microtiter plate format. The antigen is bound to the solid phase of the microtiter plate. The acylated analyte concentrations in the standards, controls and samples and the solid phase bound analyte compete for a fixed number of antibody binding sites. After the system is in equilibrium, free antigen and free antigen-antibody complexes are removed by washing. The antibody bound to the solid phase is detected by an anti-goat IgG-peroxidase conjugate using TMB as a substrate. The reaction is monitored at 450 nm.

Determination of unknown samples is achieved by comparing their absorbance with a standard curve prepared with known standard concentrations.

1.2 Background

Histamine belongs to the biogenic amines and is synthesized by decarboxylation from the amino acid histidine. It is synthesized by mast cells, basophils, platelets, histaminergic neurons, and enterochromaffine cells, where it is stored intracellularly in vesicles and released on stimulation.

Histamine acts by binding to its 4 receptors (H1R, H2R, H3R and H4R) on target cells in various tissues. It causes smooth muscle cell contraction, vasodilatation, increased vascular permeability and mucus secretion, tachycardia, alterations of blood pressure, and arrhythmias.

In humans, histamine is one of the most important mediators and takes part in the initial phase of an anaphylactic reaction ("immediate type" allergy).

Of interest is also the determination of the histamine *release* from basophilic leucocytes in allergies.

2. Procedural cautions, guidelines, warnings and limitations

2.1 Procedural cautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for certain types of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- (4) The principles of Good Laboratory Practice (GLP) have to be followed.
- (5) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (6) All kit reagents and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (7) For dilution or reconstitution purposes, use deionized, distilled, or ultra-pure water.
- (8) The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided.
- (9) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (10) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (11) Incubation times do influence the results. All wells should be handled in the same order and time intervals.
- (12) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (13) A standard curve must be established for each run.
- (14) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report provided with the kit.
- (15) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (16) Avoid contact with Stop Solution containing 0.25 M H₂SO₄. It may cause skin irritation and burns. In case of contact with eyes or skin, rinse off immediately with water.

- (17) Some reagents contain sodium azide (NaN₃) as preservatives. In case of contact with eyes or skin, rinse off immediately with water. NaN₃ may react with lead and copper plumbing to form explosive metal azides. When disposing reagents, flush with a large volume of water to avoid azide build-up.
- (18) TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them.
- (19) For information on hazardous substances included in the kit please refer to Material Safety Data Sheet (MSDS). The Material Safety Data Sheet for this product is made available directly on the website of the manufacturer or upon request.
- (20) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (21) The results obtained with this test kit are for research use only, not for use in diagnostic procedures.
- (22) Kit reagents must be regarded as hazardous waste and disposed of according to national regulations.

2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

2.2.1 Interfering substances

24-hour urine

Please note the sample preparation! If the percentage of the final concentration of acid is too high, the buffer capacity of the Acylation Buffer is insufficient. As a consequence histamine will not be acylated quantitatively.

2.2.2 Drug interferences

There are no known substances (drugs) which ingestion interferes with the measurement of histamine level in the sample.

2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

4. Materials

4.1 Contents of the kit

BA D-0024	REAC-PLATE	Reaction Plate - Ready to use
Contents:	1 x 96 well plate, empty in a resealable pouch	
BA E-0030	WASH-CONC 50x	Wash Buffer Concentrate - Concentrated 50x
Contents:	Buffer with a non-ionic detergent and physiological pH	
Volume:	1 x 20 ml/vial, light purple cap	
BA E-1240	CONJUGATE	Enzyme Conjugate - Ready to use
Contents:	Donkey anti-goat immunoglobulins conjugated with peroxidase	
Volume:	1 x 12 ml/vial, red cap	
BA E-0055	SUBSTRATE	Substrate - Ready to use
Contents:	Chromogenic substrate containing tetramethylbenzidine, substrate buffer and hydrogen peroxide	
Volume:	1 x 12 ml/black vial, black cap	
BA E-0080	STOP-SOLN	Stop Solution - Ready to use
Contents:	0.25 M sulfuric acid	
Volume:	1 x 12 ml/vial, light grey cap	
BA E-1031	HIS	Histamine Microtiter Strips - Ready to use

Contents: 1 x 96 well (12x8) antigen precoated microwell plate in a resealable foil pouch with desiccant

BA E-1210 **HIS-AS** **Histamine Antiserum** - Ready to use

Contents: Goat anti-histamine antibody, blue coloured

Volume: 1 x 12 ml/vial, blue cap

Standards and Controls - Ready to use

Cat. no.	Component	Colour/Cap	Concentration ng/ml	Concentration nmol/l	Volume/Vial
BA E-1001	STANDARD A	white	0	0	4 ml
BA E-1002	STANDARD B	light yellow	0.5	4.5	4 ml
BA E-1003	STANDARD C	orange	1.5	13.5	4 ml
BA E-1004	STANDARD D	dark blue	5	45	4 ml
BA E-1005	STANDARD E	light grey	15	135	4 ml
BA E-1006	STANDARD F	black	50	450	4 ml
BA E-1051	CONTROL 1	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
BA E-1052	CONTROL 2	dark red			4 ml

Conversion: Histamine (ng/ml) x 9 = Histamine (nmol/l)

Contents: Acidic buffer spiked with defined quantity of Histamine

BA E-1211 **ACYL-BUFF** **Acylation Buffer** - Ready to use

Contents: TRIS-buffer

Volume: 2 x 12 ml/vial, brown cap

BA E-1212 **ACYL-REAG** **Acylation Reagent** - Ready to use

Contents: Acylation reagent containing DMSO

Volume: 2 x 1.5 ml/vial, green cap


BA E-0041 **DILUENT** **Diluent** - Ready to use


Contents: Acidic buffer with non-mercury preservatives

Volume: 1 x 22 ml/vial, white cap

4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 20 - 200 µl
- Microtiter plate washing device (manual, semi-automated or automated)
- ELISA reader capable of reading absorbance at 450 nm and if possible 620 - 650 nm
- Absorbent material (paper towel)
- Water (deionized, distilled, or ultra-pure)
- Vortex mixer

 The assay can be performed with or without the usage of a microtiter plate shaker. If a shaker is used, it should have the following characteristics: shaking amplitude 3 mm, capable of approx. 600 rpm.

 For the determination of total histamine in whole blood, a **Release Buffer** is necessary! (Catalog #BA E-1726)

BA E-1726 **RELEASE-BUFF** **Release Buffer** - Ready to use

Contents: Buffer with physiological pH

Volume: 1 x 250 ml/vial, white cap

5. Sample collection and storage

Cell culture

Samples can be stored at 2 – 8 °C for up to 6 hours. For longer periods (≤ 6 months) the samples should be stored between -20 °C and -80 °C.

Urine

Spontaneous urine or 24-hour urine, collected in a bottle containing 10 - 15 ml of 6 M HCl, may be used. *If 24-hour urine is being used please record the total volume of the collected urine.* If the percentage of the final concentration of acid is too high, the buffer capacity of the Diluent is insufficient. As a consequence interfering factors are not extracted quantitatively.

Storage: up to 6 hours at 2 - 8 °C; for longer periods (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

Avoid exposure to direct sunlight.

Whole Blood

Whole blood is collected into a tube (e.g. Monovette™ or Vacuette™) containing heparin as anti-coagulant. The samples can be stored for up to 24 hours at room temperature. Please do not keep the samples refrigerated, this will lead to clotting of the leucocytes. Avoid direct sunlight.

6. Test procedure

Allow all reagents and samples to reach room temperature and mix thoroughly by gentle inversion before use. Duplicate determinations are recommended.

The binding of the antisera and of the enzyme conjugate and the activity of the enzyme are temperature dependent, and the absorption values may vary if a thermostat is not used. The higher the temperature, the higher the absorption values will be. Varying incubation times will have similar influences on the absorbance. The optimal temperature during the Enzyme Immunoassay is between 20 - 25 °C.

6.1 Preparation of reagents

Wash Buffer

Dilute the 20 ml Wash Buffer Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 1000 ml.

Storage: 1 month at 2 - 8 °C

Acylation Reagent

The Acylation Reagent has a freezing point of 18.5 °C. To ensure that the Acylation Reagent is liquid when being used, it must be ensured that the Acylation Reagent has reached room temperature and forms a homogeneous, crystal-free solution.


Whole Blood

For the determination of total histamine in whole blood, dilute heparinized whole blood 1 + 20 with **Release Buffer** (catalog # BA E-1726) and incubate for **10 min** at **90 °C** (e.g. 50 µl whole blood plus 1 ml Release Buffer)


Cool down the samples for **10 min** at **2 - 8°C**

Centrifuge for 10 min at 700 x g (*the brake should be switched-off!*)

6.2 Sample preparation and acylation

1.	Pipette 100 µl of standards, controls and 20 µl of cell culture samples or urine samples and 100 µl of perpetrated whole blood samples (<i>refer to 6.1.</i>) into the respective wells of the Reaction Plate .
2.	Add 80 µl of Diluent to the wells with cell culture samples and urine samples .
3.	Add 25 µl of Acylation Reagent (<i>refer to 6.1</i>) to all wells.
4.	Pipette 200 µl of Acylation Buffer into all wells.
5.	Incubate 15 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Alternatively without shaker: <i>shake the Reaction Plate shortly by hand and incubate 15 min at RT (20 – 25 °C).</i>
	Take 25 µl for the subsequent ELISA

6.3 Histamine ELISA

1.	Pipette 25 µl of the acylated standards, controls and samples into the wells of the Histamine Microtiter Strips .
2.	Pipette 100 µl of the Histamine Antiserum into all wells.
3.	Incubate 30 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Alternatively without shaker: <i>shake the Histamine Microtiter Strips shortly by hand and incubate for 40 min at RT (20 – 25 °C).</i>
4.	Discard or aspirate the contents of the wells and wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
5.	Pipette 100 µl of the Enzyme Conjugate into all wells.
6.	Incubate for 10 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Alternatively without shaker: <i>incubate for 20 min at RT (20 – 25 °C).</i>
7.	Discard or aspirate the contents of the wells and wash the plate 3 x by adding 300 µl of Wash Buffer , discarding the content and blotting dry each time by tapping the inverted plate on absorbent material.
8.	Pipette 100 µl of the Substrate into all wells.
9.	Incubate for 15 ± 2 min at RT (20 – 25 °C) on a shaker (approx. 600 rpm). Alternatively without shaker: <i>incubate for 15 ± 2 min at RT (20 – 25 °C).</i>
	Avoid exposure to direct sunlight!
10.	Add 100 µl of the Stop Solution to each well and shake the microtiter plate shortly by hand to ensure a homogeneous distribution of the solution.
11.	Read the absorbance of the solution in the wells within 10 minutes, using a microplate reader set to 450 nm (if available a reference wavelength between 620 nm and 650 nm is recommended).

7. Calculation of results

Measuring range	Histamine	
	Controls	0.12 - 50 ng/ml
	Cell culture and urine	0.6 - 250 ng/ml
	Whole blood	2.52 - 1050 ng/ml

The standard curve is obtained by plotting the absorbance readings (calculate the mean absorbance) of the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis). Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

Note: This assay is a competitive assay. This means: the OD-values are decreasing with increasing concentrations of the analyte. OD-values found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.

Controls

The concentrations of the **controls** can be read directly from the standard curve.

Cell culture samples

The read concentrations of **histamine in cell culture samples** have to be **multiplied by 5**.

Whole blood samples

The read concentrations of **histamine in whole blood** have to be **multiplied by 21**.

Urine samples

The read concentrations of **histamine in urine** have to be **multiplied by 5**.

The total amount of Histamine excreted in urine during 24 h is calculated as following:

$$\mu\text{g}/24\text{h} = \mu\text{g}/\text{l} \times \text{l}/24\text{h}$$

Conversion

$$\text{Histamine (ng/ml)} \times 9 = \text{Histamine (nmol/l)}$$

Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.

In a study conducted with apparently normal healthy adults, using the Histamine ELISA ^{Fast Track} the following values were observed:

Whole blood (total Histamine)	24 hour-urine	Spontaneous urine
20 - 200 ng/ml	< 45 $\mu\text{g}/\text{d}$	< 45 $\mu\text{g}/\text{g}$ creatinine


7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and elevated levels. The kit or other commercially available controls should fall within established confidence limits. The confidence limits of the kit controls are listed in the QC-Report.

7.2 Typical standard curve



Example, do not use for calculation!



8. Assay characteristics

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Histamine
	3-Methyl-Histamine	0.1
	Tyramine	0.01
	L-Phenylalanine	< 0.001
	L-Histidine	< 0.001
	L-Tyrosine	< 0.001

	Tryptamine	< 0.001
	5-Hydroxy-Indole-Acetic Acid	< 0.001
	Serotonine	< 0.001

Analytical Sensitivity (Limit of Detection)	Histamine	0.2 ng/ml	Mean signal (Zero-Standard) + 2SD
---	-----------	-----------	-----------------------------------

Precision					
Inter-Assay Variation, n = 13			Intra-Assay Variation, n = 39		
Sample	Mean ± SD (ng/ml)	CV (%)	Sample	Mean ± SD (ng/ml)	CV (%)
1	2.03 ± 0.16	8	1	0.6 ± 0.1	12
2	6.74 ± 0.37	5.6	2	4.6 ± 0.3	6.3

Linearity		Range (ng/ml)	Range (%)	Mean (%)
	Histamine	0.74 - 8.48	85 - 106	100

Recovery		Serial dilution up to	Range (%)	Mean (%)
	Histamine	1:16	92 - 120	103

9. References/Literature










- (1) Yagci et al. TCTP/HRF pathway and angiogenesis: A feasible intercourse in chronic lymphocytic leukaemia. *Leukemia Research*, 37:665-670 (2013)
- (2) Coulson et al. Paracetamol (acetaminophen) attenuates in vitro mast cell and peripheral blood mononucleocyte cell histamine release induced by N-acetylcysteine. *Clinical Toxicology*, 48(2):111-114 (2010)
- (3) Rovere et al. Histamine and Selenium in Lung Cancer. *Anticancer Research*, 26: 2937-2942 (2006)

For orders, please contact:

Immuno-Biological Laboratories, Inc. (IBL-America)
 8201 Central Ave NE, Suite P
 Minneapolis, MN 55432
 Email: info@ibl-america.com
 Web: www.ibl-america.com
 Phone: (888) 523-1246
 Fax: (763) 780-2988

 **For updated literature or any other information please contact your local supplier.**

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled



Caution

REF

Catalogue
number

RUO

For research use
only!