

Code No. 54071

Fibronectin Neosilk®, Plasma

Volume : 1 mg

Lot No. :

Expiration date :

Introduction : A plasma fibronectin is present in plasma at high concentrations and plays roles in blood coagulation and wound healing. The fibronectin from plasma has been used as a cell-adhesive substrate (a coating material) for cultured cells¹. However, the fibronectin is derived from animal or human plasma, therefore, has risks of blood-borne infections.

Fibronectin Neosilk®, Plasma is a recombinant human fibronectin which is produced by transgenic silkworms as a homodimer of a plasma-abundant splicing variant, EDA-, EDB- and IIICS+². This recombinant fibronectin has the cell-adhesive activity comparable with or more than that of natural plasma fibronectin, and has a specific feature of free from animal derived materials. This product can be used as a coating material for various cultured cells including mesenchymal stem cells.

Content : Recombinant human fibronectin (EDA-, EDB-, IIICS+) homodimers.

Source : Extracted from cocoons of transgenic silkworms.

Purification : Purified by two types of multimodal chromatography (not by gelatin affinity chromatography).

Form : Lyophilized from 2 mL of 10 mM Tris-HCl, pH8.0 and 2% sucrose.

Storage : Lyophilized product should be stored at 2–10°C.

Reconstitution method : When dissolved in 2 mL of sterile purified water, the concentration of the fibronectin will be 0.5 mg/mL. After reconstitution, white coagulates may sometimes be observed in the solution, but these do not affect the protein concentration.

Storage after reconstitution : The reconstituted solution may be stored at 2–10°C for up to two months. For longer storage, aliquot into small volume, and store below –20°C. Avoid repeated freeze-thaw cycles. The solution is stable below –20°C for at least 6 months.

Example of coating method : Dilute the reconstituted solution with PBS (-) or other solution to appropriate concentrations, and coat the culture surface with the diluted solution. The typical coating concentration is 1–5 µg/cm². In the case of coating onto a 6 well plate at 2.5 µg/cm² (24 µg/well), 2 mL of the diluted fibronectin solution to 12 µg/mL is added to the wells. After incubation at 37°C for 1 hr, remove the solution, wash with PBS (-), and seed your cells.

References : 1. Parisi L. et. al., A glance on the role of fibronectin in controlling cell response at biomaterial interface. *Jpn. Dent. Sci. Rev.* 56, 50 (2020)
2. To WS, Midwood KS. Plasma and cellular fibronectin: distinct and independent functions during tissue repair. *Fibrogenesis Tissue Repair* 4, 21 (2011)

For Non-Clinical Research Use Only

Version 2

Distributed By:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com