

# Product information

Information about other products is available at: [www.IBL-America.com](http://www.IBL-America.com)



# Urine Reagent Set for Aldosterone



For Research Use Only – Not for Use in Diagnostic Procedures



DE5298URIN



96

***Please use only the valid version of the Instructions for Use provided with the kit.***

**Table of Contents**

1	INTENDED USE .....	3
2	PRINCIPLE OF THE FUNCTION.....	3
3	WARNINGS AND PRECAUTIONS .....	3
4	MATERIALS .....	4
5	SAMPLE COLLECTION, STORAGE AND PREPARATION.....	5
6	ASSAY PROCEDURE.....	5
7	QUALITY CONTROL.....	6
8	PERFORMANCE CHARACTERISTICS .....	6
9	LIMITATIONS OF THE PROCEDURE.....	6
10	LEGAL ASPECTS .....	6

## 1 INTENDED USE

The **IBL-America Urine Reagent Set for Aldosterone** is a reagent set intended to be used for the manual preparation of urine samples to be used with

- Aldosterone ELISA (IB79134; manual ELISA)

**For research use only – Not for use in diagnostic procedures.**

## 2 PRINCIPLE OF THE TEST

The IBL-America Urine Reagent Set for Aldosterone releases the binding of glucuronic acid from the aldosterone molecule using the Release Reagent to allow precise and true determinations. The addition of Neutralization Buffer is then used to neutralize the solution in order to allow for the assay to work as intended.

## 3 WARNINGS AND PRECAUTIONS

- For research use only – Not for use in diagnostic procedures.
- Before starting the usage, read the instructions for use completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact IBL-America.
- Microbial contamination of reagents or samples may give false results.
- Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
- All indicated volumes must be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into original vials as reagent contamination may occur.

### General precautions

- Follow good laboratory practice and safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.


### Biohazard information

- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

**Information to chemical hazards and hazard classification**

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Chemicals and prepared or used reagents must be treated as hazardous waste according to the national safety guideline or regulation.
- This product does not contain substances which have carcinogenic, mutagenic or toxic for reproduction (CMR) properties.

The following kit component(s) is(are) classified as hazardous:

<div style="display: flex; justify-content: space-between;"> <div style="border: 1px solid black; padding: 2px;">RELEASE</div> <div style="border: 1px solid black; padding: 2px;">REAG</div> </div> <p>Release Reagent</p>  <p>Warning</p>	<p><b>Hazard Statement(s):</b> H290 May be corrosive to metals.</p>
	<p><b>Precautionary statement(s):</b> P234 Keep only in original container P390 Absorb spillage to prevent material damage</p>

For detailed information, please refer to the Safety Data Sheet, which is available upon request from IBL-America.

**4 MATERIALS**

**4.1 Materials provided with the kit**

Symbol	Quantity	Description	Preparation
<div style="display: flex; justify-content: space-between;"> <div style="border: 1px solid black; padding: 2px;">RELEASE</div> <div style="border: 1px solid black; padding: 2px;">REAG</div> </div> <p>Release Reagent</p>	1 x 3 mL	<p><b>Release Reagent</b> Contains 1 M HCL <i>Avoid Contact with Release Reagent. It may cause skin irritation.</i></p>	Ready to use
<div style="display: flex; justify-content: space-between;"> <div style="border: 1px solid black; padding: 2px;">NEUTRA</div> <div style="border: 1px solid black; padding: 2px;">BUF</div> </div> <p>Neutralization Buffer</p>	1 x 3 mL	<p><b>Neutralization Buffer</b> Contains Tris buffer, pH 8.5</p>	Ready to use
<div style="display: flex; justify-content: space-between;"> <div style="border: 1px solid black; padding: 2px;">DIL</div> <div style="border: 1px solid black; padding: 2px;">BUF</div> </div> <p>Dilution Buffer</p>	2 x 25 mL	<p><b>Dilution Buffer</b> Contains PBS</p>	Ready to use
-	1 x	<b>Instructions for Use</b>	
-	1 x	<b>Certificate of Analysis (CoA)</b>	

**4.2 Materials required but not provided**

- Calibrated variable precision micropipettes
- Refrigerator for storage at 2 °C - 8 °C
- Distilled water
- Timer

### 4.3 Storage and Stability of the Kit

**Unopened kits and reagents** as well as **opened reagents** must be stored at 2 °C to 8 °C. Once opened, reagent vials must be closed tightly again.

	Storage Temperature	Stability
Unopened kits and unopened reagents	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!

### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

### 4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

### 4.6 Damaged Test Kits

In case of any damage to the test kit or components, IBL-America must be informed in writing, at the latest one week after receiving the kit. Damaged single components must not be used for a test run. They have to be stored until a final solution has been found. After this, they must be disposed of according to the official regulations.

## 5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

**urine**

Samples containing sodium azide should not be used in the assay.  
Please refer to the instructions for Use of Aldosterone ELISA (IB79134)

## 6 ASSAY PROCEDURE

### 6.1 Procedural Notes

- All reagents and samples must be allowed to come to room temperature (20 °C to 26 °C) before use.
- All reagents must be mixed without foaming.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- Once the test has been started, all steps must be completed without interruption and in the same sequence for each step.
- Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- **Test performance using fully automated analysis devices:**  
Automated test performance using fully automated, open-system analysis devices is possible. However, the combination must be validated by the user.

### 6.2 Test Procedure

Please refer to IFU of *Aldosterone ELISA (IB79134)*

### 6.3 Calculation of Results

Please refer to IFU of *Aldosterone ELISA (IB79134)*

## 7 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels. The controls and the corresponding results of the Quality Control Laboratory are stated in the Certificate of Analyses (CoA) added to the kit. The values and ranges stated on the CoA always refer to the current kit lot and must be used for direct comparison of the results. If available, it is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results. Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, results should be considered invalid. In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods. After checking the above-mentioned items without finding any error contact your distributor or IBL-America directly.

## 8 PERFORMANCE CHARACTERISTICS

Performance is solely depending on the performance of the respective device (*Aldosterone ELISA IB79134*)

## 9 LIMITATIONS OF THE PROCEDURE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice. Any improper handling of samples or modification of this test might influence the results.

## 10 LEGAL ASPECTS

### 10.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact IBL-America.






### 10.2 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.



**SYMBOLS USED WITH IBL-AMERICA ASSAYS**

Symbol	English	Deutsch	Française	Espanol	Italiano
	European Conformity	CE-Konformitäts-kennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	utilisation Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di catalogo
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributed by	Vertrieb durch	Distribution par	Distribución por	Distribuzione da parte di
	Version	Version	Version	Versión	Versione
	Single-use	Einmalverwendung	À usage unique	Uso único	Uso una volta