

Manufactured for Immuno-Biological Laboratories Inc.
(IBL-America)
8201 Central Avenue, NE, Suite P
Minneapolis, MN 55432
Tel: 763-780-2955
Toll Free: 1-888-523-1246



Instructions for use

Melatonin (Research) RIA

Radioimmunoassay for the determination of Melatonin in biological samples.

For research use only, not for use in diagnostic procedures.

REF

IB88139



RUO

For Research use only-
Not for use in diagnostic
procedures

200 kBq

Melatonin Research RIA

1. Intended use and principle of the test

Radioimmunoassay for the determination of Melatonin in biological samples. For research use only, not for use in diagnostic procedures.

In principle the Melatonin Research RIA kit will work for all kinds of biological samples. As a starting point it is recommended to perform a proof-of-principle (please contact our technical service for details).

As a reference a proof of principle for rat plasma samples was performed using 100 µl of sample volume (please refer to chapter 8.).

Melatonin - the major hormone secreted by the pineal gland - is a key modulator of annual and circadian biorhythms. Its circadian profile in body fluids is an excellent marker for the setting of the endogenous clock. Daytime plasma melatonin levels are low and rise in the evening (onset). Night-time levels peak at around 03.00 hrs. (acrophase) in most healthy humans.

Onset, acrophase and offset have a stable phase relationship even when the phase of the melatonin profile is shifted.

The assay procedure follows the basic principle of radioimmunoassay, involving competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of ¹²⁵I-labelled antigen bound to the antibody is inversely proportional to the analyte concentration of the sample. When the system is in equilibrium, the antibody bound radioactivity is precipitated with a second antibody in the presence of polyethylene glycol. The precipitate is counted in a gamma counter. Determination of unknown samples is achieved by comparing their activity with a reference curve prepared with known standards.

2. Precautions, Guidelines and Warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) The principles of Good Laboratory Practice (GLP) have to be followed.
- (3) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.
- (4) The radioactive material (¹²⁵Iodine, half life 60 days, emitting ionizing X-radiation with 28 keV and G-radiation with 35.5 keV) may be received, acquired, possessed and used only by physicians, laboratories or hospitals. In compliance with regulations, a copy of the customer's current radioisotope license must be on file with the supplier. Orders cannot be shipped until the license is received by the supplier (Radiation Protection Act of June 30, 1989).
- (5) For the dilution or reconstitution purposes use deionized, distilled, or ultra-pure water.
- (6) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (7) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (8) Incubation times do influence the results. All tubes should be handled in the same order and time intervals.
- (9) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (10) A calibrator curve must be established for each run.
- (11) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report.
- (12) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (13) For information on hazardous substances included in the kit please refer to Material Safety Data Sheets (MSDS). The Material Safety Data Sheet for this product is available directly on the website of the manufacturer or upon request.
- (14) If expected reference values are reported in this test instruction they are only indicative. It is recommended that each laboratory establishes its own reference intervals.

3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 - 8 °C.

Melatonin is sensitive to light-exposure. To avoid photo-oxidative reduction of melatonin, it is necessary to keep it away from direct sunlight and from heat.

4. Materials

4.1 Content of the kit

BA R-0030 PREC-REAG **Precipitating Reagent** - Ready to use

Content: Goat anti-rabbit serum in PEG phosphate buffer

Volume: 2 x 55 ml/vial, yellow cap

BA R-3310 AS MEL **Melatonin Antiserum** - Ready to use

Content: Rabbit anti-melatonin antibody, blue coloured

Volume: 1 x 5.25 ml/vial, blue cap

BA R-3315 ENZYME **Enzyme** - Lyophilized

Content: Digestive enzyme

Volume: 4 vials, pink cap

BA R-3313 ASSAY-BUFF **Assay Buffer** - Ready to use

Content: TRIS buffer

Volume: 1 x 15 ml/vial, light purple cap

Hazards identification:



H315 Causes skin irritation.

H319 Causes serious eye irritation.

BA R-3316 ENZYME-BUFF **Enzyme Buffer** - Ready to use

Content: 1 M hydrochloric acid

Volume: 1 x 15 ml/vial, orange cap

BA R-3320 ¹²⁵I-MEL **¹²⁵I – Melatonin** - Ready to use

Content: ¹²⁵I labeled Melatonin, red coloured

Volume: 1 x 3 ml/vial, red cap

Hazards identification:



Radioactive, activity < 200 kBq

Standards and Controls - Ready to use

Cat. no.	Component	Colour/Cap	Concentration pg/ml	Concentration pmol/l	Volume/ Vial
BA R-3301	STANDARD A	white	0	0	4 ml
BA R-3302	STANDARD B	light yellow	30	129	4 ml
BA R-3303	STANDARD C	orange	100	430	4 ml
BA R-3304	STANDARD D	dark blue	300	1 290	4 ml
BA R-3305	STANDARD E	light grey	1 000	4 300	4 ml
BA R-3306	STANDARD F	black	3 000	12 900	4 ml
BA R-3307	STANDARD G	brown	10 000	43 000	4 ml
BA R-3351	CONTROL 1	light green	Refer to QC-Report for expected value and acceptable range!		4 ml
BA R-3352	CONTROL 2	dark red			4 ml

Conversion: Melatonin (pg/ml) x 4.30 = Melatonin (pmol/l)

Content: TRIS buffer with non-mercury preservatives, spiked with defined quantity of melatonin

BA R-3961 SYRINGE **Syringe** - Ready to use, single-use

Volume: 1 x 3 pcs.

BA R-3962 FILTER UNITS **Filter Units** - Ready to use, single-use

Volume: 1 x 3 pcs.

BA R-3963 PREP TUBES **Preparation Tubes** - Ready to use, single-use

Volume: 1 x 3 pcs.

BA R-3964 CHARCOAL **Charcoal Suspension** - Ready to use

Content: activated charcoal

Volume: 1 x 25 ml/vials, light purple cap

4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 20 – 1000 µl; 1,8 ml, 3 ml, 3,6 ml, 7,2 ml
- Conical plastic tubes (polypropylene, polystyrene) and suitable rack
- Centrifuge (preferable refrigerated) capable of at least 3 000 x g
- Suitable device for aspirating or decanting the tubes
- Vortex mixer
- Gamma counter
- Water (deionized, distilled, or ultra-pure)

5. Sample storage

The samples have to be stored deep frozen until use.

6. Test procedure

Allow all reagents – with the exception of Precipitating Reagent - to reach room temperature and mix thoroughly by gentle inversion before use. Number the assay tubes accordingly. Duplicate determinations are recommended.



Pipetted liquids should not adhere to the wall of the RIA tubes. If necessary please centrifuge the tubes for 1 minute at 500 x g to spin down adhering liquids.



For the assay the use of conical tubes is highly recommended.

6.1 Preparation of reagents**Equalizing Reagent:**

The preparation of a species and sample specific Equalizing Reagent is mandatory to avoid any matrix effects and false results.

For the preparation of this specific reagent a pool of exact the composition as the sample should be used (e.g. for the measurement of Melatonin in rat plasma, rat plasma has to be used as a starting material for the preparation of the Equalizing Reagent). Through the addition of a charcoal suspension the endogenous melatonin is removed from the specific biological fluid and a Melatonin-free (= "stripped") Equalizing Reagent is produced.

The needed amount of Equalizing Reagent depends on how much sample volume is used for the experiment and this has a strong influence on the limit of detection (please refer to chapter 8.). The following examples demonstrate the relation between number of runs, sample volume, and needed Equalizing Reagent:

Number of runs (standard curves)	Sample volume (µl)	Tubes (N)*	Volume of Equalizing Reagent needed**
1	100	18	1.8 ml
2	100	36	3.6 ml
2	50	36	1.8 ml
2	200	36	7.2 ml

* The number of tubes needed to prepare a standard curve in duplicates including NSB, 7 standards and 2 controls.

** The charcoal treatment of biological fluids in combination with the ultrafiltration step always leads to a loss of some material. Therefore it is recommended to take more sample volume to prepare the Equalizing Reagent than is needed in the assay (e.g. twice as much).

Example: Preparation of Equalizing Reagent for 100 µl sample volume; run in duplicates

1.	Pipette 7.2 ml of Charcoal Suspension into the Preparation Tubes .
2.	Spin down the charcoal.
3.	Decant the supernatant and remove carefully remaining liquid from the wall of the tubes.
4.	Pipette 7.2 ml of the biological fluid (serum, tissue extract etc.) onto the charcoal pellet.
5.	Close the tube carefully and mix the suspension for 30 min at RT (20 - 25 °C) on a rotating mixer .
6.	Spin down the charcoal.
7.	Remove charcoal fines by filtering the supernatant. One Filter Unit is capable of filtering 4 ml of supernatant so that two Filter Units (and two Syringes) have to be used.
8.	Mix both filtrates and store deep-frozen.

Always check the volumes of the melatonin free Equalizing Reagent prior to its use in the RIA to ensure that the prepared volume is sufficient.

Enzyme

Reconstitute the content of the vial with 3 ml of Enzyme Buffer prior to use. Mix carefully (30 minutes on a rotating mixer). The reconstituted Enzyme cannot be stored and can only be used once. Upon request additional Enzyme vials are provided.

6.2 Melatonin RIA

(this exemplary protocol refers to 100 µl sample volume; for other sample volumes please refer to the chart 6.2.1 below)

1.	Pipette 20 µl of Standard A into the tubes for the NSB and into the tubes for the samples .
2.	Pipette 20 µl of standards and controls into the respective tubes.
3.	Add 100 µl of Equalizing Reagent to the tubes with NSB, standards and controls .
4.	Pipette 100 µl of samples into the respective tubes.
5.	Pipette 25 µl of Enzyme in all tubes (except totals) and vortex.
6.	Incubate for 1 h at RT (20 – 25 °C).
7.	Pipette 50 µl of Assay Buffer into all tubes (except totals) and mix shortly.
8.	Pipette 25 µl of the ¹²⁵I Melatonin into all tubes.
9.	Pipette 50 µl of Melatonin Antiserum into all tubes (except totals and NSB); mix thoroughly.
10.	Cover tubes and incubate for 20 - 24 h at RT (20 – 25 °C).
11.	Mix the chilled (2 - 8 °C) Precipitating Reagent thoroughly, pipette each 1000 µl into all tubes (except totals), and mix on a vortex.
12.	Incubate for 20 min at 2 - 8 °C .
13.	Centrifuge for 20 min at 3 000 x g , if possible in a refrigerated centrifuge.
14.	Decant or aspirate the supernatant carefully (except totals) . Blot the tubes dry and leave them upside down for 2 minutes.
15.	Count all tubes for 1 min in a gamma counter.

6.2.1 For sample volumes differing from 100 µl, the following modifications for the RIA (6.2) in reagent volumes are recommended:

Equalizing Reagent	Sample volume	Enzyme solution	Assay Buffer	Antiserum
100 µl	100 µl	25 µl	50 µl	50 µl
25 µl	25 µl	25 µl	50 µl	50 µl
50 µl	50 µl	25 µl	50 µl	50 µl
200 µl	200 µl	50 µl	100 µl	50 µl

7. Calculation of results

Subtract the mean cpm of the non-specific binding NSB from the mean cpm of standards, controls and samples.

The standard curve from which the concentrations in the samples can be read off, is obtained by plotting the percentage of (B-NSB)/ (B₀-NSB) measured for the standards (linear, y-axis) against the corresponding standard concentrations (logarithmic, x-axis).

Use a non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).



This assay is a competitive assay. This means: the counts are decreasing with increasing concentrations of the analyte. Counts found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.

The concentrations of the **controls** can be read directly from the respective standard curve.

The concentrations of the **samples** are depending on the sample volume which is used for the assay and concentrations read from the standard curve have to be **multiplied with a volume-factor**:

$$\text{Volume factor} = \frac{20}{\text{sample volume used for the assay}}$$

Conversion

Melatonin (pg/ml) x 4.30 = Melatonin (pmol/l)

7.1 Quality control

The confidence limits of the kit controls are indicated on the QC-Report.

8. Assay characteristics (for rat plasma, 100 µl sample volume)

Analytical Sensitivity (Limit of Detection)	5.9 pg/ml
--	-----------

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)
		Melatonin
	Melatonin	100
	N-Acetylserotonin	0.98
	5-Methoxytryptophol	0.11
	5-Methoxytryptamine	0.07
	6-Methoxytryptamine	< 0.01
	5-Methoxyindol-3-acetic acid	< 0.01
	Serotonin	< 0.01
	DL-Tryptophan	< 0.01
	DL-5-Methoxytryptophan	< 0.01
	5-Hydroxy-L-Tryptophan	< 0.01

Linearity		Range	Serial dilution up to	Range (%)
	Melatonin Rat	126 – 2 128 pg/ ml	1:128	83 – 98

Recovery		Mean (%)	Range (%)	% Recovery after spiking
	Melatonin Rat	94	82 – 103	

Stability (Freeze and Thaw Stability)	Stability (comparison Fresh Sample with Freeze and Thaw Stability)						
	Sample	Fresh Sample	Freeze and Thaw Stability	Mean (pg/ml)	Deviation (%)	SD (%)	CV (%)
	1	39.5 pg/ml	40 pg/ml	39.8	1.2	0.4	0.9
	2	87.1 pg/ml	89.3 pg/ml	88.2	2.5	1.6	1.8
	3	248.6 pg/ml	240.4 pg/ml	244.5	3.3	5.8	2.4
	Stability (comparison Fresh Sample with Short-Term Temperature Stability)						
	Sample	Fresh Sample	Freeze and Thaw Stability	Mean (pg/ml)	Deviation (%)	SD (%)	CV (%)
	1	39.5 pg/ml	46.4 pg/ml	43	17.5	4.9	11.3
	2	87.1 pg/ml	85.9 pg/ml	86.5	1.4	0.8	1.0
	3	248.6 pg/ml	210.8 pg/ml	229.7	15.2	26.7	11.6

For orders, please contact:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave NE, Suite P
Minneapolis, MN 55432
Email: info@ibl-america.com
Web: www.ibl-america.com
Phone: (888) 523-1246
Fax: (763) 780-2988



For literature or any other information please contact your local supplier.

Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled
	Caution		Catalogue number		For research use only!