# **CA 19-9-IRMA**

***KIP0311***

For Informational/Research Purposes Only

**Distributed By:**

Immuno-Biological Laboratories, Inc. (IBL-America)  
8201 Central Ave. NE, Suite P  
Minneapolis, Minnesota 55432, USA  
Phone: (888) 523-1246  
Fax.: (763) 780-2988  
Email: [info@ibl-america.com](mailto:info@ibl-america.com)  
Web: [www.ibl-america.com](http://www.ibl-america.com)

For Informational/Research Purposes Only

---

**LOT** : 110218/1

Read entire protocol before use.

## CA 19-9-IRMA

### *I. INTENDED USE*

Immunoradiometric assay kit for the in vitro quantitative determination of the cancer associated antigen CA 19-9 in serum.

### *II. GENERAL INFORMATION*

- A. **Proprietary name :** DIAsource CA 19-9-IRMA Kit
- B. **Catalog number :** KIP0311 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.90

### *III. CLINICAL BACKGROUND*

CA 19-9, equivalent to Sialyl Lewis<sup>a</sup> antigen, is a mucin type glycoprotein found in the blood repetitively bearing a carbohydrate epitope which is specifically recognized by the mouse monoclonal antibody, C192.

CA 19-9 is an oncofetal antigen, expressed by several different cancers, but especially carcinomas of the gastrointestinal tract.

It is present in the fetal epithelium of the stomach, intestine, liver and pancreas.

CA 19-9 is a tumor marker for pancreatic, gallbladder, gastric and colorectal cancers.

CA 19-9 can be used to track treatment progress of previously diagnosed cancer. This blood test in conjunction with periodic CT scans will show whether the cancer is in remission or is continuing to grow.

However, it is more useful in measuring the effectiveness of cancer treatment by studying the patient's CA 19-9 levels over time. There are benign diseases of the bile ducts and pancreas that occasionally result in elevated CA 19-9 levels and a low CA 19-9 level does not exclude the possible presence of a malignant tumor. For this reason, the determination of CA 19-9 concentrations is not considered as a diagnostic test.

#### IV. PRINCIPLES OF THE METHOD

The DIASource CA 19-9-IRMA is a two-step immunoradiometric assay based on coated-tube separation. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Add calibrators or samples to the tubes. After incubation, washing removes the occasional excess of antigen. Addition of Mab2, the signal antibody labelled with  $^{125}\text{I}$ , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration. .

#### V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with Avidine	2 x 48	silver	Ready for use			
<table border="1" data-bbox="119 571 255 622"> <tr> <td>Ab</td> <td><math>^{125}\text{I}</math></td> </tr> </table> Anti-CA 19-9- $^{125}\text{I}$ (monoclonal antibody) in phosphate-citrate buffer with bovine serum albumin, bovine serum, azide (<0.1%) and inert red dye	Ab	$^{125}\text{I}$	1 vial 10,5 ml 760 kBq	red	Ready for use	
Ab	$^{125}\text{I}$					
<table border="1" data-bbox="119 772 255 801"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Calibrators 0 in Tris Maleate buffer with bovine serum albumin and azide (<0.1%).	CAL	0	1 vials 10 ml	yellow	Ready for use	
CAL	0					
<table border="1" data-bbox="119 907 255 936"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators 1-5 in Tris Maleate buffer with bovine serum albumin, bovine serum and thymol (<0.1%). See exact value on vial labels.	CAL	N	5 vials lyophilised	yellow	Add 0.5 ml distilled water	
CAL	N					
<table border="1" data-bbox="119 1064 279 1093"> <tr> <td>Ab</td> <td>BIOT</td> </tr> </table> Biotin Anti-CA 19-9 (monoclonal antibody) : phosphate buffer with bovine serum albumin, bovine serum, azide (<0.1%) and inert blue dye	Ab	BIOT	2 vials 3 ml	bleu	Ready for use	
Ab	BIOT					
<table border="1" data-bbox="71 1254 319 1288"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (NaCl, Tween 20)	WASH	SOLN	CONC	1 vial 40 ml	green	Dilute 20 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="71 1355 279 1388"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum and thymol	CONTROL	N	2 vials lyophilised	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note: use the zero calibrator for sera dilutions.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 25  $\mu\text{l}$ , 50  $\mu\text{l}$ , 100  $\mu\text{l}$  and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Pipette for delivery of 1 to 5 ml of distilled water
- Vortex mixer
- Magnetic stirrer
- Tube shaker (400rpm)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- Calibrators:** Reconstitute the calibrator 1-5 with 0.5 ml distilled water.
- Controls:** Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 19 volumes of distilled water to 1 volume of Wash Solution (20x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- The calibrators and controls are very unstable, use them immediately after reconstitution, freeze immediately in aliquots and keep them at -20°C for maximum 3 months.
- Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Serum samples must be kept at 2-8°C.
- If the test is not run within 24 hrs, storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.

#### X. PROCEDURE

##### A. Handling notes

- Do not use the kit or components beyond expiry date.
- Do not mix materials from different kit lots.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure.

- Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes.
- Briefly vortex calibrators, samples and controls and dispense 25  $\mu\text{l}$  of each into the respective tubes.
- Dispense 50  $\mu\text{l}$  Biotin Anti-CA 19-9 in each tube except those for total counts.
- Incubate for 30 minutes at room temperature on a tube shaker (400 rpm).
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Dispense 100  $\mu\text{l}$  of  $^{125}\text{I}$  Iodine labelled anti CA19-9 into each tube, including the uncoated tubes for total counts.
- Incubate for 60 minutes at room temperature on a tube shaker(400 rpm).
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CA 19-9 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended

## XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

CA 19-9-IRMA		cpm	B/T (%)
Total count		285932	100
Calibrator	0.0 U/ml	212	0.1
	7.8 U/ml	1459	0.5
	26.4 U/ml	5196	1.8
	100.6 U/ml	17489	6.1
	310.5 U/ml	45032	15.7
	1061.1 U/ml	87314	30.5

Since no international reference material is available for CA 19-9 antigen, DIAsource CA 19-9-IRMA calibrator values are assigned against a set of in-house reference standards.

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.7 U/ml

### B. Precision

#### INTRA-ASSAY PRECISION

#### INTER-ASSAY PRECISION

Serum	N	<X> ± SD (U/ml)	CV (%)	Serum	N	<X> ± SD (U/ml)	CV (%)
A	20	23.5 ± 0.6	2.7	A	20	57.7 ± 2.9	5.0
B	20	194.3 ± 5.3	2.7	B	20	201.6 ± 7.4	3.7

SD: Standard Deviation; CV: Coefficient of variation

### C. Accuracy

#### DILUTION TEST

Sample	Dilution	Theoretical Concent. (U/ml)	Measured Concent. (U/ml)
A	1/1	-	661.1
	1/2	330.6	338.3
	1/4	165.3	169.2
	1/8	82.6	91.3
	1/16	41.3	46.5
	1/32	20.7	21.8
	1/64	10.3	10.2
	1/128	5.2	5.1

Samples were diluted with calibrator 0.

#### RECOVERY TEST

Added CA 19-9 (U/ml)	Recovered CA 19-9 (U/ml)	Recovery (%)
62.50	58.6	93.8
125.00	137.0	109.6
250.00	255.9	102.4
500.00	498.1	99.6

### D. Hook-effect

Samples containing 100.000 U/ml CA 19-9 give a result higher than the last calibration point.

### E. Specificity

The DIAsource CA 19-9 IRMA is based on one mouse monoclonal antibody, C192, highly specific for the sialyl Lewis epitope.

## XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely

elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.

- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

## XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices.

## XVI. REFERENCE INTERVALS

65 apparently healthy individuals, 97% were found to have serum CA 19-9 values of less than 37.0 U/ml.

**For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.**

## XVII. PRECAUTIONS AND WARNINGS

### Safety

For *in vitro* diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days), emitting ionizing X (28 keV) and  $\gamma$  (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

## XVIII. BIBLIOGRAPHY

1. KOURI M., PYRHONON S., KUUSELA P.  
**Elevated CA 19-9 as the most significant prognostic factor in advanced colorectal carcinoma.**  
J Surg Oncol 1992;49(2):78-85.
2. KOURI M., NORLING S., PYRHONON S., KUUSELA P.  
**Poor prognosis associated with elevated serum CA 19-9 level in advanced colorectal carcinoma, independent of DNA ploidy or SPF.**  
Eur J Cancer 1993;29(A):1691-1696.
3. FREBOURG T., BERCOFF E., MANCHON N., SENANT J., BASUYAU J-P., BRETON P., JANVRESSE A., BRUNELLE P., BOURREILLE J.  
**The evaluation of CA 19-9 antigen level in the early detection of pancreatic cancer. A prospective study of 866 patients.**  
Cancer 1988;62(11):2287-2290.

4. SAFI F., SCHLOSSER W., FALKENRECK S., BERGER HG.  
**CA 19-9 serum course and prognosis of pancreatic cancer.**  
Int J Pancreatol 1996;20(3):155-161.
5. MO?TGOMERY RC., HOFFMAN JP., RILEY LB., ROGATKO A., RIDGE JA., EISENBERG BL.  
**Prediction of recurrence and survival by post-resection CA 19-9 values in patients with adenocarcinoma of the pancreas.**  
Ann Surg Oncol 1997;4(7):551-556.
6. NAKAO A., OSHIMA K., NOMOTO S., TAKEDA S., KANEKO T., ICHIHARA T., KUROKAWA T., NONAMI T., TAKAGI H.  
**Clinical usefulness of CA 19-9 in pancreatic carcinoma.**  
Semin Surg Oncol 1998;15(1):15-22.
7. Mc LAUGHLIN R., O'HANLON D., KERIN M., KENNY P., GRIMES H., GIVEN HF.  
**Are elevated levels of the tumour marker CA 19-9 of any clinical significance ?**  
Ir J Med Sci 1999;168(2):124-126.
8. KODERA Y., YAMAMURA Y., TORII A., UESAKA K., HIRAI T., YASUI K., MORIMOTO T., KATO T., KITO T.  
**The prognostic value of preoperative serum levels of CEA and CA 19-9 in patients with gastric cancer.**  
Am J Gastroenterol 1996;91(1):49-53.
9. GOETZ M., STEEN PD.  
**Flase elevation of CA 19-9 levels in a patient with a history of pancreatic cancer.**  
Am J Gastroenterol 1997;98(8):1390-1391

**XIX. SUMMARY OF THE PROTOCOL**

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
Calibrators (0 to 5)	-	25	-
Sample(s), controls	-	-	25
Biotin Anti-CA19-9	-	50	50
Incubation	30 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Tracer	100	100	100
Incubation	60 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Working Wash solution		2 ml	
Separation		Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP0311	P.I. Number : 1700462/en	Revision nr : 110218/1
-------------------------------------	-----------------------------	---------------------------

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	<b>Used symbols</b>
	Consult instructions for use
	Storage temperature
	Use by
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>CONTROL</b>	Control
<b>I V D</b>	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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