CA 15-3-IRMA

KIP0321

For Informational/Research Purposes Only

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

For Informational/Research Purposes Only

LOT : 110218/1

Read entire protocol before use.

CA 15-3-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative determination of the cancer associated antigen MUC-1 (CA 15-3) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource CA 15-3-IRMA Kit
- B. **Catalog number :** KIP0321 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.90

III. CLINICAL BACKGROUND

A. Biological activities

The CA 15-3 antigen is a high-molecular-weight mucin glycoprotein termed MUC-1. It has been characterized with monoclonal antibodies developed against purified extracts of a membrane fraction isolated from a human breast carcinoma (clone DF3) and against human milk fat globule membrane (clone 115D8).

Other names for the mucin include polymorphic epithelial mucin (PEM), epithelial membrane antigen (MEA) or episalin.

B. Clinical applications

This marker is most useful in evaluating the effect of treatment for women with advanced breast cancer. Elevated levels of CA15-3 are also associated with cancers of the ovary, lung and prostate, as well as non-cancerous conditions such as benign breast or ovarian disease, endometriosis, pelvic inflammatory disease and hepatitis. Pregnancy and lactation also can raise CA 15-3 levels.

IV. PRINCIPLES OF THE METHOD

The DIAsource CA 15-3-IRMA is a two-step immunoradiometric assay based on coated-tube separation. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Add calibrators or samples to the tubes. After incubation, washing removes the occasional excess of antigen. Addition of Mab2, the signal antibody labelled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration. .

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with anti CA 15-3(monoclonal antibody)	2 x 48	ivory	Ready for use			
<table border="1" data-bbox="124 577 258 629"> <tr> <td>Ab</td> <td>^{125}I</td> </tr> </table> Anti-CA 15-3- ^{125}I (monoclonal antibody) in phosphate buffer with bovine serum, azide (<0.1%) and inert red dye	Ab	^{125}I	1 vial 6.0 ml 870 kBq	red	Ready for use	
Ab	^{125}I					
<table border="1" data-bbox="124 779 258 831"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators 1-5 in phosphate buffer with bovine serum albumin and thymol (<0.1%). See exact value on vial labels. Calibrators are prediluted	CAL	N	5 vials lyophilised	yellow	Add 0.5 ml dilution buffer	
CAL	N					
<table border="1" data-bbox="113 992 280 1043"> <tr> <td>DIL</td> <td>BUF</td> </tr> </table> Dilution buffer: phosphate buffer with bovine serum and azide (<0.1%)	DIL	BUF	1 vial 26 ml	black	Ready for use	
DIL	BUF					
<table border="1" data-bbox="76 1155 323 1207"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (TRIS-HCl)	WASH	SOLN	CONC	2 vials 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="82 1267 280 1319"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum and thymol	CONTROL	N	2 vials lyophilised	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note: use the dilution buffer as zero calibrator.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 μl , 200 μl and 500 μl (the use of accurate pipettes with disposable plastic tips is recommended)
3. Plastic tubes for dilution of samples
4. Vortex mixer
5. Tube shaker (400 rpm)
6. Magnetic stirrer
7. 5 ml automatic syringe (Cornwall type) for washing
8. Aspiration system (optional)
9. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. Calibrators:** Reconstitute the calibrators with 0.5ml dilution buffer.
- B. Controls:** Reconstitute the controls with 0.5 ml distilled water.
- C. Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- The calibrators and controls are very unstable, use them immediately after reconstitution, freeze immediately in aliquots and keep them at -20°C for 3 months.
- Avoid subsequent freeze-thaw cycles.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum samples must be kept at 2-8°C.
- If the test is not run within 24 h., storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.

X. PROCEDURE

A. Handling notes

- Do not use the kit or components beyond expiry date.
- Do not mix materials from different kit lots.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label one plain plastic tube for each sample and control.
2. Dispense 500 μl of Dilution Buffer into each tube.
3. Add 20 μl of sample and controls into these tubes.
4. Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes.
5. Briefly vortex calibrators, pre-diluted samples and controls and dispense 50 μl of each into the respective tubes.
6. Incubate for 90 minutes at room temperature on a tube shaker (400 rpm).
7. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
8. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
9. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
10. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
11. Dispense 50 μl of ^{125}I Iodine labelled anti CA15-3 into each tube, including the uncoated tubes for total counts.
12. Incubate for 90 minutes at room temperature on a tube shaker(400 rpm).
13. Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
14. Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
15. Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
16. Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
17. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CA 15-3 (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

CA 15-3-IRMA		cpm	B/T (%)
Total count		258080	100
Calibrator	0.0 U/ml	415	0.16
	2.86 U/ml	2440	0.95
	9.10 U/ml	6957	2.70
	28.30 U/ml	18367	7.12
	104.00 U/ml	48032	18.61
	331.00 U/ml	89798	34.79

Since no international reference material is available for CA 15-3 antigen, DIAsource CA 15-3-IRMA calibrator values are assigned against a set of in-house reference standards.

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.33 U/ml.

B. Precision

INTRA-ASSAY PRECISION

INTER-ASSAY PRECISION

Serum	N	<X> ± SD (U/ml)	CV (%)	Serum	N	<X> ± SD (U/ml)	CV (%)
A	20	14.2 ± 0.4	2.7	A	20	12.8 ± 0.8	6.4
B	20	94.1 ± 2.8	3.0	B	20	86.9 ± 5.4	6.2

SD: Standard Deviation; CV: Coefficient of variation

C. Accuracy

DILUTION TEST

Sample	Dilution	Theoretical Concent. (U/ml)	Measured Concent. (U/ml)
A	1/1	-	184.61
	1/2	92.31	91.79
	1/4	46.15	46.70
	1/8	23.08	22.33
	1/16	11.54	11.34
	1/32	5.77	5.61

Samples were diluted with Dilution Buffer.

RECOVERY TEST

Added CA 15-3 (U/ml)	Recovered CA 15-3 (U/ml)	Recovery (%)
21.80	20.48	93.9
43.75	43.34	99.1
80.91	81.59	100.8
169.92	164.18	96.6
307.87	299.60	97.3

D. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY				
	0'	10'	20'	30'
S 1 (U/ml)	11.4	11.4	11.4	11.2
S 2 (U/ml)	187.7	183.3	187.2	187.2

E. Hook-effect

Samples containing 100.000 U/ml CA15-3 give a result higher than the last calibration point.

F. Specificity

The DIAsource CA 15-3 IRMA is based on two mouse monoclonal antibodies, the catching MAb Ma695 targeting a sialylated carbohydrate epitope and the detecting MAb Ma 552 targeting the PDTRPAG hydrophilic peptide of the protein core.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices.

XVI. REFERENCE INTERVALS

Among 232 apparently healthy individuals, 99% of the results were below 30.0 U/ml.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in All

animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

- GENDLER SJ, SPICER AP. +
Epithelial mucin genes.
Ann Rev Physiol 1995;57:112-120.
- KUFE D, INGHIRAMI G, ABE M, et al.
Differential reactivity of a novel monoclonal antibody with human malignant versus benign breast tumors.
Hybridoma 1984;3:223-232.
- DUFFY MJ
CA 15-3 and related mucins as circulating markers in breast cancer.
Ann Clin Biochem 1999;36:579-586.
- COVENEY EC, GERAGHTY JG, et al.
The clinical value of CEA and CA 15-3 I breast cancer management.
Int J Biol Markers 1995;10:35-41.
- SOLETORMOS G, NIELSEN D, et al.
CA 15-3, CEA and TPA for monitoring metastatic breast cancer during first line chemotherapy and follow-up.
Clin Chem 1996;42:564-575.
- CHEUNG K, GRAVES CR, et al.
Tumor marker measurements in the diagnosis and monitoring of breast cancer.
Cancer Treat Rev 2000;26:91-102.
- MOLINA R, DUFFY MJ, et al.
Tumor markers in breast cancer : EGTM Recommendations.
Anticancer Res 1999;19:2803-2805.
- GION M, PELOSO L, et al.
Tumor markers in breast cancer monitoring should be scheduled according to initial stage and follow-up time.
Cancer J 2001;7:181-190.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	DILUTED SAMPLE CONTROLS µl
Calibrators (0 to 5)	-	50	-
Diluted Samples, controls	-	-	50
Incubation	90 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant) 2 ml	
Working Wash solution		Aspirate (or decant) 2 ml	
Separation		Aspirate (or decant) 2 ml	
Working Wash solution		Aspirate (or decant) 2 ml	
Separation		Aspirate (or decant)	
Tracer	50	50	50
Incubation	90 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant) 2 ml	
Working Wash solution		Aspirate (or decant) 2 ml	
Separation		Aspirate (or decant) 2 ml	
Working Wash solution		Aspirate (or decant) 2 ml	
Separation		Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP0321	P.I. Number : 1700864/en	Revision nr : 110218/1
-------------------------------------	-----------------------------	---------------------------

Revision date : 2011-02-18

For Informational/Research Purposes Only

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XVIII. BIBLIOGRAPHIE

- GENDLER SJ, SPICER AP.+
Epithelial mucin genes.
Ann Rev Physiol 1995;57:112-120.
- KUFE D, INGHIRAMI G, ABE M, et al.
Differential reactivity of a novel monoclonal antibody with human malignant versus benign breast tumors.
Hybridoma 1984;3:223-232.
- DUFFY MJ
CA 15-3 and related mucins as circulating markers in breast cancer.
Ann Clin Biochem 1999;36:579-586.
- COVENEY EC, GERAGHTY JG, et al.
The clinical value of CEA and CA 15-3 I breast cancer management.
Int J Biol Markers 1995;10:35-41.
- SOLETORMOS G, NIELSEN D, et al.
CA 15-3, CEA and TPA for monitoring metastatic breast cancer during first line chemotherapy and follow-up.
Clin Chem 1996;42:564-575.
- CHEUNG K, GRAVES CR, et al.
Tumor marker measurements in the diagnosis and monitoring of breast cancer.
Cancer Treat Rev 2000;26:91-102.
- MOLINA R, DUFFY MJ, et al.
Tumor markers in breast cancer : EGTM Recommendations.
Anticancer Res 1999;19:2803-2805.

- GION M, PELOSO L, et al.

Tumor markers in breast cancer monitoring should be scheduled according to initial stage and follow-up-time.

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XVIII. LITERATUR

- GENDLER SJ, SPICER AP.+
Epithelial mucin genes.
Ann Rev Physiol 1995;57:112-120.
- KUFE D, INGHIRAMI G, ABE M, et al.
Differential reactivity of a novel monoclonal antibody with human malignant versus benign breast tumors.
Hybridoma 1984;3:223-232.
- DUFFY MJ
CA 15-3 and related mucins as circulating markers in breast cancer.
Ann Clin Biochem 1999;36:579-586.
- COVENEY EC, GERAGHTY JG, et al.
The clinical value of CEA and CA 15-3 I breast cancer management.
Int J Biol Markers 1995;10:35-41.
- SOLETORMOS G, NIELSEN D, et al.
CA 15-3, CEA and TPA for monitoring metastatic breast cancer during first line chemotherapy and follow-up.
Clin Chem 1996;42:564-575.
- CHEUNG K, GRAVES CR, et al.
Tumor marker measurements in the diagnosis and monitoring of breast cancer.
Cancer Treat Rev 2000;26:91-102.
- MOLINA R, DUFFY MJ, et al.
Tumor markers in breast cancer : EGTM Recommendations.
Anticancer Res 1999;19:2803-2805.
- GION M, PELOSO L, et al.
Tumor markers in breast cancer monitoring should be scheduled according to initial stage and follow-up time.
Cancer J 2001;7:181-190.

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1. GENDLER SJ, SPICER AP.+
Epithelial mucin genes.
Ann Rev Physiol 1995;57:112-120.
2. KUFU D, INGHIRAMI G, ABE M, et al.
Differential reactivity of a novel monoclonal antibody with human malignant versus benign breast tumors.
Hybridoma 1984;3:223-232.
3. DUFFY MJ
CA 15-3 and related mucins as circulating markers in breast cancer.
Ann Clin Biochem 1999;36:579-586.
4. COVENEY EC, GERAGHTY JG, et al.
The clinical value of CEA and CA 15-3 I breast cancer management.
Int J Biol Markers 1995;10:35-41.
5. SOLETORMOS G, NIELSEN D, et al.
CA 15-3, CEA and TPA for monitoring metastatic breast cancer during first line chemotherapy and follow-up.
Clin Chem 1996;42:564-575.
6. CHEUNG K, GRAVES CR, et al.
Tumor marker measurements in the diagnosis and monitoring of breast cancer.
Cancer Treat Rev 2000;26:91-102.
7. MOLINA R, DUFFY MJ, et al.
Tumor markers in breast cancer : EGTM Recommendations.
Anticancer Res 1999;19:2803-2805.
8. GION M, PELOSO L, et al.
Tumor markers in breast cancer monitoring should be scheduled according to initial stage and follow-up time.
Cancer J 2001;7:181-190.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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