CT-U.S.-Irma

KIP0429

For Informational/Research Purposes Only

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

For Informational/Research Purposes Only

LOT : 120914/2

Read entire protocol before use.

CT-U.S. - IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative measurement of human Calcitonin in serum.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource CT-U.S. -IRMA Kit
- B. Catalog number :** KIP0429: 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

Calcitonin(CT) is a 32 amino acid peptide hormone secreted by the para-follicular C-cells of the thyroid gland under serum calcium control. After acute administration this peptide acts as a potent hypocalcemic and hypophosphatemic hormone by increasing renal calcium clearance and reducing bone resorption. However its precise physiological role in bone metabolism is not yet fully understood.

Various forms of CT may be detected in blood samples, including a CT monomer, an oxidized monomer, a dimer, higher molecular weight forms, and possibly precursor of CT. The concentrations of these peptides vary with clinical status, renal function and tissular origin of CT (normal or ectopic production).

Medullary thyroid carcinoma (MTC) is a malignant tumor, developed from the C-cells, secreting calcitonin in large excess. This disease occurs either as a sporadic (80%) or a familial (20%) form, which is transmitted as an autosomal dominant gene or as a component of multiple endocrine neoplasia (IIb).

Moderate hypercalcitoninemia is also observed in pregnancy, pernicious anemia, renal failure and during the neonatal period. Preferably, monomeric form of CT is detected in this assay.




The measurement of CT by the present IRMA is used for :

- diagnosis of medullary thyroid carcinoma (MTC)
- follow up of malignant tumors, to check the success of surgery and to monitor for recurrence
- diagnosis of the preclinical cases of the familial forms of MTC (MEN II or Sipple syndrome) by the use of stimulation tests (calcium or pentagastrin)
- study of the pathophysiology of the calcium-phosphate and bone metabolism.

IV. PRINCIPLES OF THE METHOD

The DIASource Calcitonin-U.S. IRMA is an immunoradiometric assay based on coated tube separation. It uses monoclonal antibodies directed against distinct epitopes of human Calcitonin. The capture antibody is attached to the lower inner surface of a plastic tube. Calibrators and samples are dispensed into the tubes and bind to the capture antibody. The signal antibody (radiolabelled) is added immediately and after an incubation which allows the immunologic reaction, the tubes are emptied and washed to remove the excess unbound antibody. The radioactivity bound to the tubes is directly proportional to the initial antigen concentration in the calibrators and samples.

V. REAGENTS PROVIDED

Reagents	Quantity	Colour Code	Reconstitution
 Tubes coated with anti CT (monoclonal antibodies)	2 x 48	blue	Ready for use
 Anti-CT- ¹²⁵ I (monoclonal antibodies) in TRIS Buffer with bovine serum albumin, azide (< 0.1%) and inert red dye	1 vial 5.5 ml 720 kBq	red	Ready for use
 Calibrators 0-5 in Calcitonin free human serum with gentamycin and thymol (see exact values on vial labels)	6 vials lyophil.	yellow	Add 1 ml distilled water
 CT free human serum (to be used for samples dilution) with gentamycin and thymol	1 vial lyophil.	black	Add distilled water (see the volume on the label)
 Wash solution (TRIS-HCl)	1 vial 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
 Controls 1 and 2 in human serum with gentamycin and thymol	2 vials lyophil.	silver	Add 1 ml distilled water

Note : 1. CT free human serum is to be used for samples dilution.
2. 1 pg of our reference preparation is equivalent to 0.19 μ IU 2nd IS 89/620.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 μ l, 200 μ l and 1 ml. (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. 5 ml automatic syringe (Cornwall type) for washing
6. Aspiration system (optional).
7. Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. Calibrators :** Reconstitute the calibrators with 1 ml distilled water.
- B. Controls :** Reconstitute the controls with 1 ml distilled water.
- C. CT free serum :** Reconstitute the CT free serum with distilled water. (see the volume on the label)
- D. Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators, controls and CT free serum should be frozen immediately after use and kept at -20°C for 3 months. Only one freeze thawing cycle is allowed, discard the calibrators, controls and CT free serum after the second use.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

Serum samples are recommended for this assay.

- Do not use hemolyzed samples.
- Do not use lipemic samples.
- If a specimen is expected or known to have a concentration above the highest calibrator, it has to be diluted with the CT free serum to fall within the measuring interval.
- If samples are not assayed the same day as the blood collection, then it is advisable to freeze them until the assay.
- Samples can only be thawed once.
- For repeat testing, freeze them in aliquots and discard each sample after first thawing.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

1. Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes
2. Homogenize calibrators, controls, specimens and dispense 200 μ l of each into the respective tubes.
3. Add 50 μ l of anti-CT-¹²⁵I (tracer) to all tubes, including the uncoated tubes for total counts.
4. Shake the tube rack gently by hand.
5. Incubate for 18 \pm 1 hours at 2-8°C.
6. Take the tubes for the total counts apart and aspirate the contents of the coated tubes. Be sure to remove all the liquid, remaining droplets will increase the background c.p.m..
7. Wash tubes twice with 2 ml Wash Solution and aspirate. Avoid foaming during the addition of the wash solution. For manual washing procedure first aspirate the foam layer and then the liquid. For automated wash cycles continuous aspiration is recommended. The precision is improved by aspirating the tubes a second time two minutes after emptying the last tube.
8. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of CT (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve. Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4 parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

CT-U.S.-IRMA		cpm	B/T x 100 (%)
Total count		306174	100
Calibrator	0 pg/ml	259	0.08
	8.2 pg/ml	1872	0.61
	17.0 pg/ml	4467	1.46
	44.8 pg/ml	11091	3.62
	154.0 pg/ml	39847	13.01
	674.0 pg/ml	148740	48.58

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.9 pg/ml.

B. Specificity

Some potentially interfering hormones have been tested in this assay. At concentrations up to 100 ng/ml, none of the following hormones showed significant interference :

- CGRP
- Salmon-calcitonin
- PDN 21
- Procalcitonin N terminal.

C. Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	X ± S.D. (pg/ml)	CV %	Serum	N	X ± S.D. (pg/ml)	CV %
A	10	76.1 ± 2.1	2.7	A	10	26.9 ± 2.5	9.3
B	10	337.4 ± 6.5	1.9	B	10	128.8 ± 6.0	4.6

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Added Calcitonin (pg/ml)	CT Measured Concent. (pg/ml)	Recovery (%)
0		
6	6.1	101.7
15.5	16.0	103.2
48	48.5	101.0
142.5	138.7	97.3
500	517.7	103.5

DILUTION TEST

Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)	Recovery (%)
1/1		396.0	-
1/2	198.0	193.6	97.8
1/4	99.0	92.3	93.2
1/8	49.5	49.3	99.6
1/16	24.8	23.6	95.4
1/32	12.4	12.1	97.8
1/64	6.2	6.4	103.4

E. Hook effect

A serum sample with a spiked concentration of 500.000 pg/ml of synthetic Calcitonin gives a signal above the highest calibrator concentration.

No Hook effect was observed on serum samples containing up to 16500 pg/mL of native Calcitonin.

It is recommended to dilute samples above 100 pg/mL with CT free human serum provided in the kit. Samples have to be sequentially diluted 1:10.

Healthy subjects, patients having medullary thyroid carcinoma or other tumors present different immunochemical forms of calcitonin.

The size of these forms may vary from type to type and may therefore respond differently in immunoassays and dilution tests.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make their own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than once.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

Normal values


These values are given only for guidance; each laboratory should establish its own normal range of values.

The ranges are expressed as 2.5% to 97.5% percentiles.

Identification	N	Mean (pg/ml)	SD (pg/ml)	Range (pg/ml)
Men	94	5.7	2.3	1.9 – 9.6
Women	98	3.5	1.9	0.5 – 7.8
All	192	4.6	2.4	0.9 – 9.5

Pentagastrin test:

Patients with low base levels of CT:

Patients with medullary thyroid carcinoma, and patients with other malignancies can have a non linear dilution test for calcitonin sample; this could be due to a modified form of calcitonin (polymers and/or glycosylated forms) that is not completely recognized by our antibodies.

XVII. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiology procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and Foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.

- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in nodular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) µl
Calibrators (0-5)	-	200	-
Samples	-	-	200
Tracer	50	50	50
Incubation	18 ± 1 hours at 2-8°C		
Separation	-	aspirate	
Washing solution	-	2.0 ml	
Separation	-	aspirate	
Washing solution	-	2.0 ml	
Separation	-	aspirate	
Counting	Count tubes for 60 seconds		

DIASource Catalogue Nr : KIP0429	P.I. Number : 1700541/en	Revision nr : 120914/1
-------------------------------------	-----------------------------	---------------------------

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the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.

Cancer, 51,5:856-862.

- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and Foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in nodular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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XVIII. BIBLIOGRAFIE

- GRAZE K., SPILER I.J., TASHIAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978) **Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.** Engl. J. Med., 299,18:980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974) **A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.** J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983) **the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.** Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981) **Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.** J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978) **Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.** Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985) **Parathyroid hormone, calcitonin, and the calciferols.** In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987) **SCC antigen and other tumor markers in lung cancer: preliminary results.** Excerpta Medica, 162-170.
- ELIARD, P.H. (1989) **Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.** 71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995) **Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.** Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994) **Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.** J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994) **Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.** J. Bone Miner Res. 9, 1, 53-57.

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- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and Foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in nodular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18:980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)

SCC antigen and other tumor markers in lung cancer: preliminary results.

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XVIII. BIBLIOGRAFIA

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.

- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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XVIII. BIBLIOGRAFIA

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18:980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.

- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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XVIII. REFERENSER

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18:980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.

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XVIII. BIBLIOGRAFIA

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18;980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.

Ann. Surg., 188,2:139-141.

- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and Foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in nodular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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1. GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18:980-985.
2. HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
3. ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
4. WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
5. WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
6. AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
7. BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
8. ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.
71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.
9. NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
10. PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
11. QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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XVIII. IRODALOM

- GRAZE K., SPILER I.J., TASHIJAN A.H., MELVIN K.E.W., CERVI-SKINNER S., GAGEL R.F., MILLER H.H., WOLFE H.J., DELELLIS R.A., LEAPE L., FELDMAN Z.T. and REICHLIN S. (1978)
Natural history of familial medullary thyroid carcinoma; Effect of a program for early diagnosis.
Engl. J. Med., 299,18:980-985.
- HENNESSY J.F., WELLS S.A., ONTJES D.A. and COOPER C.W. (1974)
A comparison of pentagastrin injection and calcium infusion as provocative agents for the detection of medullary carcinoma of the thyroid.
J. Clin. Endocrinol. Metab., 39:487-495.
- ROUGIER Ph., CALMETTES C., LAPLANCHE A., TRAVAGLI J.P., LEFEVRE M., PARMENTIER C., MILHAUD G. and TUBIANA M. (1983)
the values of calcitonin and carcinoembryonic antigen in the treatment and management of nonfamilial medullary thyroid carcinoma.
Cancer, 51,5:856-862.
- WALLACH S.R., ROYSTON I., TAETLE R., WOHL H. and DEFTOS L. (1981)
Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung.
J. Clin. Endocrinol. Metab., 53,3:602-606.
- WELLS S.A., BAYLIN S.B., LINEHAN W.M., FARRELL R.E., COX E.B. and COOPER C.W. (1978)
Provocative agents and the diagnosis of medullary carcinoma of the thyroid gland.
Ann. Surg., 188,2:139-141.
- AURBACH G.D., MARX S.J. and SPIEGEL A.M. (1985)
Parathyroid hormone, calcitonin, and the calciferols.
In: Williams Textbook of endocrinology (7th edition; Wilson J.D. and foster D.W. eds) W.B. Saunders Company, Philadelphia, 1137-1217.
- BODY J.J. et al. (1987)
SCC antigen and other tumor markers in lung cancer: preliminary results.
Excerpta Medica, 162-170.
- ELIARD, P.H. (1989)
Evaluation of a highly sensitive two-site immunoradiometric assay (IRMA) for human calcitonin (hCT): comparison with the RIA's for

hCT and for the carboxyl-terminal flanking peptide (PDN-21) of the hCT gene.

71th Annual meeting of the Endocrine Society, Seattle, Washington, Abst. N° 1800 p472.

- NICOLI P. et al. (1995)
Abnormal calcitonin basal levels and pentagastrin response in patients with chronic renal failure on maintenance hemodialysis.
Eur. J. Endocrinol. 132, 1, 75-81.
- PACINI F. et al. (1994)
Routine measurement of serum calcitonin in modular disease allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma.
J. Clin. Endocrinol. Metab. Excerpta Medica, 78, 4, 824-9.
- QUESADA J. M. et al. (1994)
Calcitriol corrects deficient calcitonin secretion in the Vit. D deficient elderly.
J. Bone Miner Res. 9, 1, 53-57.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACER BUF	Tracer buffer
µTIT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
CART	Cartridge
WASH SOLN	Wash buffer

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