hCG+ β -IRMA

KIP0981 - KIP0984

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P
Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246
Fax.: (763) 780-2988
Email: info@ibl-america.com
Web: www.ibl-america.com

Read entire protocol before use.

hCG+ β IRMA

I. INTENDED USE

Immunoradiometric assay for the in vitro quantitative measurement of human chorionic gonadotropin in serum and plasma.

This kit is not intended to be used for the risk evaluation of trisomy 21.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource hCG+ β IRMA kit
- B. **Catalog number :** KIP0981 : 96 tests
KIP0984 : 4 x 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.90

III. CLINICAL BACKGROUND

- A. **Biological activities**
hCG is a glycoprotein synthesised by the syncytiotrophoblast of the placenta throughout pregnancy. hCG-molecular weight 37.9 kDa - comprises two subunits. The hCG α subunit -molecular weight 14.9 Kda - is chemically similar to the α subunits of FSH, LH and TSH hormones. The hCG β subunit - molecular weight 23.0 kDa - has a structure similar to that of the LH β subunit, differing by only a few epitopes. hCG has biological characteristics similar to LH. During pregnancy, hCG stimulates the remaining corpus luteum and the placental tissue to secrete the various steroid hormones. In addition to its stimulating action on the luteal and placental tissue, hCG, by crossing the placenta, is essential to differentiate the genital tract of the fetus, which occurs around the 7th week of pregnancy.
- B. **Clinical application of hCG+ β - IRMA**
- **Diagnostic and monitoring test in pregnancy**
hCG and its free subunits α and β appear in the serum and urine of pregnant women about 9 days following ovulation. The hCG level then increases rapidly to reach a peak between the 8th and the 12th week.
 - **Tumour marker test in trophoblastic tumours**
Hydatiform moles and choriocarcinomas may secrete large amounts of native hCG and its two free subunits α and β into the peripheral blood circulation.
 - **Tumour marker test in non-trophoblastic cancers**
10 to 15 % of the breast, lung, and digestive tract cancers release hCG and/or either of its two constitutive subunits α and β .

IV. PRINCIPLES OF THE METHOD

The DIAsource hCG+β - IRMA is an immunoradiometric assay based on coated-tube separation in which several monoclonal antibodies directed against distinct epitopes of hCG have been used. Mabs 1 - the capture antibodies are attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mabs 1. Addition of Mabs 2, the signal antibody labelled with ¹²⁵I, will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

V. REAGENTS PROVIDED

Reagents	Quantity 96 tests	Quantity 4x96 tests	Colour Code	Reconstitution			
 Tubes coated with anti hCG+β (monoclonal antibodies)	2 x 48	8 x 48	grey	Ready for use			
<table border="1" data-bbox="108 667 242 712"> <tr> <td>Ab</td> <td>¹²⁵I</td> </tr> </table> Anti- hCG+β - ¹²⁵ I (monoclonal antibodies) in TRIS Buffer with bovine serum albumin, sodium azide (< 0.1 %) and inert red dye	Ab	¹²⁵ I	1 vial 5.5 ml 450 kBq	4 vials 5.5 ml 4 x 450 kBq	red	Ready for use	
Ab	¹²⁵ I						
<table border="1" data-bbox="108 869 242 913"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators 1-6 in bovine serum and thymol (see exact value on vial labels)	CAL	N	6 vials lyophil.	2 x 6 vials lyophil.	yellow	Add 0.5 ml distilled water	
CAL	N						
<table border="1" data-bbox="108 1003 242 1048"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Calibrator 0 in bovine serum and azide (0,5 %)	CAL	0	3 vials 10 ml	8 vials 10 ml	yellow	Ready for use	
CAL	0						
<table border="1" data-bbox="76 1126 290 1171"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (TRIS-HCl)	WASH	SOLN	CONC	1 vial 10 ml	3 vials 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC					
<table border="1" data-bbox="76 1238 290 1283"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls 1 and 2 in human serum with thymol	CONTROL	N	2 vials lyophil.	2 x 2 vials lyophil.	silver	Add 0.5 ml distilled water	
CONTROL	N						

Note: Also use the content of the zero calibrator for dilutions of samples. 1 mIU of the calibrator preparation is equivalent to 1 mIU of ⁴hIS 75/589

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit :

- Distilled water
- Pipettes for delivery of: 50 µl and 500 µl (the use of accurate pipettes with disposable plastic tips is recommended)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Vortex mixer
- Magnetic stirrer
- Tube shaker (400 rpm)
- Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators 1-6:** Reconstitute the calibrators with 0.5 ml distilled water.
- Controls :** Reconstitute the controls with 0.5 ml distilled water.
- Wash solution :** Prepare an adequate volume of Working Wash Solution by adding 69 volumes of distilled water to 1 volume of wash solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash Solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for 7 days at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for 3 months. Avoid subsequent freeze-thaw cycles.

- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 24h., storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- Serum and heparinized or EDTA plasma provide similar results.
 $y(\text{serum}) = 1.03x(\text{HEP. plasma}) + 0.07$ $r = 0.99$ $n = 30$
 $y(\text{serum}) = 1.00x(\text{EDTA plasma}) - 0.47$ $r = 0.99$ $n = 30$

X. PROCEDURE

A. Handling notes

- Do not use the kit or components beyond expiration date.
- Do not mix materials from different kit lots.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.
- High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- Prepare a calibrator curve for each run, do not use data from previous runs.

B. Procedure

- Label coated tubes in duplicate for each calibrator, control, specimen. For the determination of the total activity label two plain plastic tubes.
- Briefly vortex calibrators, controls, specimens and dispense 50 µl of each into the respective tubes.
- Dispense 50 µl of anti-hCG+β-¹²⁵I tracer to all tubes, including the uncoated tubes for total counts.
- Shake the rack containing the tubes gently by hand to liberate any trapped bubbles.
- Incubate for 30 min. at room temperature on a tube shaker (400 rpm) . **Respect carefully the incubation time.**
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated-tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
- Aspirate (or decant) the content of each tube (except total counts).
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of hCG+β (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used in place of the real time calibrator curve.

hCG+β IRMA		cpm	B/T (%)
Total count		293692	
Calibrator	0.0 mIU/ml	194	0.07
	5.5 mIU/ml	1127	0.38
	16.6 mIU/ml	2496	0.85
	54.0 mIU/ml	8095	2.76
	153.0 mIU/ml	25435	8.66
	307.0 mIU/ml	49646	16.90
	585.0 mIU/ml	93037	31.68

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero standards were assayed along with a set of other calibrators. The sensitivity, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 1.6 mIU/ml.

B. Specificity

Cross-reactive hormones were added to a low hCG+β value serum and to a high value calibrator (150 mIU/ml). The apparent hCG+β response was measured.

Added hormones		CAL 1 Observed hCG values	CAL 2 Observed hCG values
FSH	-	9.8 mIU/ml	284 mIU/ml
LH	250 mIU/ml	11.3 mIU/ml	292 mIU/ml
TSH	250 μIU/ml	10.9 mIU/ml	262 mIU/ml
αhCG	34000 fmol/ml	10.3 mIU/ml	290 mIU/ml
βhCG	1000 fmol/ml	11.0 mIU/ml	252 mIU/ml
		361.0 mIU/ml	699 mIU/ml

This demonstrates that the hCG+β-IRMA does not cross react with FSH, LH, TSH, αhCG and presents a 100 % cross-reaction with the free β-chain of hCG.

Note : Using the following molecular weights 37.9 kDa (hCG), 23.0 kDa (βhCG) and 14.9 kDa (αhCG), the molar equivalents are :

- hCG : 1000 fmol = 37.9 ng - 325 mIU 1st IRP 75/537
- βhCG : 1000 fmol = 23.0 ng
- αhCG : 1000 fmol = 14.9 ng

C. Precision

INTRA ASSAY

INTER ASSAY

Serum	N	<X> ± S.D. (mIU/ml)	CV (%)	Serum	N	<X> ± S.D. (mIU/ml)	CV (%)
A	20	67.3 ± 1.2	1.8	C	10	35.2 ± 2.2	6.2
B	20	245.9 ± 3.0	1.2	D	10	149.6 ± 12.6	8.4

D. Accuracy

RECOVERY TEST

Sample	Added hCG (mIU/ml)	Recovered (mIU/ml)	Recovery (%)
Serum 1	50.0	52.1	104
	100.0	101.7	104
	200.0	209.0	102
	400.0	417.1	104

DILUTION TEST

Sample	Dilution	Theoretical Concent. (mIU/ml)	Measured Concent. (mIU/ml)
Serum 1	1/1	563.6	563.6
	1/2	281.8	260.2
	1/4	140.9	126.4
	1/8	70.5	62.6
	1/16	35.2	34.1
	1/32	17.6	18.1
	1/64	8.8	8.7
	1/128	4.4	4.8

Samples were diluted with zero calibrator.

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY

	0'	10'	20'	30'
Serum 1 (mIU/ml)	31.2	31.1	30.9	31.4
Serum 2 (mIU/ml)	144.1	142.4	140.5	143.5

F. Hook effect

A Hook effect can be observed above 150000 mIU/ml. Samples from pregnant women can yield such high concentrations. Therefore, it is advisable, in cases where an accurate quantitative result of the hCG concentration should be obtained (e.g. follow up of pregnancies), to test 3 dilutions as indicated in the scheme below.

- A: undiluted
- B: 17x diluted (25 μl of A + 400 μl zero calibrator)
- C: 289x diluted (25 μl of B + 400 μl zero calibrator)

Screening for pregnancies should be tested undiluted.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

The values provided below are given only for guidance; each laboratory should establish its own normal range of values.

Pregnancy cut-off:

The cut-off for pregnancy determination, defined as the apparent concentration three standard deviations above the average counts of 97 samples from non-pregnant women, was 4.6 mIU/ml.

It is recommended to retest samples between the cut-off and 10 mIU/ml.

XVII. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. BRAUNSTEIN G.D., RASOR J., ADLER D., DANZER H., WADE H.E., (1976)
Serum chorionic gonadotropin levels throughout normal pregnancy.
Am. J. Obstet. Gynecol., 126:678.
2. CHEN F., SETSUKO G., YOSHIHITO F., YUTAKA T., (1987)
Radioimmunoassay of the serum free β -subunit of human chorionic gonadotropin in trophoblastic disease.
J. Clin. Endocrinol. Metab., 64:313.
3. DAWOOD M.Y., SASCENA B.B., LANDESMAN R., (1977)
Human chorionic gonadotropin and its subunit in hydatidiform mole and choriocarcinoma.
Obstet. Gynecol., 50:172.
4. HAY D.L., (1982)
Chorionic gonadotropin.
Clin. Biochem., 3:35.
5. GASPARD V.J., REUTER A.M., DEVILLE J.L., VRINDTS-GEVAERT Y., BAGSHAW K.D., FRANCHIMONT P., (1980)
Serum concentration of human chorionic gonadotropin and its alpha and beta subunit II Trophoblastic tumours.
Clin. Endocrinol. (Oxf.), 13:219.
6. PIERCE J.G., PARSONS T.F., (1981)
Glycoprotein hormones : structure and function.
Annu. Rev. Biochem., 50:465.

7. MANCINI G. et al., (1992)
hCG, AFP and V E₃ pattern in the 14-20th weeks of Down's syndrome pregnancies.
Prenatal Diagnost., 12(7):619-24.
8. WHITEHEAD N. et al., (1993)
Elevated material serum human chorionic gonadotropin increases the chance of adverse pregnancy outcome.
Am. J. Obstet. Gynecol., 169(5):1359-60.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) ml
Calibrators (0-5)	-	0.05	-
Samples	-	-	0.05
Tracer	0.05	0.05	0.05
Incubation	30 minutes at room temperature with shaking at 400 rpm		
Separation	-	aspirate (or decant)	
Washing solution	-	2.0	
Separation	-	aspirate (or decant)	
Washing solution	-	2.0	
Separation	-	aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP0981 - KIP0984	P.I. Number : 1700467/en	Date of issue : 110218/1
---	-----------------------------	-----------------------------

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Data wydania : 2011-02-18

	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µP	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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