# **hGH-IRMA**

***KIP1081***

**Distributed By:**

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: [info@ibl-america.com](mailto:info@ibl-america.com)

Web: [www.ibl-america.com](http://www.ibl-america.com)

---

**LOT** : 110218/2

en

Read entire protocol before use.

## hGH-IRMA

### I. INTENDED USE

Immunoradiometric assay kit for the *in vitro* quantitative measurement of human growth hormone (hGH) in serum.

### II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource hGH-IRMA
- B. **Catalogue number :** KIP1081 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.  
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

**For technical assistance or ordering information contact :**  
**Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91**

### III. CLINICAL BACKGROUND

#### A. Biological activities of hGH

hGH is a polypeptide hormone (molecular weight 21,500 Da) produced by the acidophil cells of the anterior pituitary under the control of two main substances from the median eminence : Growth-hormone Releasing Factor (GRF) and an inhibitory agent, somatostatin. Dopaminergic, adrenergic and serotonergic neuroendocrine pathways also play an important role in the control of hGH secretion. Excitatory stimuli of hGH secretion include hypoglycemia, exercise, fasting, meals with a high protein content, deep sleep, stress, glucagon, L Dopa, amino acids, etc. Inhibitory stimuli include glucose, cortisol, hGH and free fatty acids. Because of its short plasma half life ( $\pm 25$  minutes) and of the frequent excitatory or inhibitory stimuli, hGH displays frequent and large variations of concentration in serum.

One of the main physiological functions of hGH is to act on the liver and other tissues to produce Somatomedins, which in turn induce growth by direct action on target tissues. In contrast to hGH, the concentration of somatomedin in serum is kept stable by virtue of being largely bound to circulating plasma proteins.

#### B. Clinical application of hGH-IRMA

· **Growth retardation**

hGH hyposecretion is one of the various causes of small stature in children. Serum hGH measurement with a highly sensitive assay, especially following a provocative stimulus (absence of response), is an important way to establish this diagnosis because this group of patients can be treated by administration of hGH.

· **Hypopituitarism**

Serum hGH measurement is also an index of pituitary function when hypopituitarism (either idiopathic or due to tumour and surgery) is suspected.

· **Gigantism and acromegaly**

Serum hGH measurement, especially following a provocative inhibitory test (absence of response), is an important way to establish the diagnosis of hGH hypersecretion due to acidophilic pituitary tumour. This results in gigantism in children and acromegaly in adults. Both of these disorders may be treated by surgery or radiation.

#### IV. PRINCIPLES OF THE METHOD

The DIASource hGH-IRMA is an immunoradiometric assay based on coated-tube separation. Mab 1 - the capture antibody is attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mab 1. Addition of Mab 2, the signal antibody labeled with  $^{125}\text{I}$ , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

#### V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with anti hGH (monoclonal antibodies)	2 x 48	white	Ready for use			
<table border="1" data-bbox="119 616 255 660"> <tr> <td>Ab</td> <td><math>^{125}\text{I}</math></td> </tr> </table> TRACER: $^{125}\text{I}$ Iodine labeled anti-hGH (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide (< 0.1%) and an inert red dye.	Ab	$^{125}\text{I}$	1 vial 5.5 ml 350 kBq	red	Ready for use	
Ab	$^{125}\text{I}$					
<table border="1" data-bbox="119 806 255 851"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero calibrator (0 $\mu\text{IU/ml}$ ) in human serum with thymol	CAL	0	1 vial lyophilized	yellow	Add 2 ml distilled water	
CAL	0					
<table border="1" data-bbox="119 907 255 952"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators hGH N = 1 to 6 (see exact values on vial labels) in human serum with thymol	CAL	N	6 vials lyophilized	yellow	Add 0.5 ml distilled water	
CAL	N					
<table border="1" data-bbox="71 1052 319 1097"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (NaCl, Tween 20)	WASH	SOLN	CONC	1 vial 40 ml	green	Dilute 20 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="71 1153 279 1198"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum with thymol	CONTROL	N	2 vials lyophilized	silver	Add 0.5 ml distilled water	
CONTROL	N					

**Note :** Use the zero calibrator for sample dilutions.  
1  $\mu\text{IU}$  of the calibrator preparation is equivalent to 1  $\mu\text{IU}$  NIBSC 2<sup>nd</sup> IS 98/574.

#### VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 50  $\mu\text{l}$ , 500  $\mu\text{l}$  and 2 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring  $^{125}\text{I}$  may be used (minimal yield 70%).

#### VII. REAGENT PREPARATION

- Calibrators :** Reconstitute the zero calibrator with 2.0 ml distilled water and other calibrators with 0.5 ml distilled water.
- Controls :** Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 19 volumes of distilled water to 1 volume of Wash Solution (20x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

#### VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for 3 months. Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

#### IX. SPECIMEN COLLECTION AND PREPARATION

- Serum samples must be kept at 2-8°C.
- If the test is not run within 24 hrs., storage at -20°C is recommended.
- Avoid successive freezing and thawing.

#### X. PROCEDURE

##### A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

##### B. Procedure

- Label coated tubes in duplicate for each calibrator, sample, and control. For the determination of total counts, label 2 normal tubes
- Briefly vortex calibrators, controls and samples and dispense 50 $\mu\text{l}$  of each into respective tubes.
- Dispense 50  $\mu\text{l}$  of  $^{125}\text{I}$  Iodine labeled anti-hGH into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand.
- Incubate for 2 hours at room temperature.
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 3 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
- Aspirate (or decant) the content of each tube (except total counts).
- Wash tubes again with 3 ml Wash solution (except total counts) and aspirate (or decant).
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

#### XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- On semilogarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of hGH (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

## XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

hGH-IRMA		cpm	B/T (%)
Total count		126065	
Calibrator	0 $\mu$ IU/ml	241	0.19
	1.0 $\mu$ IU/ml	1254	0.99
	2.4 $\mu$ IU/ml	2743	2.18
	8.0 $\mu$ IU/ml	7817	6.20
	25.5 $\mu$ IU/ml	19515	15.48
	59.1 $\mu$ IU/ml	34310	27.22
	119.9 $\mu$ IU/ml	49965	39.63

## XIII. PERFORMANCE AND LIMITATIONS

### A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.04  $\mu$ IU/ml.

### B. Specificity

Cross-reactive hormones were added to a low and to a high hGH value calibrator. The apparent hGH response was measured.

added Hormone	hGH CAL 1 $\mu$ IU/ml	hGH CAL 5 $\mu$ IU/ml
-	0.60	65
hCG 300000 mIU/ml	0.63	63
hPL 50000 ng/ml	0.69	68
PRL 2000 ng/ml	1.35	63

**Note :** this table shows the cross-reactivity for the hGH-IRMA.

### C. Precision

INTRA-ASSAY				INTER-ASSAY			
Serum	N	<X> $\pm$ SD ( $\mu$ IU/ml)	CV (%)	Serum	N	<X> $\pm$ SD ( $\mu$ IU/ml)	CV (%)
A	6	6.7 $\pm$ 0.3	4.4%	A	20	3.8 $\pm$ 0.3	8.1%
B	6	19.0 $\pm$ 0.2	1.0%	B	20	17.9 $\pm$ 1.2	6.7%

SD : Standard Deviation; CV: Coefficient of variation

### D. Accuracy

DILUTION TEST			
Sample	Dilution	Theoretical Concent. ( $\mu$ IU/ml)	Measured Concent. ( $\mu$ IU/ml)
Serum 1	1/1	-	35.6
	1/2	17.8	17.4
	1/4	8.9	8.8
	1/8	4.5	5.0
	1/16	2.2	2.4
	1/32	1.1	1.2

Sample was diluted with the zero calibrator.

### Conversion factor:

1 $\mu$ IU hGH-IRMA Calibrator = 0.33 ng

## RECOVERY TEST

Sample	added hGH ( $\mu$ IU/ml)	Recovered hGH ( $\mu$ IU/ml)	Recovered (%)
Serum 1	5	5.1	102%
	10	9.5	95%
	25	24.7	99%
	50	53.7	107%
Serum 2	60	63	105%
	40	46.8	117%
	20	22.1	111%
	10	10.2	102%
	5	4.8	96%

### E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to coated tubes.

## TIME DELAY

Serum ( $\mu$ IU/ml)	0'	15'	30'
Serum 1	4.4	4.5	4.8
Serum 2	18.8	18.1	19.3
Serum 3	27.6	26.2	27.0
Serum 4	36.5	36.8	36.6

### F. Hook effect

A sample spiked with hGH up to 4000  $\mu$ IU/ml gives a signal above the highest calibrator concentration.

## XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

## XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices.

## XVI. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

In normal subjects, growth hormone (hGH) secretion is pulsatile. During day-time, hGH concentrations range from <0.2-10  $\mu$ IU/ml. During sleep, hGH concentrations increase consistently ( $\pm$  30  $\mu$ IU/ml).

hGH secretion is greatly stimulated by exercise and stress (venous puncture, hypoglycemia, ...), but is decreased by hyperglycemia.

In normal subjects (n=34), two hours after oral glucose load (75 g in adults), hGH levels were lower than 10  $\mu$ IU/ml and the hGH response to stimulation tests (insulin, arginine, glucagon administration) exceeded 20  $\mu$ IU/ml.

hGH levels were elevated (even after glucose load) in acromegaly (> 10  $\mu$ IU/ml). In hGH deficiency, the response to stimulation test is absent or blunted (short stature by hGH deficiency; hypopituitarism from various origins).

**XVII. PRECAUTIONS AND WARNINGS**

**Safety**

For in vitro diagnostic use only.

This kit contains <sup>125</sup>I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents or serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves. All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosity procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

**XVIII. BIBLIOGRAPHY**

1. CAMANNI F., MASSARA F., BELFORTE L., ROSATELLO A., MOLENATTI G.M. (1977)  
**Effect of dopamine on plasma growth hormone and prolactin levels in normal and acromegalic subjects.**  
J. Clin. Endo. and Metab., 44:465.
2. DAUGHADAY W.H. (1981)  
**The adenohipophysis**  
In : "Textbook of Endocrinology", R.H. Williams ed., W.B. Saunders, Philadelphia, p:73.
3. DROBNY E.C., AMBURN K., BAUMANN G. (1983)  
**Circadian variation of basal plasma growth hormone in man.**  
J. Clin. Endo. and Metab., 57:524-528.
4. HASHIDA S., NAKAWAGA K., ISHIKAWA E., OHTAKI S. (1985)  
**Basal level of human growth hormone (hGH) in normal serum.**  
Clin. Chim. Acta, 151:185-186.
5. LEWIS V.J. et al. (1980)  
**Human growth hormone : A complex of proteins.**  
Rec. Progr. Horm. Res., 36:488.

6. LI C.H. (1974)  
**Human growth hormone : perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
7. REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
8. CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
9. KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

**XIX. SUMMARY OF THE PROTOCOL**

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Calibrators (0 to 6)	-	0.05	-
Samples, controls	-	-	0.05
Tracer	0.05	0.05	0.05
Incubation	2 hours at room temperature		
Separation	-	aspirate (or decant)	
Working Wash solution	-	3.0	
Separation	-	aspirate (or decant)	
Working Wash solution	-	3.0	
Separation	-	aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP1081	P.I. Number : 1700469/en	Revision nr : 110218/1
-------------------------------------	-----------------------------	---------------------------

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- LEWIS V.J. et al. (1980)  
**Human growth hormone : A complex of proteins.**  
Rec. Progr. Horm. Res., 36:488.
- LI C.H. (1974)  
**Human growth hormone : perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
- REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
- CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
- KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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- LI C.H. (1974)  
**Human growth hormone : perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
- REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
- CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
- KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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5. LEWIS V.J. et al. (1980)  
**Human growth hormone : A complex of proteins.**  
Rec. Progr. Horm. Res., 36:488.
6. LI C.H. (1974)  
**Human growth hormone : perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
7. REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
8. CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
9. KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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- LI C.H. (1974)  
**Human growth hormone: perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
- REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
- CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
- KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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- HASHIDA S., NAKAWAGA K., ISHIKAWA E., OHTAKI S. (1985)  
**Basal level of human growth hormone (hGH) in normal serum.**  
Clin. Chim. Acta, 151:185-186.
- LEWIS V.J. et al. (1980)  
**Human growth hormone : A complex of proteins.**  
Rec. Progr. Horm. Res., 36:488.
- LI C.H. (1974)  
**Human growth hormone : perspective on its chemistry and physiology.**  
In : "Advances in human growth hormone research", S. Raiti ed. U.S. Department of Health, Education and Welfare Publ. DMEW publ. no. (NJH) 74-612, p. 321.
- REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
- CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
- KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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- REUTENS AT. (1995)  
**Evaluation and application of a highly sensitive assay for serum growth hormone (GH) in the study of adult GH deficiency.**  
J. Clin. Endocrinol. Metab., 80(2): 480-485
- CROWLEY S. et al. (1995)  
**Growth and the growth hormone axis in prepubertal children with asthma.**  
J. Pediatr. 126(2): 297-303.
- KELLY JT. et al. (1993)  
**Effect of different oral oestrogen formulations on insulin-like growth factor I, growth hormone and growth hormone binding protein in post-menopausal women.**  
Clin. Endocrinol. 39(5) : 561-567

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	<b>Used symbols</b>
	Consult instructions for use
	Storage temperature
	Use by
<b>LOT</b>	Batch code
<b>REF</b>	Catalogue number
<b>CONTROL</b>	Control
<b>IVD</b>	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
<b>µP</b>	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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