Urinary hGH-IRMA

KIP1131

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

LOT : 121023/1

Read entire protocol before use.

Urinary hGH-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the *in vitro* quantitative measurement of human Growth Hormone (hGH) in urine.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource Urinary hGH-IRMA Kit
- B. **Catalog number :** KIP1131: 25 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological Activity

hGH is a polypeptide hormone (molecular weight 21,500 Da) produced by the acidophil cells of the anterior pituitary under the control of two main substances from the median eminence : Growth-hormone Releasing Factor (GRF) and an inhibitory agent, somatostatin. Dopaminergic, adrenergic and serotonergic neuroendocrine pathways also play an important role in the control of hGH secretion. Excitatory stimuli of hGH secretion include hypoglycemia, exercise, fasting, meals with a high protein content, deep sleep, stress, glucagon, L Dopa, amino acids, etc. Inhibitory stimuli include glucose, cortisol, hGH and free fatty acids.

Because of its short plasma half life (± 25 minutes) and of the frequent excitatory or inhibitory stimuli, hGH displays frequent and large variations of concentration in serum. One of the main physiological functions of hGH is to act on the liver and other tissues to produce somatomedins, which in turn induce growth by direct action on target tissues. In contrast to hGH, the concentration of somatomedin in serum is kept stable by virtue of being largely bound to circulating plasma proteins.

B. Clinical Application

It is known that serum growth hormone (GH) concentrations can fluctuate greatly due to large day-to-day variability of GH secretion. The most predictable GH peaks occurs about one hour after onset of sleep. GH release during the day is very pulsative: for this reason a simple blood collection provides no reliable information on GH secretion. It has been suggested that sensitive urinary Growth Hormone determination, based on 24 hours urine collection might reflect central GH secretion and therefore can help in the reliable measurement of GH secretion which is necessary in several pathological settings as. Overnight collection reflects the major episodes of secretion of growth hormone during slow wave sleep and avoids the effect of physiological during the day. There is a good correlation between the u-HGH values and the integrated serum HGH values during a 12-h or 24-h profile (ref 5,7) and between the serum and u-HGH values during provocation tests (ref 10,11) . Also when serum values are high as in gigantism or acromegaly u-HGH values are high too (ref 10,15) The concentration of HGH in urine is only about 0.01% of that in serum but such levels can now be measured with the use of monoclonal antibodies bound to a solid phase as demonstrated first by Girard et al. in Basle (ref 6) and it seems that these very low concentrations do, in fact, have physiological significance.

- Growth retardation

hGH hyposecretion is one of the various causes of small stature in children.

- Gigantism and acromegaly

hGH hypersecretion due to acidophilic pituitary tumour. This results in gigantism in children and acromegaly in adults. Both of these disorders may be treated by surgery or radiation.

IV. PRINCIPLES OF THE METHOD

The Urinary-hGH-IRMA is an immunoradiometric assay based on coated-bead separation. Mab1, the capture antibodies, are attached to the plastic beads. Calibrators, controls and samples are first incubated with the beads to extract the hGH. After incubation, washing removes the occasionally excess of antigen. Addition of Mab2, the signal antibody labelled with ¹²⁵I, will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the beads reflects the antigen concentration.

V. REAGENTS PROVIDED

Reagents	Quantity 25 tests	Colour Code	Reconstitution
BEAD Beads coated with anti hGH (monoclonal antibodies)	100	white	Ready for use
Ab 125I CONC Anti-hGH- ¹²⁵ I (monoclonal antibodies) in Phosphate Buffer with bovine albumin and sodium azide (0.05 %)	1 vial 0.6 ml 380 kBq	red	Dilute 51x with Tracer Buffer (see section VII. C)
TRACER BUF Tracer Buffer: Borate Buffer with bovine casein, bovine serum, EDTA and azide (0.05%)	1 vial 26 ml	black	Ready for use
CAL N Calibrators 1-5 in human serum with thymol (see exact values on vial labels)	5 vials lyophil.	yellow	Add 1 ml distilled water
DIL CAL Calibrator Diluent in Phosphate Buffer with bovine serum and sodium azide (< 0.1%)	1 vial 100 ml	blue	Ready for use
DIL BUF Diluent Buffer: Phosphate Buffer Saline with azide (0.05%)	1 vial powder	white	Dissolve in 300 ml distilled water (use a magnetic stirrer).
WASH SOLN CONC Wash solution (TRIS-HCl)	1 vial 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
CONTROL N Controls 1 and 2 in Phosphate buffer with bovine albumin and thymol	2 vials lyophil.	silver	Add 11 ml Calibrator Diluent
Glass vials (26x45 mm) with caps for extraction	25 vials	white	Ready for extraction

Note: 1 pg of the calibrator preparation is equivalent to 1 pg of NIBSC 2nd IS98/574

VI. SUPPLIES NOT PROVIDED

The following material is required and available at DIAsource :

1. Vial shaker (IKA Vibrax 300 rpm) N° art. 4300004.
2. Vial extraction insert N° art. 4300606 to use in combination with the vial shaker

The following material is required but not provided in the kit :

1. Distilled water
2. Pipettes for delivery of: 100 µl, 200 µl, 1 ml and 10 ml. (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Plastic tube (12x75)
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration system

8. Any gamma counter capable of measuring ¹²⁵I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. Calibrators** : Reconstitute the calibrators with 1 ml distilled water.
- B. Controls** : Reconstitute the controls with 11 ml Calibrator Diluent.
- C. Tracer** : Prepare an adequate volume of Tracer solution by adding 20 µl of Anti-hGH-¹²⁵I to 1 ml of Tracer Buffer. Use a vortex to homogenize.
- D. Diluent Buffer** : Transfer quantitatively the Diluent Buffer powder into 300 ml distilled water Use a magnetic stirrer to homogenize..
- E. Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C; except the Diluent Buffer, it must be stored at room temperature.
- After reconstitution, the calibrators are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, Tracer solution is stable until expiry date, if kept in the original well-closed vial at -20°C.
- After reconstitution, the Diluent Buffer is stable for one month at room temperature.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Urine samples must be kept at 2-8°C.
- If the test is not run within 24 hours, storage in aliquots at -20°C is recommended.
- Avoid successive freezing and thawing

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

I. Extraction step :

1. Label glass vial (26x45 mm) in single for extraction : 1 NSB, 5 calibrators, 2 controls and up to 17 samples.
2. Add 10 ml of Calibrator Diluent into the NSB vial.
3. Add 10 ml of Calibrator Diluent and 100 µl of each calibrator into the respective vials
4. Add 10 ml control or sample in the respective vial.
5. Dispense 10 ml of Diluent Buffer in each vial (NSB, calibrators, controls and samples).
6. Put four beads in each vial.
7. Close the vials with a cap and incubate overnight on a vial shaker at 300 rpm (supplied by DIAsource under the art. 4300004). The foam including the extraction vials can be directly inserted into the vial extraction insert (available from DIAsource under the art. 4300606) and used on the vial shaker.
8. Aspirate the content of each vial. Be sure that the plastic tip of the aspirator reaches the bottom of the vial in order to remove all the liquid.
9. Wash the vials with 6 ml Working Wash solution.
10. Aspirate the content of each vial.
11. Wipe gently the beads on absorbing paper and transfer the beads in the corresponding plastic tubes.

II. Incubation step :

12. Label plastic tubes (12 x 75 mm) in duplicate for each calibrator. For the determination of total count, label 2 tubes.
13. Put two beads in each tube for NSB, calibrators, controls and samples.
14. Dispense 500 µl of the Tracer solution in each tube (calibrators, controls and samples) including the tubes for total counts.
15. Shake the tube rack gently by hand to liberate any trapped air bubbles.
16. Incubation 18 to 24 hours at room temperature.

17. Aspirate the content of each tube. Be sure that the plastic tip of the aspirator reaches the bottom of the tube in order to remove all the liquid.
18. Wash tube with 2 ml Working Wash solution (except total counts) and aspirate. Avoid foaming during the addition of Working Wash solution.
19. Aspirate the content of each tube (except total counts).
20. Wash tubes again with 2 ml Working Wash solution (except total counts) and aspirate.
21. After the washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
22. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. On semi-logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of Urinary hGH (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
3. Read the concentration for each control and sample by interpolation on the calibration curve.
4. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

Urinary hGH-IRMA		cpm	B/T (%)
Total count		273476	100
NSB		549	0.2
Calibrator	1.65 ng/L	1735	0.6
	4.95 ng/L	3350	1.2
	16.00 ng/L	10924	3.9
	49.50 ng/L	26732	9.5
	215.00 ng/L	89325	31.7

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations above the average counts at zero binding, was 0.42 ng/L.

B. Specificity

Some potentially interfering hormones have been tested in this assay. None of the following hormones showed significant interference :

- hCG 300000 mIU/ml
- hPL 50000 ng/ml
- PRL 2000 ng/ml

C. Precision

INTRA ASSAY

INTER ASSAY

Urine	N	<X> ± SD (ng/L)	CV (%)	Urine	N	<X> ± SD (ng/L)	CV (%)
A	20	9.79 ± 0.86	8.7	D	10	9.36 ± 0.81	8.6
B	20	6.15 ± 0.54	8.8	E	10	103.27 ± 7.49	7.3

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST

Sample	Added hGH (ng/L)	Recovered hGH (ng/L)	Recovery (%)
1	50.00	46.35	92.5
	16.00	14.87	92.9
	4.95	4.58	92.7

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/L)	Measured Concent. (ng/L)
1	1/1	-	135.03
	1/2	67.52	68.78
	1/5	27.01	25.69

Samples were diluted with calibrator diluent.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practices

XVI. REFERENCE INTERVALS

The values provided below are given only for guidance; each laboratory should establish its own normal range of values.

Normal subjects	Mean (ng/24h)	Range (ng/24h)	N
All subjects	1.91	0.42 – 4.61	50

Remark : the range is based on 2.5 % and 97.5 % percentiles.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. BUTLER J.
Biochemical tests of growth hormone status in short children.
Ann. Clin. Biochem. 2001; 38:1-2.
2. EVANS A J.
Screening tests for growth hormone deficiency.
Journal of the Royal society of medicine. 1995; 88:161-165.
3. BERCU BB, SHULMAN D, ROOT AW, SPILLOTIS BE.
Growth hormone (GH) provocative testing frequently does not reflect endogenous GH secretion.
J. Clin. Endocrinol. Metab. 1986; 63:709-716.
4. ZADIK Z, CHALEW MD, GILULA Z, KOWARSKI AA.
Reproducibility of growth hormone testing procedures : a comparison between 24-hour integrated concentration and pharmacological stimulation.
J. Clin. Endocrinol. Metab. 1990; 71:1127-1130.
5. LARCEY KA, HEWISON A, PARKIN JM
Exercise as a screening test for growth hormone deficiency in children.
Arch. Dis. Child. 1973; 48:508-512.
6. GIRARD J, ERB T, PAMPALONE A, EBERLE AN, BAUMANN JB.
Growth hormone in urine : development of an ultrasensitive assay applicable to plasma and urine.
Horm. Res.. 1987; 28:71-80.
7. ALBINI CH, SOTOS J, SHERMAN B, JOHNSON, et al.
Diagnostic significance of urinary growth hormone measurements in children with growth failure : correlation between serum and urine growth hormone.
Pediatr. Res. 1991; 29:619-622.
8. EDGE JA, HOURD P, EDWARDS R, DUNGER DB.
Urinary growth hormone during puberty in normal and diabetic children.
Clin. Endocrinol. 1989; 30:413-420.
9. OKUNO A, YANO K, ITOH T, et al.
Urine growth hormone determinations compared with other methods in the assessment of growth hormone secretion.
Acta. Paed. Scand. 1987; 337(suppl):74-81.
10. SUKEGAWA I, HIZUKA N, TAKANO K, et al.
Measurement of nocturnal urinary growth hormone values.
Acta. Endocrinol. 1989; 121:290-296.
11. EVANS AJ, WOOD PJ.
Development of an assay for human growth hormone in urine using commercially available reagents.
Ann. Clin. Biochem. 1989; 26:353-357.
12. HOURD P, EDWARDS R.
Measurement of human growth hormone in urine : development and validation of a sensitive and specific assay.
J. Endocrinol. 1989; 121:167-175.
13. HASHIDA S, ISHIKAWA E, KATO T, et al.
Human growth hormone (hGH) in urine and its correlation to serum hGH examined by a highly sensitive sandwich enzyme immunoassay.
Clin. Chim. Acta. 1987; 162:229-235.
14. GIRARD J, FISHER-WASELS T.
Measurement of urinary growth hormone.
Horm. Res. 1990; 33(suppl4):12-18.
15. EVANS AJ, WILLIS DS, WOOD PJ.
The assay of urinary growth hormone in normal and acromegalic adults.
Clin. Endocrinol. 1991; 35:413-418.
16. WALKER JM, WILLIAMSON S, BETTS PR, EVANS AJ, WOOD PJ.
Urinary growth hormone excretion as a screening test for human growth deficiency.
Arch. Dis. Child. 1990; 65:89-92.
17. TOMITA H, OGAWA M, KOMIJO T, et al.
A highly sensitive sandwich enzyme immunoassay of urinary growth hormone in children with short stature, Turner's syndrome and simple obesity.
Acta. Endocrinol. 1989; 121:513-519.

18. HASHIDA S, ISHIKAWA E, NAKAGAWA K, et al.
Level of human growth hormone (hGH) in urine determined by a highly specific and sensitive sandwich enzyme immunoassay.
Anal. Leut. 1986; 19:625-638.
19. CROWNE EC, WALLACE WHB, SHALET SM, ADDISON GM, PRICE DA.
Relationship between urinary and serum growth hormone and pubertal status.
Arch. Dis. Child. 1991; 67:91-95.
20. PRICE DA, ADDISON GM, HERBERT ED.
Increase in urinary growth hormone excretion in puberty.
Arch. Dis. Child. 1990; 65:1203-1204.
21. HATTORI N, KATO Y, MURAKAMI Y, et al.
Urinary growth hormone levels measured by ultrasensitive enzyme immunoassay in patients with renal insufficiency.
J. Clin. Endocrinol. Metab. 1988; 66:727-732.
22. TASSONI P, CACCIARI E, CAU M, et al.
Variability of growth hormone response to pharmacological and sleep tests performed twice in short children.
J. Clin. Endocrinol. Metab. 1990; 71:230-234.
23. FORNITO MC, CALOGERO AE, MANGIOI A, et al.
Intra- and inter-individual variability in growth hormone responses to growth hormone releasing 164P hormone.
J. Neuroendocrinol. 1990; 2:87-90.
24. MAIN K, JARDEN M, LOTTE A, DINESEN B, HERTEL T, JUUL A, MIJLLER J, SKAKKEBAEK N.
The impact of gender and puberty on reference values for urinary growth hormone excretion : a study of 3 morning urine samples in 517 children and adults
J. Clin. Endocrinol. Metab.. 1994; 79:865-871.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	NSB ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Extraction				
Calibrator Diluent	-	10 ml	10 ml	-
Calibrators (1-5)	-	-	100 µl	-
Controls, samples	-	-	-	10 ml
Diluent Buffer	-	10 ml	10 ml	10 ml
Beads	-	4 beads	4 beads	4 beads
Incubation	Overnight at room temperature with shaking at 300 rpm			
Separation	-		aspirate	
Washing solution	-		6.0 ml	
Separation	-		aspirate	
Incubation				
Extracted NSB	-	2 beads	2 beads	2 beads
Calibrators, controls and samples	0.5	0.5	0.5	0.5
Tracer				
Incubation	18 - 24 hours at room temperature			
Separation	-		aspirate	
Washing solution	-		2.0 ml	
Separation	-		aspirate	
Washing solution	-		2.0 ml	
Separation	-		aspirate	
Counting	Count tubes for 60 seconds			

DIAsource Catalogue Nr : KIP1131	P.I. Number : 1700440/en	Date of issue : 121023/1
-------------------------------------	-----------------------------	-----------------------------

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XVIII. BIBLIOGRAPHIE

1. BUTLER J. **Biochemical tests of growth hormone status in short children.** Ann. Clin. Biochem. 2001; 38:1-2.
2. EVANS A J. **Screening tests for growth hormone deficiency.** Journal of the Royal society of medicine. 1995; 88:161-165.
3. BERCU BB, SHULMAN D, ROOT AW, SPILLOTIS BE. **Growth hormone (GH) provocative testing frequently does not reflect endogenous GH secretion.** J. Clin. Endocrinol. Metab. 1986; 63:709-716.
4. ZADIK Z, CHALEW MD, GILULA Z, KOWARSKI AA. **Reproducibility of growth hormone testing procedures : a comparaison between 24-hour integrated concentration and pharmacological stimulation.** J. Clin. Endocrinol. Metab. 1990; 71:1127-1130.
5. LARCEY KA, HEWISON A, PARKIN JM **Exercise as a screening test for growth hormone deficiency in children.** Arch. Dis. Child. 1973; 48:508-512.
6. GIRARD J, ERB T, PAMPALONE A, EBERLE AN, BAUMANN JB. **Growth hormone in urine : development of an ultrasensitive assay applicable to plasma and urine.** Horm. Res.. 1987; 28:71-80.
7. ALBINI CH, SOTOS J, SHERMAN B, JOHNSON, et al. **Diagnostic significance of urinary growth hormone measurements in children with growth failure : correlation between serum and urine growth hormone.** Pediatr. Res. 1991; 29:619-622.
8. EDGE JA, HOURD P, EDWARDS R, DUNGER DB. **Urinary growth hormone during puberty in normal and diabetic children.** Clin. Endocrinol. 1989; 30:413-420.
9. OKUNO A, YANO K, ITOH T, et al. **Urine growth hormone determinations compared with other methods in the assessment of growth hormone secretion.** Acta. Paed. Scand. 1987; 337(suppl):74-81.
10. SUKEGAWA I, HIZUKA N, TAKANO K, et al. **Measurement of nocturnal urinary growth hormone values.** Acta. Endocrinol. 1989; 121:290-296.
11. EVANS AJ, WOOD PJ. **Development of an assay for human growth hormone in urine using commercially available reagents.** Ann. Clin. Biochem. 1989; 26:353-357.
12. HOURD P, EDWARDS R. **Measurement of human growth hormone in urine : development and validation of a sensitive and specific assay.** J. Endocrinol. 1989; 121:167-175.
13. HASHIDA S, ISHIKAWA E, KATO T, et al. **Human growth hormone (hGH) in urine and its correlation to serum hGH examined by a highly sensitive sandwich enzyme immunoassay.** Clin. Chim. Acta. 1987; 162:229-235.
14. GIRARD J, FISHER-WASELS T. **Measurement of urinary growth hormone.** Horm. Res. 1990; 33(suppl4):12-18.
15. EVANS AJ, WILLIS DS, WOOD PJ. **The assay of urinary growth hormone in normal and acromegalic adults.** Clin. Endocrinol. 1991; 35:413-418.
16. WALKER JM, WILLIAMSON S, BETTS PR, EVANS AJ, WOOD PJ. **Urinary growth hormone excretion as a screening test for human growth deficiency.** Arch. Dis. Child. 1990; 65:89-92.
17. TOMITA H, OGAWA M, KOMIJO T, et al. **A highly sensitive sandwich enzyme immunoassay of urinary growth hormone in children with short stature, Turner's syndrome and simple obesity.** Acta. Endocrinol. 1989; 121:513-519.
18. HASHIDA S, ISHIKAWA E, NAKAGAWA K, et al. **Level of human growth hormone (hGH) in urine determined by a highly specific and sensitive sandwich enzyme immunoassay.** Anal. Leut. 1986; 19:625-638.
19. CROWNE EC, WALLACE WHB, SHALET SM, ADDISON GM, PRICE DA. **Relationship between urinary and serum growth hormone and pubertal status.** Arch. Dis. Child. 1991; 67:91-95.
20. PRICE DA, ADDISON GM, HERBERT ED. **Increase in urinary growth hormone excretion in puberty.** Arch. Dis. Child. 1990; 65:1203-1204.
21. HATTORI N, KATO Y, MURAKAMI Y, et al. **Urinary growth hormone levels measured by ultrasensitive enzyme immunoassay in patients with renal insufficiency.** J. Clin. Endocrinol. Metab. 1988; 66:727-732.
22. TASSONI P, CACCIARI E, CAU M, et al. **Variability of growth hormone response to pharmacological and sleep tests performed twice in short children.** J. Clin. Endocrinol. Metab. 1990; 71:230-234.
23. FORNITO MC, CALOGERO AE, MANGIOI A, et al. **Intra- and inter-individual variability in growth hormone responses to growth hormone releasing 164P hormone.** J. Neuroendocrinol. 1990; 2:87-90.
24. MAIN K, JARDEN M, LOTTE A, DINESEN B, HERTEL T, JUUL A, MIJLLER J, SKAKKEBAEK N. **The impact of gender and puberty on reference values for urinary growth hormone excretion : a study of 3 morning urine samples in 517 children and adults.** J. Clin. Endocrinol. Metab.. 1994; 79:865-871.

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Date de révision: 2012-10-23

	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
CART	Cartridge
WASH SOLN	Wash buffer

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