17OH-RIA-CT

KIP1409

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P
Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246
Fax.: (763) 780-2988
Email: info@ibl-america.com
Web: www.ibl-america.com

LOT : 110217/2

Read entire protocol before use.

17OH-RIA-CT

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 17- α -hydroxyprogesterone (17OH) in serum and plasma.

II. GENERAL INFORMATION

- A. Proprietary name :** DIAsource 17OH-RIA-CT Kit
- B. Catalog number :** KIP1409 : 96 tests
- C. Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :

Tel: +32 (0) 10 84.99.11

Fax: +32 (0) 10 84.99.91

III. CLINICAL BACKGROUND

A Biological activities of 17- α -hydroxyprogesterone

17- α -hydroxyprogesterone (17OH) is a C-21 steroid hormone (molecular weight 330.3) which is produced from 17- α -hydroxy pregnenolone in the adrenals and also in the ovaries, testes and placenta. 17OH is hydroxylated at 11 and 21 positions to produce cortisol via 11-deoxycortisol.

B. Clinical applications of 17- α -hydroxyprogesterone determination

As a rule, serum or amniotic fluid 17OH dosages are relevant to diagnose congenital adrenal hyperplasia (CAH). This CAH is due to a specific enzyme defect (six distinct enzyme deficiencies have been described). As a result of these deficiencies, ACTH increases and produces adrenal hyperplasia and the raise of many steroid precursors. But it is also very interesting to know the value of 17OH in patients with varicocele. (17OH and testosterone represent markers of Leydig cell function) and in ageing male patients to detect Benign Prostatic Hypertrophy (BPH) and carcinoma of the prostate (PCA) (plasma 17OH is significantly lower in PCA and BPH groups than in normal men).

There are other domains for 17OH investigations as : male infertility, girls with peripubertal virilization, children with premature adrenache (in these cases, the values of 17OH are increased without or after ACTH stimulation).

IV. PRINCIPLES OF THE METHOD

A fixed amount of ^{125}I labelled steroid competes with the steroid to be measured present in the sample or in the calibrator for a fixed amount of antibody sites immobilized on the wall of a polystyrene tube. Neither extraction nor chromatography is required because of the high specificity of the coated antibodies. After 3 hours incubation at 37°C , an aspiration step terminates the competition reaction. The tubes are then washed with 3 ml of washing solution and aspirated. A calibration curve is plotted and the 17OH concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with anti $17\text{-}\alpha\text{-OH-progesterone}$	2 x 48	purple	Ready for use			
<table border="1" data-bbox="119 593 255 638"> <tr> <td>Ag</td> <td>^{125}I</td> </tr> </table> TRACER: ^{125}I Iodine labelled $17\text{-}\alpha\text{-OH-progesterone}$ (HPLC grade) in phosphate/citric acid buffer with bovine serum albumin and azide (<0.1%)	Ag	^{125}I	1 vial 55 ml 190 kBq	Red	Ready for use	
Ag	^{125}I					
<table border="1" data-bbox="119 795 255 840"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero calibrator in human serum and azide (0.5%)	CAL	0	1 vial 3 ml	yellow	Ready for use	
CAL	0					
<table border="1" data-bbox="119 896 255 940"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators $17\text{-}\alpha\text{-OH-progesterone}$ N = 1 to 5 (see exact values on vial labels) in human serum and azide (0.5%)	CAL	N	5 vials 0.5 ml	yellow	Ready for use	
CAL	N					
<table border="1" data-bbox="71 1064 327 1108"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (TRIS-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="87 1176 279 1220"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum and thymol	CONTROL	N	2 vials lyophilized	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note : Use the zero calibrator for sample dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 25 μl , 500 μl and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Water bath at 37°C
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Controls :** Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C .
- After reconstitution, controls are stable for one week at 2 to 8°C . For longer storage periods, aliquots should be made and kept at -20°C for 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original

well-closed vial at 2 to 8°C .

- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at $2\text{-}8^\circ\text{C}$.
- If the test is not run within 48 hrs., storage at -20°C is recommended.
- Avoid successive freezing and thawing.
- Serum and heparinized plasma provide similar results (avoid citrate samples) : $Y(\text{Serum}) = 0.92 \times (\text{hep. plasma}) + 0.10$ $r = 1.0$ $n = 18$
- EDTA plasma provides 15% lower results than serum and heparinized plasma : $Y(\text{Serum}) = 1.16 \times (\text{EDTA plasma}) - 0.10$ $r = 0.99$ $n = 18$

X. PROCEDURE

A. Handling notes

- Do not use the kit or components beyond expiry date.
- Do not mix materials from different kit lots.
- Bring all the reagents to room temperature prior to use.
- Thoroughly mix all reagents and samples by gentle agitation or swirling.
- Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
- Respect the incubation times.
- Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

- Label coated tubes in duplicate for each calibrator, sample and control. For the determination of total counts, label 2 normal tubes
- Briefly vortex calibrators, samples and controls and dispense 25 μl of each into respective tubes.
- Dispense 0.5 ml of ^{125}I Iodine labelled $17\text{-}\alpha\text{-OH-progesterone}$ into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand, to liberate any trapped bubbles.
- Incubate for 3 hours at 37°C .
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 3 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

- Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B₀(%)) values for each calibrator point as a function of the $17\text{-}\alpha\text{-OH-progesterone}$ concentration of each calibrator point, reject obvious outliers.
- Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.
- By interpolation of the sample (B/B₀(%)) values, determine the $17\text{-}\alpha\text{-OH-progesterone}$ concentrations of the samples from the reference curve.
- For each assay, the percentage of total tracer bound in the absence of unlabelled $17\text{-}\alpha\text{-OH-progesterone}$ (B₀/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

17OH-RIA-CT		cpm	B/Bo (%)
Total count		48579	
Calibrator	0.00 ng/ml	19831	100.0
	0.15 ng/ml	15587	78.6
	0.53 ng/ml	10398	52.4
	1.02 ng/ml	7074	35.7
	3.10 ng/ml	2930	14.8
	11.10 ng/ml	1028	5.2

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 0.02 ng/ml.

B. Specificity

The specificity was estimated by spiking a pool of 17OH samples (± 0.5 ng/ml) with steroids which might be present in patient samples.

Compound	Added amount (ng/ml)	Cross-Reactivity (%)
17OH-Progesterone	-	100
Progesterone	500	0.7
17- α -hydroxypregnenolone	1000	0.3
21-deoxycortisol	1000	0.5
Pregnenolone	500	0.022
11-deoxycortisol	500	0.5
Corticosterone	500	ND
11-deoxycorticosterone	500	0.035
Cortisol	2500	0.008
Testosterone	5000	ND
Androstenedione	5000	0.001
Estradiol	5000	0.001

C. Precision

INTRA-ASSAY

INTER-ASSAY

Serum	N	<X> \pm SD (ng/ml)	CV (%)	Serum	N	<X> \pm SD (ng/ml)	CV (%)
A	10	0.43 \pm 0.03	6.2	A	20	0.77 \pm 0.07	9.2
B	10	2.35 \pm 0.13	5.6	B	20	1.97 \pm 0.10	5.2
C	10	7.55 \pm 0.38	5.1				

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST

Serum	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
1	1/1	-	3.83
	1/2	1.92	2.16
	1/4	0.96	1.13
	1/8	0.48	0.53
	1/16	0.24	0.25
2	1/1	-	6.65
	1/2	3.33	3.32
	1/4	1.66	2.12
	1/8	0.83	0.88
	1/16	0.42	0.39

Samples were diluted with the zero calibrator.

RECOVERY TEST

Sample	added 17OH (ng/ml)	Recovered 17OH (ng/ml)	Recovered (%)
1	6.65	6.10	92
	5.17	4.32	84
	2.93	2.69	92
	1.25	1.22	98
	0.92	0.87	95
2	7.37	7.23	98
	5.90	5.78	98
	3.60	3.57	99
	2.01	2.14	107
	1.69	1.65	98

Conversion factor :

From ng/ml to nmol/L : $\times 3.03$

From nmol/L to ng/ml : $\times 0.33$

To the best of our knowledge, no international reference material exists for this parameter.

Reference method: 50 patient serum samples (values from 0.3 ng/ml to 10.0 ng/ml) were assayed along with the tritiated reference method. The linear regression analysis is quite good.

17OH-RIA-CT = 0.99 (^3H reference) + 0.06 ng/ml $r = 0.96$

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 36 minutes after the calibrator has been added to coated tubes.

TIME DELAY

Serum (ng/ml)	0'	12'	24'	36'
Serum 1	0.50	0.58	0.56	0.56
Serum 2	0.86	0.86	0.93	0.83
Serum 3	1.22	1.25	1.24	1.26
Serum 4	1.86	1.90	1.93	1.93
Serum 5	3.29	3.70	3.46	3.48
Serum 6	4.52	4.76	5.18	4.65

XIV. INTERNAL QUALITY CONTROL

- The percentage of total tracer bound in the absence of unlabelled 17- α -OH-progesterone (BO/T) must be $> 25\%$.
- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

	Concentration range (2.5 to 97.5% percentiles) (ng/ml)	Number of subjects
Normal males	0.59 – 3.44	72
Normal Females		
. Follicular phase	0.11 - 1.08	48
. Luteal phase	0.95 - 5.00	49
Pregnancy		
. First trimester	2.50 – 9.78	43
. Second trimester	3.40 – 8.50	31
. Third trimester	4.53 – 18.86	35

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons. Purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.

3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
Calibrators (0-5)	-	25	-
Samples, controls	-	-	25
Tracer	500	500	500
Incubation	3 hours at 37°C		
Separation	-	aspirate	
Working Wash solution		3.0 ml	
Separation		aspirate carefully	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP1409	P.I. Number : 1700475/en	Revision nr : 110217/1
-------------------------------------	-----------------------------	---------------------------

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1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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XVII. IRODALOM

- ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
- BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
- CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.
- DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
- GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenance.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
- NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
- SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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XVII. BIBLIOGRAFIA

1. ANDO S. et al (1983)
Plasma levels of 17-OH-progesterone and testosterone in patients with varicoceles.
Acta Endocrinologica, 102,463-469.
2. BONACCORSI A. et al (1987)
Male infertility due to congenital adrenal hyperplasia: testicular biopsy findings, hormonal evaluation, and therapeutic results in three patients.
Fertility and Sterility, 47, 4, 664-670.
3. CACCIARE et al (1987)
Neonatal screening program for congenital adrenal hyperplasia in a homogeneous Caucasian population.
Annals New York Academy of sciences 85-88.

4. DRAFTA D. et al (1982)
Plasma steroids in benign prostatic hypertrophy and carcinoma of the prostate.
J. Steroid Biochem. 17, 689-693.
5. GRANOFF A. et al (1985)
17-hydroxyprogesterone responses to adrenocorticotropin in children with premature adrenence.
Journal of Clinical Endocrinology and Metabolism, 3, 409-415.
6. NEW M. et al (1983)
Genotyping steroid 21-hydroxylase deficiency hormonal reference data.
Journal of Clinical Endocrinology and Metabolism 57, 320-326.
7. SOLYOM J. et al (1987)
Detection of late-onset adrenal hyperplasia in girls with peripubertal virilization.
Acta Endocrinologica (Copenh.) 115, 413-418.

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
	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
U L I	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
WASH SOLN	Wash buffer

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