IGF-1-RIA-CT

KIP1588

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

Read entire protocol before use.

IGF-1-RIA-CT

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human Insulin-like Growth Factor-1 (IGF-1) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource IGF-1-RIA-CT Kit
- B. **Catalog number :** KIP1588 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. Biological activities

Insulin-like growth factor 1 (IGF-1) or Somatomedin-C (SM-C) is a basic 70 amino acid single chain polypeptide (MW : 7649 Da) similar to proinsulin (50% sequence homology), and to the other well-characterized member of the somatomedin family : IGF II (67AA, 70 % sequence homology with IGF-1). IGF-1 is the most important factor, which mediates the growth promoting actions of growth hormone, a pituitary hormone with highly fluctuating blood levels due to pulsatile release. The blood concentration of IGF-1 is more stable due to the binding to carrier proteins. The concentration of the predominant binding protein (MW 53000) as well as the production of IGF-1, are regulated by growth hormone. IGF-1 is produced by the liver, and other tissues, and it has endocrine, paracrine and autocrine activities. It stimulates growth and regulates differentiation of various tissues, displays insulin-like activities and promotes cartilage growth. Although GH is the most important factor controlling IGF-1 secretion and concentration, other factors are also determinant: the age (with a peak at adolescence), the sex, the nutritional status, and other hormones (oestrogen, thyroxin, prolactin, ...). Specific trophic stimuli mainly control IGF-1 secretion in the local microenvironment of a particular organ (paracrine activities), while blood IGF-1 concentration is the most important variable for balanced systemic growth (endocrine activities).

B. Clinical applications

- **Growth retardation:** Growth retardation may be due to several causes, among which deficient GH production (hypopituitarism), which is associated with low IGF-1 blood levels. Because of the difficulties to get interpretable results from GH measurements (by dynamic multiple or stimulation tests), the determination of the stable IGF-1 concentration in plasma is often considered as a simple screening test to evaluation "GH impregnation" of the patient before deciding more extensive investigations. In several clinical situations with impaired growth, low IGF-1 levels may be observed despite normal or high GH production (i.e. malnutrition, chronic diseases states, some genetic dwarfs like Pygmies, ...). Interestingly, children with discrete GH neuro-secretory dysfunction may display low IGF-1 values despite normal GH levels by conventional testing. The results of IGF-1 assay must be interpreted cautiously by considering the normal variations of IGF-1 during childhood and adolescence (see Rosenfeld et al).
- **Acromegaly:** IGF-1 levels are elevated in acromegaly (excess production of GH) and may serve as an indicator of disease severity. Results are more readily interpreted because the normal values are more easily defined in adults. IGF-1 measurements are also useful to monitor treatment.
- **Research:** The IGF-1 RIA kit is an invaluable tool to study the modifications of this growth factor during physiologic (i.e. pregnancy) or pathologic (i.e. diabetes) situations, and the local regulation of IGF-1 production in relation to its paracrine and autocrine activities (wound healing, organ regeneration, neoplastic growth, foetal development, gonadal regulation, etc).

IV. PRINCIPLES OF THE METHOD

In the present kit, DIAsource has introduced a pre-treatment step in order to improve the clinical performance of the assay. It is well established that the binding proteins interfere with the radioimmunoassay for IGF-1. A fixed amount of ^{125}I labelled IGF-1 compete with the IGF-1 to be measured present in the sample or in the calibrator for a fixed amount of antibody sites immobilized on the wall of a polystyrene tube. After 2 hours incubation at room temperature with shaking, an aspiration step terminates the competition reaction. The tubes are then washed with 2 ml of washing solution and aspirated. A calibration curve is plotted and the IGF-1 concentrations of the samples are determined by dose interpolation from the calibration curve.

V. REAGENTS PROVIDED

Reagents	96 Test Kit	Colour Code	Reconstitution			
 Tubes coated with anti IGF-1	2 x 48	green	Ready for use			
<table border="1" data-bbox="71 649 271 694"> <tr> <td>Ag</td> <td>125I</td> <td>CONC</td> </tr> </table> TRACER: ^{125}I Iodine labelled IGF-1 (HPLC grade) in phosphate buffer with bovine albumin and azide (<0.1%)	Ag	125I	CONC	1 vial 0.4 ml 300 kBq	red	Dilute 101x with Tracer Buffer (see section VII. C)
Ag	125I	CONC				
<table border="1" data-bbox="95 817 287 862"> <tr> <td>TRACER</td> <td>BUF</td> </tr> </table> Tracer Buffer in phosphate buffer with bovine casein and azide (<0.1%)	TRACER	BUF	1 vial 25 ml	black	Ready for use	
TRACER	BUF					
<table border="1" data-bbox="127 963 255 1008"> <tr> <td>CAL</td> <td>0</td> </tr> </table> Zero calibrator in phosphate buffer with ovalbumin and gentamycin	CAL	0	1 vial lyophilised	yellow	Add 1 ml distilled water	
CAL	0					
<table border="1" data-bbox="119 1086 263 1131"> <tr> <td>CAL</td> <td>N</td> </tr> </table> Calibrators - N = 1 to 5 (see exact values on vial labels) in phosphate buffer with ovalbumin and gentamycin Calibrators are prediluted	CAL	N	5 vials lyophilised	yellow	Add 1 ml dilution buffer	
CAL	N					
<table border="1" data-bbox="71 1265 319 1310"> <tr> <td>PRE</td> <td>SOLN</td> </tr> </table> Pre-treatment Solution containing HCl 0.1N	PRE	SOLN	1 vial 5 ml	black	Ready for use	
PRE	SOLN					
<table border="1" data-bbox="71 1366 319 1411"> <tr> <td>DIL</td> <td>BUF</td> </tr> </table> Dilution buffer containing Tris-HCl buffer with bovine casein and azide (<0.1%)	DIL	BUF	1 vial 50 ml	green	Ready for use	
DIL	BUF					
<table border="1" data-bbox="63 1500 327 1556"> <tr> <td>WASH</td> <td>SOLN</td> <td>CONC</td> </tr> </table> Wash solution (TRIS-HCl)	WASH	SOLN	CONC	1 vial 10 ml	brown	Dilute 70 x with distilled water (use a magnetic stirrer).
WASH	SOLN	CONC				
<table border="1" data-bbox="79 1601 279 1657"> <tr> <td>CONTROL</td> <td>N</td> </tr> </table> Controls - N = 1 or 2 in human serum with thymol	CONTROL	N	2 vials lyophilised	silver	Add 0.5 ml distilled water	
CONTROL	N					

Note : 1. Use the dilution buffer for sample dilutions.
2. 1 ng of the calibrator preparation is equivalent to 1 ng of the NIBSC 1st IRR 87/518.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 50 μl , 100 μl , 200 μl , 500 μl , 1 ml and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Plastic tubes for pre-treatment of samples
- Vortex mixer
- Magnetic stirrer
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Tube shaker (400 rpm)
- Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators:** Reconstitute the zero calibrator with 1.0 ml distilled water and other calibrators with 1.0 ml dilution buffer.
- Controls:** Reconstitute the controls with 0.5 ml distilled water.
- Tracer solution :** Prepare an adequate volume of Tracer solution by adding 20 μl of Tracer to 2 ml of Tracer Buffer. Use a vortex to homogenize.

Example Tracer dilution

Number of tubes	Tracer (μl)	Tracer buffer (ml)	Tracer solution volume (ml)
24	50	5.0	5.050
48	100	10.0	10.100
72	150	15.0	15.150
96	200	20.0	20.200

Extemporaneous dilution of the tracer is required

- Working Wash solution:** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months. Avoid successive freezing and thawing.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Freshly prepared Tracer solution should be used on the same day.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATIONS

- Serum samples must be kept at 2 to 8°C.
- If the test is not run within 48 hrs, storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use. Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times. Prepare a calibration curve for each run, do not use data from previous runs.

B. Pre-treatment step

- Label one plastic tube for each sample and control.
- Dispense 50 μl of each sample and control into the tube.
- Add 50 μl of pre-treatment solution into this tube.
- Vortex each tube during 5 seconds.
- Incubate 30 minutes at room temperature.
- Add 1 ml of the dilution buffer into the tube.
- Vortex each tube.

C. Procedure

- Label coated tubes in duplicate for each calibrator, control and sample. For the determination of total counts, label 2 normal tubes.
- Briefly vortex calibrators, controls and samples and dispense 100 μl of each into respective tubes.
- Dispense 200 μl of Tracer solution into each tube, including the uncoated tubes for total counts.
- Shake the tube rack gently by hand to liberate any trapped air bubbles.
- Incubate for 120 minutes at room temperature on a tube shaker (400 rpm).
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts) and aspirate (or decant). Avoid foaming during the addition of the Working Wash solution.
- Let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0(\%) = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

3. Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B₀(%)) values for each calibrator point as a function of the IGF-1 concentration of each calibrator point. Reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4-parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B₀ (%)) values, determine the IGF-1 concentrations of the samples from the calibration curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled IGF-1 (B₀/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

IGF-1	cpm	B/Bo (%)
Total count	71038	
Calibrator		
0.0 ng/ml	38117	100.0
33.0 ng/ml	35753	93.8
81.4 ng/ml	30816	80.8
228.8 ng/ml	19914	52.2
640.2 ng/ml	10525	27.6
1529.0 ng/ml	5331	14.0

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 3.4 ng/ml.

B. Specificity

The percentages of cross-reaction estimated by comparison of the concentration yielding a 50% inhibition are respectively:

Compound	Cross-Reactivity (%)
IGF-1	100.0
IGF-II	0.7
Insuline	ND
GH	ND

Note: this table shows the cross-reactivity for the anti IGF-1.

C. Precision

INTRA-ASSAY PRECISION

INTER-ASSAY PRECISION

Serum	N	<X> ± SD (ng/ml)	CV (%)	Serum	N	<X> ± SD (ng/ml)	CV (%)
A	20	36.1 ± 3.3	9.1	A	17	129.1 ± 11.6	9.0
B	20	81.4 ± 1.9	1.9	B	17	362.5 ± 14.9	4.1
C	20	402.8 ± 6.8	1.7				

SD: Standard Deviation; CV: Coefficient of variation

D. Accuracy

DILUTION TEST

Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
Serum A	1/1	-	818
	1/2	409	418
	1/4	205	216
	1/8	102	110
	1/16	51	54
Serum B	1/1	-	501
	1/2	251	250
	1/4	125	130
	1/8	63	62
	1/16	31	26

Samples were diluted with the dilution buffer.

RECOVERY TEST

Sample	Added IGF-1 (ng/ml)	Recovered IGF-1 (ng/ml)	Recovered (%)
C1	17	18	108
C2	53	55	104
C3	160	164	103
C4	382	411	107

Conversion factor:

From ng/ml to nmol/L: x 0.1307

From nmol/L to ng/ml: x 7.649

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

Normal subjects

Age Group	MALES (ng/ml)			FEMALES (ng/ml)		
	Mean	Range	N	Mean	Range	N
0 - 2 years	102	73-184	14	98	59-143	3
3 - 5 years	124	103-189	7	159	84-447	9
6 - 8 years	177	115-249	9	276	79-432	10
9 - 11 years	362	181-656	9	272	175-445	8
12 - 14 years	315	168-557	10	418	202-1101	19
15 - 17 years	409	224-592	9	414	138-658	25
18 - 20 years	330	190-390	4	355	144-519	10
21 - 30 years	313	235-408	6	310	191-478	52
31 - 40 years	225	154-270	10	279	180-437	53
41 - 50 years	222	160-318	15	233	123-406	39
51 - 60 years	195	144-286	26	217	122-327	36
> 60 years	171	94-245	23	186	91-320	27

Rem: ranges are based on min - max values.

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.

8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS µl	CALIBRATORS µl	SAMPLE(S) CONTROLS µl
PRE-TREATMENT			
Samples, controls	-	-	50
Pre-treatment solution	-	-	50
Incubation	Vortex 5 seconds and incubate 30 minutes at room temperature		
Dilution buffer	-	-	1000
Shaking	Vortex		
INCUBATION			
Calibrators (0 to 5)	-	100	-
Pre-treated Samples, controls	-	-	100
Tracer solution	200	200	200
Incubation	120 minutes at room temperature with shaking at 400 rpm		
Separation	-	Aspirate (or decant)	
Working Wash solution		2.0 ml	
Separation		Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP1588	P.I. Number : 1700868/en	Revision nr : 120823/1
-------------------------------------	-----------------------------	---------------------------

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8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research, 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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XVII. LITERATUR

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.

11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.

9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IGF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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XVII. BIBLIOGRAFIA

1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOEGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.

12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research , 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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1. DAUGHADAY W.H. and ROTWEIN P. (1989)
Insulin-like growth factors I and II. peptide, messenger ribonucleic acid and gene structures, serum and tissue concentrations.
Endocrine Rev., 10(1) : 68-91.
2. FROESCH E.R. and ZAPF J. (1985)
Insulin-like growth factors and insulin : comparative aspects.
Diabetologia, 28 : 485-493.
3. FURLANETTO R.W., UNDERWOOD L., VAN WYK J.J., D'ERCOLE A.J. (1977)
Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay.
J. Clin. Invest. 60 : 648.
4. ROSENFELD R.G., WILSON D.M., LEE P.D.K. and HINTZ R.L. (1986)
Insulin-like growth factors I and II in evaluation of growth retardation.
J. Pediatr., 109 : 428.
5. RUDMAN D., KUTNER M.H. and CHAWLA R.K. (1985)
The short child with subnormal plasma somatomedin-C.
Pediatric Res., 19(10) : 975-980.
6. ALBERTSSON-WIKLAND K. and HALL K. (1987)
Growth hormone treatment in short children : relationship between growth and serum insulin-like growth factor I and II levels.
J. Clin. Endocrinol Metab., 65 : 671.
7. WASS J.A.H., CLEMMONS D.R., UNDERWOOD L.E., BARROW L., BESSER G.M. and VAN WYK J.J. (1982)
Changes in circulating somatomedin-C levels in bromocriptine treated acromegaly.
Clin. Endocrinol., 17 : 369-377.
8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and

radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.

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8. DAUGHADAY W.H., MARIZ I.K. and BLETHEN S.L. (1980)
Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites - a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol extracted serum.
J. Clin. Endocrinol. Metab., 51 : 781.
9. BREIERB.H., GALLAHER B.W. and GLUCKMAN P.D. (1991)
Radioimmunoassay for insulin-like growth factor-I : solutions to some potential problems and pitfalls.
Journal of Endocrinology, 128 : 347-357.
10. SCHOUTEN J.S.A.G. and al. (1993)
IGF1 : a prognostic factor of knee osteoarthritis.
British J. Rheumatol., 32 : 274-280.
11. GRONBAEK H., SKJAERBAEK C., NIELSEN B., FRYSTYK J., FOGH M.L., FLYVBJERG A., ORSKOV H. (1995)
Growth hormone and insulin-like growth factor-I: a suggested role in renal transplantation and graft vessel disease.
Transplant. Proceed., 27/3 : 2133-2136.
12. TSITOURAS P.D., ZHONG Y.G., SPUNGEN A.M., BAUMAN W.A. (1995)
Serum testosterone and growth hormone insulin-like growth factor-I in adults with spinal cord injury.
Hormone and metabolic Research, 27/6 : 287-292.
13. KOCH A., DORR H.G., GERLING S., BEHRENS R., BOHLES H.J. (1995)
Effect of growth hormone on IDF-I levels in patients with growth hormone deficiency and Wilson disease.
Hormone Research, 44/1 : 40-44.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACER BUF	Tracer buffer
µT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
CART	Cartridge
WASH SOLN	Wash buffer

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