25OH-VIT.D3-RIA-CT

KIP1961

For Informational/Research Purposes Only

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

For Informational/Research Purposes Only

LOT : 130729/1

Read entire protocol before use.

25OH-VIT.D3-RIA-CT

I. INTENDED USE

Radioimmunoassay for the in vitro quantitative measurement of 25OH-Vit.D3 in human serum and heparin plasma.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource 25OH-Vit.D3-Ria-CT Kit
- B. **Catalog number :** KIP1961 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.91

III. CLINICAL BACKGROUND

A. *Physiological function of 25OH-Vit.D₃*

25-Hydroxyvitamin D₃ (25OHD₃) is the trivial name of 9, 10-secocholesta-5, 7, 10(19)-triene-3 β, 25-diol. This secosteroid is produced in the liver by 25-hydroxylation of cholecalciferol or Vitamin D₃. 25OHD₃ is a precursor for other Vitamin D metabolites and has only a limited biological activity in itself. The most active derivative is 1 α,25-Hydroxyvitamin D₃, produced in the kidney (or placenta) by 1 α-hydroxylation of 25OHD₃. This hormonally regulated steroid stimulates the intestinal absorption of both calcium and phosphorus. It also stimulates bone resorption and mineralisation thereby preventing the development of rickets, osteoporosis or osteomalacia. This Vitamin D hormone might also be active in other tissues responsible for calcium transport (placenta, kidney, mammary gland, ...) and endocrine glands (beta-cells, parathyroid glands, ...). 25-hydroxyvitamin D₃ is a main precursor for the metabolites.

B. *Regulatory mechanism*

The production of 25-hydroxyvitamin D mainly occurs in the liver although other tissues (intestine, kidney) might perform the same hydroxylation. Although there might be some feedback inhibition of 25-hydroxyvitamin D on its own production, a higher substrate availability (Vitamin D₃) results in higher 25-hydroxyvitamin D production and higher circulating 25-hydroxyvitamin D concentrations in blood.

C. *Clinical applications*

This assay is of importance for the diagnosis of Vitamin D deficiency, insufficiency or intoxication.






IV. PRINCIPLES OF THE METHOD

At first calibrators, controls and samples (serum or heparin plasma) are extracted with acetonitrile.

A fixed amount of ^{125}I labelled 25OH Vitamin D_3 competes with the 25OH Vitamin D_3 from either extracted samples, controls or calibrators for a fixed amount of specific antibody sites immobilized to the lower and inner surface of plastic tubes.

After 2 hours incubation at room temperature, an aspiration step stops the competition reaction. The tubes are then washed with 2 ml washing solution and counted in a gamma counter.

V. REAGENTS PROVIDED

Reagents	Quantity	Colour Code	Reconstitution
 Tubes coated with anti 25OH-Vitamin. D_3	2 x 48	pink	Ready for use
 ^{125}I 25OH-Vitamin. D_3	1 vial lyophil. 160 kBq	red	Reconstitute extemporaneously by adding 6 ml of the tracer buffer
 Calibrator 0: horse serum / phosphate buffer, with gentamycin	1 vial lyophil.	yellow	Add 0.5 ml distilled water
 Calibrators 1-5 in horse serum / phosphate buffer, with gentamycin (see exact values on vial labels)	5 vials lyophil.	yellow	Add 0.5 ml distilled water
 Incubation Buffer with sodium azide (< 0.1%)	1 vial 45 ml	black	Ready for use
 Ethanol solution with sodium azide (< 0.1%)	1 vial 7 ml	red	Ready for use
 Acetonitrile	1 vial 25 ml	black	Ready for use
 Wash solution	1 vial 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
 Controls 1 and 2 in human serum with thymol	2 vials lyophil.	silver	Add 0.5 ml distilled water

Note : Use Calibrator 0 for dilution of samples with values above the highest calibrator before extraction step.
No international reference material is available.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. Distilled water
2. Pipettes for delivery of: 50 μl , 100 μl , 400 μl , 500 μl .
3. Vortex mixer
4. Magnetic stirrer
5. Glass tubes (12x75 mm) for extraction step
6. 5 ml automatic syringe (Cornwall type) for washing
7. Aspiration and washing device
8. Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).
10. Centrifuge operating at 800-1500 g

VII. REAGENT PREPARATION

- A. Calibrators :** Reconstitute the calibrators with 0.5 ml distilled water.
- B. Controls :** Reconstitute the controls with 0.5 ml distilled water.
- C. ^{125}I 25OH-Vit. D_3 :** Reconstitute with 6 ml of the tracer buffer.
- D. Working Wash solution :** Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, calibrators and controls are stable for one week at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for maximum 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- The reconstituted tracer has to be frozen after first use. Then, it is stable until the expiry date.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum and heparin plasma samples must be kept at 2-8°C.
- This kit is suitable for serum and heparin plasma samples. A correlation has been established between 23 serum and heparin plasma samples from same patients: Plasma = 0.95 Serum + 1.25, R = 0.89.
- If the test is not run within 48 h., storage at -20°C is recommended
- After thawing serum and heparin plasma samples should be mixed (Vortex), then centrifuged.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling. In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

Extraction step :

1. Label glass tubes (12x75 mm) for extraction: 6 calibrators, 2 controls and up to 40 samples in duplicate.
2. Dispense 100 μl of each calibrator, control or sample in the respective tubes.
3. Add 0.5 ml acetonitrile to each tube.
4. Mix for 7 seconds with a vortex.
5. Centrifuge for 5 minutes at room temperature (at 800-1500 g).

Incubation step :

1. Label coated tubes in duplicate for each calibrator, sample, control. For the determination of total counts, label 2 normal tubes.
2. Add 100 μl of the supernatant obtained after the extraction step in the corresponding tubes. Pipette tips have to be saturated with corresponding supernatant before the addition in the tube.
3. Dispense 400 μl Incubation Buffer in each tube, except those for total counts.
4. Add 50 μl Tracer in each tube, including total counts.
5. Incubate for 2 hours, under stirring (300-700 RPM), at room temperature.
6. Aspirate the content of each tube (except total counts).
7. Wash tubes twice with 2 ml Wash Solution and aspirate. Avoid foaming during the addition of the wash solution.
8. After the washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
9. Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

1. Calculate the mean of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0 \times 100 = \frac{\text{Counts (Calibrator or sample)}}{\text{Counts (Zero Calibrator)}} \times 100$$

3. Using a 3 cycle semi-logarithmic or logit-log graph paper plot the (B/B₀ x 100) values for each calibrator point as a function of the 25OH.D₃ concentration of each calibrator point, reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B₀ x 100) values, determine the 25OH.D₃ concentrations of the samples from the reference curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled 25OH.D₃ (B₀/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

25OH-Vit.D ₃	cpm	B/Bo x 100
Total count	31868	-
Calibrator 0 ng/ml	22218	100 %
3.5 ng/ml	17967	81 %
10.6 ng/ml	13916	63 %
21.0 ng/ml	9603	43 %
62.8 ng/ml	5113	23 %
145.2 ng/ml	4107	19 %

Note : 1 ng/ml = 2.5 pmol/ml

XIII. PERFORMANCE AND LIMITATIONS

A. Detection Limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent 25OH-Vitamin.D₃ concentration two standard deviations below the average counts at zero binding, was 1.2 ng/ml.

B. Specificity

The percentage of cross reaction estimated by comparison of the concentration yielding a 50 % inhibition are respectively :

Compound	Cross-Reactivity (%)
25OH-Vitamin D ₃	100
25OH-Vitamin D ₂	<0.3
1,25(OH) ₂ -Vitamin.D ₃	89*
1,25(OH) ₂ -Vitamin.D ₂	<0.01
Vitamin D ₃	<0.03
Vitamin D ₂	<0.03
24,25(OH) ₂ -Vitamin.D ₃	<0.8

* As 1,25(OH)₂-Vit.D₃ concentrations are practically 1000 times lower than 25-OH-Vit.D₃, this cross-reactivity is insignificant and does not interfere in this 25-OH-Vit.D₃ assay.

The assay performance is not affected by hemolysis (5 g/L haemoglobin tested) or bilirubinemia (0.25 g/L bilirubin tested).

C. Precision

INTRA-ASSAY				INTER-ASSAY			
Sample	N	<X> ± SD (ng/ml)	C.V. (%)	Sample	N	<X> ± SD (ng/ml)	C.V. (%)
A	49	19.1 ± 1.4	7.2	A	9	6.7 ± 0.5	7.2
B	43	10.4 ± 0.9	8.7	B	9	13.1 ± 1.0	7.3

SD : Standard Deviation; CV: Coefficient of variation

D. Accuracy

RECOVERY TEST			
Added 25OH-Vit D ₃ (ng/ml)	Measured 25OH-Vit D ₃ conc.		Recovery (%)
	Total (ng/ml)	Blanked (ng/ml)	
0	7,9		
8,9	17,6	9,7	109%
26,8	33,8	25,9	97%
38,6	48,6	40,7	105%

DILUTION TEST				
Sample	Dilution	Theoretical Concentration	Measured Concentration	Recovery (%)
		(ng/ml)	(ng/ml)	
Serum A	1/1	-	42,5	-
	1/2	21,3	20,9	98,4%
	1/4	10,6	11,2	105,4%
	1/8	5,3	5,7	107,3%
	1/16	2,7	2,6	97,9%
	1/32	1,3	1,3	97,9%

E. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 15 and 30 minutes after the calibrator has been added to the coated tubes.

TIME DELAY			
	0' (ng/ml)	15' (ng/ml)	30' (ng/ml)
Sample 1	9.5	9.2	10.1
Sample 2	36.8	33.9	36.2

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplo results of the samples should rely on Good Laboratory Practises

XV. EXPECTED VALUES

Dietary intake, race, season and age are known to affect the normal levels of 25OH.Vit.D₃.

The expected values given hereafter should not be considered as absolute. DIASource evaluated serum collected from 40 male and 39 female, healthy according to Ca, PTH and albumin values, from Western Europe. The age of the volunteers fell within the range of 17 - 58 years. Samples were collected during the months of December 2010 and January 2011. The mean for the population (n = 79) was 12.5 ng/mL, ranging from 4.1 to 28.7 ng/ml (based on 2.5 % to 97.5 % percentiles).

Each laboratory should establish its own range based on their local population.

Recent literature has suggested the following ranges for the classification of 25 OH Vitamin D status: Deficiency: 0-10 ng/mL; Insufficiency: 10-30 ng/mL; Sufficiency: 30 to 150 ng/mL; Toxicity: >150 ng/mL.

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves. All radioactive handling should be executed in a designated area, away from regular passage. A log book for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosafety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

XVII. BIBLIOGRAPHY

- ZERWEKH J.E. (2008)
Blood biomarkers of Vitamin D status.
Am. J. Clin. Nutr., 87(suppl):1087S-91S.
- HOLICK M.F. (2006)
Resurrection of Vitamin D deficiency and rickets.
J. Clin. Invest., 116:2062-2072.
- HEANEY R.P. (2000)
Vitamin D: how much do we need, and how much is too much.
Osteoporos. Int., 11:553-555.
- DAWSON-HUGHES B., HEANEY R.P., HOLICK M.F., LIPS P., MEUNIER P.J. (1997)
Prevalence of Vitamin D insufficiency in an adult normal population.
Osteoporos. Int., 7:439-443.
- BISCHOFF-FERRARI H.A., GIOVANNUCCI E., WILLETT W.C., DIETRICH T., DAWSON-HUGHES B. (2006)
Estimation of optimal serum concentrations of 25-hydroxyvitamin D for multiple health outcomes.
Am. J. Clin. Nutr., 84:18-28.
- HOLICK M.F. (2007)
Vitamin D deficiency.
N. Engl. J. Med., 357:266-281.

XVIII. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
EXTRACTION			
Calibrators	-	0.1	-
Samples/Controls	-	-	0.1
Acetonitrile	-	0.5	0.5
Vortex	Vortex for 7 seconds		
Centrifugation	5 minutes at 800-1500 g		
INCUBATION			
Supernatant of extraction	-	0.1	0.1
Incubation Buffer	-	0.4	0.4
Tracer	0.05	0.05	0.05
Incubation	2 hours at room temperature under stirring (300-700 RPM)		
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate	
Wash solution	-	2.0	
Separation	-	aspirate carefully	
Counting	Count tubes for 60 seconds		

DIAsource Catalogue Nr : KIP1961	P.I. Number : 1700543/en	Revision nr : 130729/1
-------------------------------------	-----------------------------	---------------------------

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
UJT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

For Information/Research Purposes Only

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