NSE-IRMA

KIP2471

For Informational/Research Purposes Only

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

For Informational/Research Purposes Only

LOT : 110218/1

Read entire protocol before use.

NSE-IRMA

I. INTENDED USE

Immunoradiometric assay kit for the in vitro quantitative measurement of human Neuron Specific Enolase (NSE) in serum.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource NSE-IRMA Kit
- B. **Catalog number :** KIP2471 : 96 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel : +32 (0)10 84.99.11 Fax : +32 (0)10 84.99.90

III. CLINICAL BACKGROUND

A. Biological activities

The glycolytic enzyme enolase (2-phospho-D-glycerate hydrolase, EC 4.2.1.11) exists as several dimeric isoenzymes ($\alpha\alpha$, $\alpha\beta$, $\alpha\gamma$, $\beta\beta$ and $\gamma\gamma$). The $\gamma\gamma$ and $\alpha\gamma$ enolase isoenzymes are also known as neuron-specific enolase (NSE) as they are produced in central and peripheral neurons and malignant tumors of neuroectodermal origin.

Elevated serum concentrations of NSE are detected in patients with medullary carcinoma of the thyroid, pancreatic islet cell tumor, pheochromocytoma carcinoid, neuroblastoma and small-cell lung carcinoma (SCLC).

B. Clinical applications

The diagnosis of lung cancer generally requires medical imaging, endoscopy, intra-operative findings and histology.





An important indication for tumor marker determinations in lung cancer is in assessing the efficacy of therapy and post-operative follow-up care.

In patients with SCLC, serum levels of NSE reflect the response to chemotherapy.

IV. PRINCIPLES OF THE METHOD

The DIASource NSE-IRMA is an immunoradiometric assay based on coated-tube separation. Mab1, the capture antibody, is attached to the lower and inner surface of the plastic tube. Calibrators or samples added to the tubes will at first show low affinity for Mab1. Addition of Mab2, the signal antibody labelled with ^{125}I , will complete the system and trigger the immunological reaction. After washing, the remaining radioactivity bound to the tube reflects the antigen concentration.

V. REAGENTS PROVIDED

Reagents	96 tests Kit	Colour Code	Reconstitution
 Tubes coated with anti NSE (monoclonal antibodies)	2 x 48	orange	Ready for use
 ^{125}I Anti-NSE- ^{125}I (monoclonal antibodies) in phosphate buffer with bovine serum albumin, azide (<0.1%) and inert red dye	1 vial 5.5 ml 700 kBq	red	Ready for use
 CAL O Calibrator 0 in bovin serum With thymol	1 vial lyophil.	yellow	Add 3ml distilled water
 CAL N Calibrators 1-5 in bovin serum with thymol (see exact value on vial labels)	5 vials lyophil.	yellow	Add 0.5 ml distilled water
 WASH SOLN CONC Wash solution (TRIS-HCl)	1 vial 10 ml	brown	Dilute 70x with distilled water (use a magnetic stirrer).
 CONTROL N Controls 1 and 2 in human serum and thymol	2 vials lyophil.	silver	Add 0.5 ml distilled water

Note: Use the zero calibrator for sera dilutions.

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 50 μl , 500 μl and 3 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Pipette for delivery of 5 to 10 ml of distilled water
- Vortex mixer
- Magnetic stirrer
- Tube shaker (400rpm)
- 5 ml automatic syringe (Cornwall type) for washing
- Aspiration system (optional)
- Any gamma counter capable of measuring ^{125}I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- Calibrators** : Reconstitute the calibrator 0 with 3 ml distilled water and the calibrators 1-5 with 0.5 ml.
- Controls** : Reconstitute the controls with 0.5 ml distilled water.
- Working Wash solution** : Prepare an adequate volume of Working Wash solution by adding 69 volumes of distilled water to 1 volume of Wash Solution (70x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the vial label, if kept at 2 to 8°C.
- After reconstitution of the calibrators and controls, aliquots should be made and kept at -20°C for maximum 3 months.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, tracer is stable until expiry date, if kept in the original well closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

NSE is also found in erythrocytes, plasma cells and platelets, it may be released into serum if separation from red cells does not occur within 60 minutes of venipuncture.

- Serum must be kept at 2-8°C.
- If the test is not run within 24 h., storage at -20°C is recommended.
- Avoid subsequent freeze-thaw cycles.
- Do not use plasma samples.

X. PROCEDURE

A. Handling notes

Do not use the kit or components beyond expiry date. Do not mix materials from different kit lots. Bring all the reagents to room temperature prior to use.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample. Thoroughly mix all reagents and samples by gentle agitation or swirling.

High precision pipettes or automated pipetting equipment will improve the precision. Respect the incubation times.

Prepare a calibration curve for each run, do not use data from previous runs.

B. Procedure

- Label coated tubes in duplicate for each calibrator, sample and control. For determination of total counts, label 2 normal tubes.
- Briefly vortex calibrators, samples, controls and dispense 50 μl of each into the respective tubes.
- Dispense 50 μl of anti-NSE- ^{125}I tracer into each tube, including the uncoated tubes for total counts.
- Shake the rack containing the tubes gently by hand to liberate any trapped air bubbles.
- Incubate for 2 hours at room temperature on a tube shaker (400 rpm).
- Aspirate (or decant) the content of each tube (except total counts). Be sure that the plastic tip of the aspirator reaches the bottom of the coated tube in order to remove all the liquid.
- Wash tubes with 2 ml Working Wash solution (except total counts). Avoid foaming during the addition of the Working Wash solution.
- Aspirate (or decant) the content of each tube (except total counts).
- Wash tubes again with 2 ml Wash solution (except total counts) and aspirate (or decant).
- After the last washing, let the tubes stand upright for two minutes and aspirate the remaining drop of liquid.
- Count tubes in a gamma counter for 60 seconds.

XI. CALCULATION OF RESULTS

- Calculate the mean of duplicate determinations.
- On semi logarithmic or linear graph paper plot the c.p.m. (ordinate) for each calibrator against the corresponding concentration of NSE (abscissa) and draw a calibration curve through the calibrator points, reject the obvious outliers.
- Read the concentration for each control and sample by interpolation on the calibration curve.
- Computer assisted data reduction will simplify these calculations. If automatic result processing is to be used, a 4-parameter logistic function curve fitting is recommended.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

NSE-IRMA		cpm	B/T (%)
Total count		275069	100
Calibrator	0.0 ng/ml	399	0.15
	2.5 ng/ml	2669	0.97
	7.8 ng/ml	6710	2.44
	26.0 ng/ml	20604	7.49
	78.0 ng/ml	57102	20.76
	270.0 ng/ml	144748	52.62

Since no international reference material is available for NSE antigen, BSE NSE-IRMA calibrator values are assigned against a set of in-house reference standards.

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twelve zero calibrators were assayed along with a set of the other calibrators. The detection limit, defined as the apparent concentration of the average count at zero binding plus two standard deviations, was 0.19 ng/ml.

B. Precision

INTRA ASSAY				INTER ASSAY			
Serum	Replicate	<X> ± SD (ng/ml)	CV (%)	Serum	Replicate	<X> ± SD (ng/ml)	CV (%)
A	20	12.9 ± 0.5	3.6	A	20	12.7 ± 0.5	3.8
B	20	104.4 ± 1.3	1.2	B	20	101.4 ± 3.8	3.8

SD : Standard Deviation; CV: Coefficient of variation

C. Accuracy

RECOVERY TEST

Added NSE (ng/ml)	Recovered NSE (ng/ml)	Recovery (%)
20.0	20.0	100.0
39.0	38.7	99.0
79.0	76.0	96.2
156.4	149.4	95.5

DILUTION TEST			
Sample	Dilution	Theoretical Concent. (ng/ml)	Measured Concent. (ng/ml)
1	1/1		164.7
	1/2	82.4	79.2
	1/4	41.2	39.4
	1/8	20.6	20.3
	1/16	10.3	10.2
	1/32	5.1	5.1
	1/64	2.6	2.5
	1/128	1.3	1.2
2	1/1		240.5
	1/2	120.3	120.4
	1/4	60.1	61.3
	1/8	30.1	29.9
	1/16	15.0	14.7
	1/32	7.5	6.6
	1/64	3.8	3.1
	1/128	1.9	1.4

Samples were diluted with zero calibrator .

D. Time delay between last calibrator and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the calibrator has been added to the coated tubes.

	TIME DELAY			
	0'	10'	20'	30'
S 1 (ng/ml)	11.5	11.8	11.4	11.9
S 2 (ng/ml)	45.1	44.7	43.6	43.1

E. Hook-effect

Samples containing 15000 ng/ml NSE give a result higher than the last calibration point.

F. Specificity

The monoclonal antibodies used are specific for the g-subunit of enolase. No measurable cross-reactions with other enolase have been observed.

XIV. LIMITATIONS

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with in vitro immunoassays. Patients routinely exposed to animals or animal serum products can be prone to this interference and anomalous values may be observed in case of the presence of heterophilic antibodies. Carefully evaluate the results of patients suspected of having these antibodies. If results are not consistent with other clinical observations, additional information should be required before diagnosis.

XV. INTERNAL QUALITY CONTROL

If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.

- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

XVI. REFERENCE INTERVALS

Among 332 apparently healthy individuals, 98% of the results were below 12.5 ng/ml.

For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history and other findings.

XVII. PRECAUTIONS AND WARNINGS

Safety

For *in vitro* diagnostic use only.

This kit contains ¹²⁵I (half-life: 60 days) ,emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HbsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections.

Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVIII. BIBLIOGRAPHY

1. LAMERZ R.
NSE (neuron-specific enolase).
Thomas L, ed. Clinical laboratory diagnostics. Verlag Frankfurt : TH-Books, 1998; 979-981.
2. LIIPPO KK, TERHO T.
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.
Acta Oncol 1991; 30:321-324.
3. OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTÉ G., NIKSSON K.
Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.
J. Pediatr Surg 1982; 17:381-385.
4. TAPIA FL., POLAK JM., BARBOSA AJ., et al.
Neuron-specific enolase is produced by neuroendocrine tumours.
Lancet 1981 ; 1 :808-811.
5. FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.
Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.
Cancer 1989; 63:1143-1149.
6. BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.
Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.
Clin Cancer Res 2000; 6:597-601.
7. EBERT W., HOPPE M., MULEY T., DRINGS P.
Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.
Anticancer RES 1997; 17:2875-2878.

XIX. SUMMARY OF THE PROTOCOL

	TOTAL COUNTS ml	CALIBRATORS ml	SAMPLE(S) CONTROLS ml
Calibrators (0-5) Samples, Controls Tracer	- - 0.05	0.05 - 0.05	- 0.05 0.05
Incubation	2 hours at room temperature with shaking at 400 rpm		
Separation Working Wash solution Separation Working Wash solution Separation	- - - - -	Aspirate (or decant) 2.0 Aspirate (or decant) 2.0 Aspirate (or decant)	
Counting	Count tubes for 60 seconds		

DIASource Catalogue No: KIP2471	P.I. Number : 1700863/en	Revision Number : 110218/1
------------------------------------	-----------------------------	-------------------------------

Revision date: 2011-02-18

For Informational/Research Purposes Only

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XVIII. BIBLIOGRAPHIE

- LAMERZ R.
NSE (neuron-specific enolase).
Thomas L, ed. Clinical laboratory diagnostics. Verlag Frankfurt : TH-Books, 1998; 979-981.
- LIIPPO KK, TERHO T.
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.
Acta Oncol 1991; 30:321-324.
- OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTÉ G., NIKSSON K.
Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.
J. Pediatr Surg 1982; 17:381-385.
- TAPIA FL., POLAK JM., BARBOSA AJ., et al.
Neuron-specific enolase is produced by neuroendocrine tumours.
Lancet 1981 ; 1 :808-811.
- FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.
Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.
Cancer 1989; 63:1143-1149.
- BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.
Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.
Clin Cancer Res 2000; 6:597-601.
- EBERT W., HOPPE M., MULEY T., DRINGS P.
Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.
Anticancer RES 1997; 17:2875-2878.

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XVIII. LITERATUR

- LAMERZ R.
NSE (neuron-specific enolase).
Thomas L, ed. Clinical laboratory diagnostics. Verlag Frankfurt : TH-Books, 1998; 979-981.
- LIIPPO KK, TERHO T.
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.
Acta Oncol 1991; 30:321-324.
- OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTJE G., NIKSSON K.
Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.
J. Pediatr Surg 1982; 17:381-385.
- TAPIA FL., POLAK JM., BARBOSA AJ., et al.
Neuron-specific enolase is produced by neuroendocrine tumours.
Lancet 1981 ; 1 :808-811.
- FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.
Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.
Cancer 1989; 63:1143-1149.
- BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.
Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.
Clin Cancer Res 2000; 6:597-601.
- EBERT W., HOPPE M., MULEY T., DRINGS P.
Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.
Anticancer RES 1997; 17:2875-2878.

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XVIII. BIBLIOGRAFIA

- LAMERZ R.
NSE (neuron-specific enolase).
Thomas L, ed. Clinical laboratory diagnostics. Verlag Frankfurt : TH-Books, 1998; 979-981.
- LIIPPO KK, TERHO T.
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.
Acta Oncol 1991; 30:321-324.
- OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTG G., NIKSSON K.
Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.
J. Pediatr Surg 1982; 17:381-385.
- TAPIA FL., POLAK JM., BARBOSA AJ., et al.
Neuron-specific enolase is produced by neuroendocrine tumours.
Lancet 1981 ; 1 :808-811.
- FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.
Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.
Cancer 1989; 63:1143-1149.
- BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.
Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.
Clin Cancer Res 2000; 6:597-601.
- EBERT W., HOPPE M., MULEY T., DRINGS P.
Monitoring of therapy in inoperable lung cancer patients by measurement of CYFRA 21-1, TPA, CEA and NSE.
Anticancer RES 1997; 17:2875-2878.

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2. LIIPPO KK, TERHO T.
Concomitant monitoring of serum neuron-specific enolase and creatine kinase BB in small cell lung cancer.
Acta Oncol 1991; 30:321-324.
3. OLDESTAD L., PAHLMAN S., LACKGREN G., LARSSON E., GROTTG G., NIKSSON K.
Neuron specific enolase : a marker for differential diagnosis of neuroblastoma and Wilms tumor.
J. Pediatr Surg 1982; 17:381-385.
4. TAPIA FL., POLAK JM., BARBOSA AJ., et al.
Neuron-specific enolase is produced by neuroendocrine tumours.
Lancet 1981 ; 1 :808-811.
5. FISCHBACH W., SCHWARZ-WALLRAUCH C., JANY B.
Neuron-specific enolase and thymidine kinase as an aid to the diagnosis and treatment monitoring of small-cell lung cancer.
Cancer 1989; 63:1143-1149.
6. BONNER JA., SLOAN JA., ROWLAND KM Jr., KLEE GG., KUGLER JW., MAILLIARD JA., et al.
Significance of neuron-specific enolase levels before and during therapy for small-cell lung cancer.
Clin Cancer Res 2000; 6:597-601.
7. EBERT W., HOPPE M., MULEY T., DRINGS P.
Monitoring of therapy in inoperable lung cancer patients by measurement of CYFERA 21-1, TPA, CEA and NSE.
Anticancer RES 1997; 17:2875-2878.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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