11-DESOXYCORTISOL- RIA-CT

KIPI20000

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

LOT : 110228/1

11-DESOXYCORTISOL-RIA-CT

Radioimmunoassay for the Quantitative Determination of
11- Desoxycortisol (11-DOC) in Human serum

KIPI20000

FOR IN VITRO DIAGNOSTIC USE ONLY

en

DIAsource ImmunoAssays SA - Rue du Bosquet 2, B-1348 Louvain-la-Neuve, Belgium - Tel: +32 10 84 99 11 - Fax : +32 10 84 99 90

1. INTENDED USE :

The 11-Desoxycortisol (Compound S) is an intermediate steroid in the glucocorticoids biosynthesis. Precursor of the cortisol, it comes from 17-hydroxyprogesterone after action of the 11 beta hydroxylase. This parameter is interesting for the diagnosis and follow up of treatment in case of surrenal enzymatic deficiency in 11 beta hydroxylase which is responsible for congenital surrenal hyperplasia in children and hyperandrogenics in adult women

Under hypothalamic-pituitary control via the Adrenocorticotrophic hormone (ACTH), the secretion of 11-Desoxycortisol follows a nycthemeral variation: it reaches a maximum in the morning (around 8 a.m.) and a minimal during the night (between 0 to 4 a.m.)


The measure of the 11-Desoxycortisol can be performed in serum using immunological competition methods (RIA).

The metopyrone inhibits the 11 beta hydroxylase and the conversion of the 11-Desoxycortisol in cortisol.

The metopyrone test is an indicator of the ACTH reserve.

2. PRINCIPLE OF THE METHOD :

The 11-Desoxycortisol RIA obeys the law of mass action according to the following equation :

Since the concentrations of ${}^{125}\text{I} - 11\text{-DOC}$ and coated antibodies are constant, the advancing state of the equation depends on the concentration of 11-DOC. The amount of ${}^{125}\text{I} - 11\text{-DOC}$ bound to the coated tube is inversely proportional to the concentration of 11 DOC in the sample.

Following the incubation, the tube is aspirated to remove excess unbound ${}^{125}\text{I}$ labelled 11-DOC.

Patient sample concentrations are read from a calibration curve.

3. MATERIAL PROVIDED AND STORAGE :

Stored at 2 - 8°C, the material can be used up to the expiration date printed on each label.

3.1.

--

2 x 48 Polystyrene tubes (12 x 75 mm)

coated with anti-11 Doc polyclonal antibodies.

Systematically, allow the coated tubes to reach room temperature before use. Store at 2-8°C.

3.2.

Ag	${}^{125}\text{I}$
----	--------------------

 ${}^{125}\text{I} - 11\text{ DOC E analog tracer}$ (yellow, 52 mL) :

1 bottle, containing < 5 μCi (185 kBq) of ${}^{125}\text{I}$ -labelled 11 DOC analog in protein-based buffer containing < 0.1% NaN_3 as a preservative. Store at 2-8°C.

3.3.

CAL	N
-----	---

11-DOC Calibrators

5 vials, 0.5 ml each, containing 0.3, 1.5, 5, 20 and 65 ng/ml 11-DOC calibrators in human serum with < 0.1% sodium azide as a preservative. Store at 2-8°C. Can be used up to the expiration date printed on each label.

3.4.

CAL	0
-----	---

11-DOC Zero Calibrator

One vial, 1 ml, containing 0 ng/ml 11-DOC zero calibrator in human serum with < 0.1% sodium azide as a preservative. Store at 2-8°C. Can be used up to the expiration date printed on the label.

3.5.

CONTROL	N
---------	---

Control 1 - 2

2 vials, lyophilised, 0.5 ml each, level I and level II containing low and high 11-DOC levels in human and horse serum with < 0.1% sodium azide as a preservative.

Store at 2-8°C. Can be used up to the expiration date printed on each label.

Before use, reconstitute the contents of the controls with 0.5 ml of distilled water.

3.6.

WASH	SOLN	CONC
------	------	------

Concentrated Wash Buffer

1 bottle concentrated buffered solution containing sodium azide (NaN_3 < 0.1 %). Poor the solution in 700 ml of distilled water.

Note: Conversion factor: 1 ng/ml = 2.886 nmol/L

1 nmol/L = 0.3465 ng/ml

4. MATERIAL REQUIRED BUT NOT PROVIDED :

- bench surfaces protected by absorbent paper to reduce the effects of radioactive spillage.
- waste disposal containers appropriately labelled and designed as suitable for solid or liquid radioactive materials.
- manual or automated precision micropipettes for dispensing samples or reagents without cross-contamination.
- absorbent paper.
- vacuum pump connected through a trap for aspiration.
- water bath
- a gamma scintillation counter.
- appropriate graph paper for plotting the results.

5. METHODOLOGY

5.1. Collection and handling of blood samples :

The blood sample may be collected into a dry tube.

After separation from the red blood cells, serum samples may be assayed immediately, within 24 hours if stored at 2 - 8°C, or later, after period up to several months if stored at -20°C. Repeating freezing and thawing must be avoided.

For samples from children of less than 2 years old, the 11-Desoxycortisol must be assayed immediately after an extraction step.

5.2. Assay procedure :

Reagents stored at 2° - 8° C. must be brought at room temperature prior to use. Do not mix reagents of different lots. Label the tubes for T (« Total Counts » do not use coated tubes) calibrators, samples and control sera.

Perform the assay in duplicate. Calibrators, controls and samples must be assayed at the same time.

1. Calibration curve:

Pipette 25 μl of each calibrator into the corresponding tubes.

2. Unknowns and control sera:

Pipette 25 μl of each sample or control sera into the corresponding tubes.

3. Add 500 μl of ${}^{125}\text{I} - 11\text{-DOC}$ tracer to each tube.

4. Vortex, cover and incubates 2.5 hours at $37 \pm 2^\circ\text{C}$.

5. Carefully aspirate or decant (before decanting, add 2 ml of washing solution to each tube.) the solution of all tubes. (Except total counts tubes).

6. Add 2 ml of washing solution to each tube. Aspirate or decant carefully.

7. Count the radioactivity fixed in each tube for at least 60 seconds.

Data processing :

Determine the mean count rate for each set of duplicate tubes.

Calculate the ratio B/B0 as follows :

$$B/B0 \% = [\text{CAL or Sample cpm} / B0 (\text{CAL } 0) \text{ cpm}] \times 100$$

Draw the calibration curve on semilogarithmic paper by plotting the ratio B/B0 % (linear scale) obtained for each calibrator versus its respective concentration expressed in ng/ml (logarithmic scale). 11-DOC concentrations in samples may be read directly from the calibration curve. If a computer is used to calculate the results, the data can be fitted to the appropriate equation : spline smoothed.

5.3. Example of a typical assay:

	Content (ng/ml)	cpm 1st duplicate	cpm 2nd duplicate	Mean cpm	B/Bo (%)	11-DOC (ng/ml)
Total counts	-	63015	62510	62762	-	-
CAL 0	0	31865	31445	31655	100	
CAL 1	0.3	27429	28004	27711	87.6	
CAL 2	1.5	20079	19803	19941	62.9	
CAL 3	5	13246	13152	13199	41.7	
CAL 4	20	6369	6509	6439	20.3	
CAL 5	65	3464	3480	3472	11	
C ₁ Low	-	22977	22569	22773	72	1
C ₂ High	-	13711	13699	13705	43.3	4.5
Sample	-	4172	4242	4214	13.3	45.2

Example of a typical assay, do not use for calculations

6. PERFORMANCE CHARACTERISTICS:

6.1. Specificity

Steroid	% Cross-reactivity
11-desoxycortisol	100
17 α Hydroxyprogesterone	5.6
Desoxycorticosterone	0.46
Progesterone	0.59
Cortisol	0.09
Estradiol-17 β	0.03
Androstenedione, Corticosterone, Dexamethasone, Cortisone, Testosterone, DHEA-S, Aldosterone	N.D.

6.2. Minimum detectable concentration of 11-DOC

The minimum detectable concentration has been assayed at 0.11 ng/ml and corresponds to the concentration given by two standards deviations below the mean cpm of 20 replicate determinations of the zero calibrators.

6.3. Recovery test:

When sera of known 11-DOC contents have their 11-DOC supplemented by addition of 11-DOC in equal volumes (1/1), a satisfactory correlation between theoretical and assayed 11-DOC is obtained.

Added 11-DOC (ng/ml) (1:1 in serum sample)	1.5	5	20	65
Theoretical (ng/ml)	0.85	2.6	10.1	32.6
Assayed 11-DOC (ng/ml)	0.70	2.5	10.2	30.2
% recovery	82.4	96.2	101	92.6

6.4. Dilution test :

The dilution test indicates that there is immunological identity between the 11-DOC present in the sample and the 11-DOC used to calibrate the test.

Dilution Factor	1	1/2	1/4	1/8	1/16	1/32	1/64
Assayed 11-DOC (ng/ml)	40.1	18.4	9.8	4.7	2.5	1.1	0.6
Expected 11-DOC (ng/ml)	-	20.05	10	5	2.5	1.25	0.63
% Recovery	-	92	98	94	100	88	96

6.5. Reproducibility:

	Within assay variation 5 replicates		Between assay variation 8 replicates	
	Mean (ng/mL)	% CV	Mean (ng/mL)	% CV
Pool 1	1.1	6.2	1.95	8.7
Pool 2	3.7	5.2	5.48	11.5
Pool 3	28.3	7.7	36.85	15.1

7. LIMITATION OF THE PROCEDURE

- The results obtained from this or any other diagnostic kit should be used and interpreted only in the context of an overall clinical picture.
- Do not use lipemic, haemolyzed, icteric or turbid specimens.

8. EXPECTED VALUES

Without Stimulation : < 7.2 ng/ml

After Stimulation with Metopyrone : 72 – 225 ng/ml

9. WARNING AND PRECAUTION

For in vitro diagnostic use only

Safety

This kit contains ¹²⁵I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons; purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiation safety procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
IVD	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 1251	Tracer
Ab 1251	Tracer
Ag 1251 CONC	Tracer concentrated
Ab 1251 CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
µP	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor

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