CE

DHT-RIA

KIP19900

For Informational Research Purposes Only

Distributed By:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P

Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

For Informational/Research Purposes Only

LOT : 120418/3

Read entire protocol before use.

DHT-RIA

I. INTENDED USE

Radioimmunoassay for the *in vitro* quantitative measurement of human 5- α -Dihydrotestosterone (DHT) in serum and plasma.

II. GENERAL INFORMATION

- A. **Proprietary name :** DIAsource DHT-RIA - Kit
- B. **Catalog number :** KIP19900 : 100 tests
- C. **Manufactured by :** DIAsource ImmunoAssays S.A.
Rue du Bosquet, 2, B-1348 Louvain-la-Neuve, Belgium.

For technical assistance or ordering information contact :
Tel: +32 (0) 10 84.99.11 Fax: +32 (0) 10 84.99.91

III. CLINICAL BACKGROUND

A. *Biological activities of 5 α -Dihydrotestosterone*

5 α -Dihydrotestosterone (DHT; 17 β -Hydroxy-5 α -androstane-3-one), the most potent naturally-occurring androgen, is produced from testosterone through the action of cholesterol 5 α -reductase [1]. The concentrations of 5 α -reductase are highest in certain peripheral tissues, including genital skin and hair follicles, and is localized intracellularly in apparent association with the nuclear membrane. DHT exerts its biological action by intracellular binding to the androgen receptor; this complex is then transferred to the nucleus where DNA-binding occurs with resultant effects on DNA transcription [2]. Most of the residual DHT undergoes intracellular metabolism to 3 α -androstane-20-one and 3 α -androstane-20-one glucuronide. Only a small proportion of DHT escapes into the peripheral circulation, where it is present primarily complexed to sex-hormone binding globulin [3].

B. *Clinical applications of 5 α -Dihydrotestosterone level measurement*

Testosterone causes virilization of the Wolffian ducts during fetal life, while DHT is responsible for the development of the male external genitalia and prostate, and is primarily responsible for the physical changes which occur during male sexual maturation. An autosomal-recessive genetic deficiency of 5 α -reductase, sometimes called male pseudohermaphroditism or pseudovaginal perineoscrotal hypospadias, leads to inadequate differentiation of DHT-dependent peripheral tissues. Male infants with this disorder have ambiguous genitalia and are often raised as females, although significant virilization may occur later in life presumably due to the natural increase in testosterone levels [2].

Measurement of DHT concentrations can be complicated by antibody cross-reactivity to testosterone. The DHT Double Antibody radioimmunoassay utilizes a sample oxydation / extraction procedure to remove most of the testosterone, coupled with a specific immunoassay for DHT.

IV. PRINCIPLES OF THE METHOD

The procedure follows the basic principle of radioimmunoassay where there is competition between the radioactive and non-radioactive antigen for a fixed number of antibody binding sites. The amount of (I-125)-labeled analyte bound to the antibody is inversely proportional to the concentration of unlabeled analyte present. The separation of free and bound antigen is achieved by using a double antibody system.

V. REAGENTS PROVIDED

Reagents	100 Tests	Colour code	Reconstitution
ANTISERUM Lyophilized anti-DHT rabbit serum with sodium azide (< 0.1%)	1 vial 10.2 ml Lyophilized	Blue	Reconstitute with 10.2 ml of buffer
Ag 125I TRACER: 125I labelled DHT in a protein based buffer with sodium azide (<0.1%)	1 vial 10.5 ml 185 kBq	Red	Ready to use
PEG Sheep anti-rabbit gamma globulin in a buffer containing polyethylene glycol and sodium azide (< 0.1%)	1 vial 105 ml	Green	Ready to use Mix well, with the use of a magnetic stirrer before use
KMnO4 Potassium permanganate solution	1 vial 50 ml	Black	Ready to use
BUF Buffer : phosphate buffered saline 0.1M with sodium azide (< 0.1%)	1 vial 100 ml	Black	Ready to use
CAL Calibrator : ethanolic stock solution, containing 100 ng/ml DHT.	1 vial	Yellow	Preparation of the standard curve : refer to VII B.
CONTROL N Controls : DHT in a buffer containing sodium azide (< 0.1%)	2 vials Lyophilized	Silver	Reconstitute with 2 ml Water. Exact concentrations are indicated on the labels

VI. SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

- Distilled water
- Pipettes for delivery of: 100 µl, 200 µl, 500 µl and 1 ml (the use of accurate pipettes with disposable plastic tips is recommended)
- Vortex mixer
- Magnetic stirrer
- Water bath at 37°C
- 12 x 75 mm conical propylene tubes
- 12 x 75 mm glass tubes.
- Aspiration system (optional)
- Refrigerated centrifuge capable of 1800 g.
- Hexane HPLC grade
- Ethanol HPLC grade
- Any gamma counter capable of measuring 125I may be used (minimal yield 70%).

VII. REAGENT PREPARATION

- A. Controls :** Reconstitute the controls with 2 ml distilled water. Allow 15 minutes to reconstitute; mix gently to avoid foaming. **Controls don't need to be extracted or oxidized.**

B. Standard curve preparation :

Prepare the following dilutions to obtain the curve:

Calibrator	Conc (pg/ml)	Calibrator (ml)	Buffer (ml)
C7	2500	Stock cal: 0.1	3.9
C6	1000	Stock cal: 0.1	9.9
C5	500	C6: 1.0	1.0
C4	200	C6: 0.2	0.8
C3	100	C6: 0.1	0.9
C2	50	C6: 0.1	1.9
C1	25	C2: 0.5	0.5

Buffer serves as "0" calibrator. Store the calibrators at 2-8°C for up to 12 weeks.

- C. Antiserum:** reconstitute the vial with 10.2 ml of buffer.

VIII. STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kit components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C.
- After reconstitution, controls are stable for 12 weeks at 2 to 8°C. For longer storage periods, aliquots should be made and kept at -20°C for 3 months.
- Freshly prepared calibrators are stable for up to 12 weeks at 2-8°C.
- After its first use, tracer is stable until expiry date, if kept in the original well-closed vial at 2 to 8°C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

IX. SPECIMEN COLLECTION AND PREPARATION

- Serum or plasma samples must be kept at 2-8°C.
- If the test is not run within 48 hrs., storage at -20°C is recommended.
- Avoid successive freezing and thawing.

X. PROCEDURE

- A. Handling notes**
Do not use the kit or components beyond expiry date.
Do not mix materials from different kit lots.
Bring all the reagents to room temperature prior to use.
Thoroughly mix all reagents and samples by gentle agitation or swirling.
The PEG reagent must be mixed at room temperature with the use of a magnetic stirrer at least 10 minutes before use.
Use a clean disposable pipette tip for addition of each different reagent and sample in order to avoid cross-contamination. High precision pipettes or automated pipetting equipment will improve the precision.
Respect the incubation times.
Prepare a calibration curve for each run, do not use data from previous runs.

- B. Extraction of the samples:**
Do not apply this procedure to calibrators or controls !

- Take **300 µl** (duplicate measurement) of serum or plasma in **glass** tubes.
- Add **600 µl** of KMnO4.
- Vortex about **2 x 10 seconds**.
- Incubate for **15 minutes** at room temperature (18-25°C).
- Add **3 ml** of a mixture of Hexane / Ethanol (98:2 vol/vol).
- Cap the tubes and shake thoroughly **3 x 30 seconds (do not use a vortex!)**.
- Centrifuge at **1800 g** for **5 minutes** at room temperature.
- For each sample, prepare and identify 2 polypropylene or glass tubes.
- Transfer **1 ml** of the upper organic layer into **each duplicate tube** and evaporate to dryness, preferentially under a stream of air or nitrogen..
- Add **200 µl** of buffer into each tube. Vortex **30 seconds**. Incubate for **20 minutes** at room temperature. Vortex again for **30 seconds**.
- The samples are now ready for the RIA determination.

C. RIA test

- Prepare and identify tubes (in duplicate) for: Total counts, non-specific binding (NSB), calibrators and controls. Use the same model of tubes than those containing the extracted samples (X.B.8.)
- Add **100 µl** of each calibrator (including 100 µl of buffer as calibrator 0) and **100 µl** of each control into their respective tubes.
- Add **100 µl** of buffer into each calibrator or control tube. (Add 300 µl of buffer into NSB tubes)
- Add **100 µl** of antiserum to each tube (calibrators, controls, samples). Do not add antiserum in Total counts and NSB tubes.

5. Add **100 µl** of 125I Tracer to all tubes.
6. Mix well, cover with a film and incubate for **30 minutes at 37°C** in a water bath.
7. Add **1 ml** PEG reagent to all tubes except Total counts. **The PEG reagent must be mixed thoroughly before use with a magnetic stirrer.**
8. Incubate at room temperature for 20 minutes.
9. Centrifuge at **1800 g, at 4-12°C, for at least 20 minutes.**
10. Aspirate or decant the supernatant by inverting the tubes. Blot the tubes to remove any remaining liquid.
11. Count the tubes for one minute in a gamma counter.

XI. CALCULATION OF RESULTS

1. Calculate the mean cpm of duplicate determinations.
2. Calculate the bound radioactivity as a percentage of the binding determined at the zero calibrator point (0) according to the following formula :

$$B/B_0 (\%) = \frac{\text{cpm (Calibrator or sample)} - \text{cpm NSB}}{\text{cpm (Zero Calibrator)} - \text{cpm NSB}} \times 100$$

3. Using a 3 cycle semi-logarithmic or logit-log graph paper, plot the (B/B₀(%)) values for each calibrator point as a function of the DHT concentration of each calibrator point, reject obvious outliers.
4. Computer assisted methods can also be used to construct the calibration curve. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.
5. By interpolation of the sample (B/B₀(%)) values, determine the DHT concentrations of the samples from the reference curve.
6. For each assay, the percentage of total tracer bound in the absence of unlabelled DHT (B₀/T) must be checked.

XII. TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

DHT-RIA	cpm	B/B ₀ (%)
Total count	72265	
NSB	2391	3.3
Calibrator		
0.0 pg/ml	30453	100.0
25.0 pg/ml	26395	86.7
50.0 pg/ml	23457	77.0
100.0 pg/ml	19476	63.9
200.0 pg/ml	14909	49.0
500.0 pg/ml	9639	31.7
1000.0 pg/ml	7180	23.6
2500.0 pg/ml	4568	15.0

XIII. PERFORMANCE AND LIMITATIONS

A. Detection limit

Twenty zero calibrators were assayed along with a set of other calibrators. The detection limit, defined as the apparent concentration two standard deviations below the average counts at zero binding, was 20 pg/ml.

B. Specificity

Compound	Cross-Reactivity (%)
DHT	100.000
Testosterone (after oxydation)	0.100
Estriol	0.030
Estradiol	0.005
Progesterone	0.004

C. Precision

INTRA-ASSAY

INTER-ASSAY

Serum	N	Mean (pg/ml)	CV (%)	Serum	N	Mean (pg/ml)	CV (%)
A	10	159.3	6.0	A	12	60.5	18.6
B	10	411.6	4.8	B	12	201.4	12.3
C	10	700.1	4.8	C	12	653.9	8.7

D. Accuracy

DILUTION TEST

Serum	Dilution	Theoretical Concent. (pg/ml)	Measured Concent. (pg/ml)
1	1/1	-	737.0
	1/2	369.0	373.4
	1/4	184.0	192.5
	1/8	92.0	88.0
	1/16	46.0	52.1

Samples were diluted with the buffer.

RECOVERY TEST

Added DHT (pg/ml)	Theoretical conc (pg/ml)	Recovered DHT (pg/ml)	Recovered (%)
0.0	-	311.1	-
50.0	180.5	185.1	102.5
100.0	205.5	210.5	102.4
200.0	255.6	260.4	101.9
500.0	405.6	429.5	105.9
1000.0	655.6	672.8	102.6

XIV. INTERNAL QUALITY CONTROL

- If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.
- If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Do not freeze-thaw more than twice.
- Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises.

XV. REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

	Concentration range (pg/ml)
Normal males	
Pre puberty	< 50
Males	250 - 800
Normal Females	
Follicular phase	50 - 200
Luteal phase	100 - 300
Menopausal	< 100

XVI. PRECAUTIONS AND WARNINGS

Safety

For in vitro diagnostic use only.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations.

This radioactive product can be transferred to and used only by authorized persons. Purchase, storage, use and exchange of radioactive products are subject to the legislation of the end user's country. In no case the product must be administered to humans or animals.

All radioactive handling should be executed in a designated area, away from regular passage. A logbook for receipt and storage of radioactive materials must be kept in the lab. Laboratory equipment and glassware, which could be contaminated with radioactive substances, should be segregated to prevent cross contamination of different radioisotopes.

Any radioactive spills must be cleaned immediately in accordance with the radiosity procedures. The radioactive waste must be disposed of following the local regulations and guidelines of the authorities holding jurisdiction over the laboratory. Adherence to the basic rules of radiation safety provides adequate protection.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with reagents (sodium azide as preservative). Azide in this kit may react with lead and copper in the plumbing and in this way form highly explosive metal azides. During the washing step, flush the drain with a large amount of water to prevent azide build-up.

Ethanol: the calibrator contains ethanol, which is highly flammable.

KMnO₄: this component is hazardous in case of skin or eye contact (irritant), of ingestion, of inhalation. In case of skin contact, flush immediately with water. In case of eye contact, check immediately for and remove any contact lenses, and flush eyes with plenty of water. Get medical attention immediately.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

XVII. BIBLIOGRAPHY

1. Miller WL: Molecular biology of steroid hormone synthesis. *Endocrin Rev* 9:295-318, 1988.
2. Pang S, Riddick L: Hirsutism. IN Lifshitz F: *Pediatric Endocrinology, A Clinical Guide, second edition*. Marcel Dekker, New York, 1990, pp. 259-291.
3. Wilson JD: Syndromes of androgen resistance. *Biol Reprod* 46:168-173, 1992.
4. Yalow R, Berson S: Introduction and general considerations. IN: Odell WD, Doughday WH (eds): *Principles of Competitive Protein Binding Assays*. J. B. Lippincott Co., Philadelphia, 1971, pp. 1-19.
5. Burtis CA, Ashwood ER: *Tietz Textbook of Clinical Chemistry, 2nd edition*. W.B. Saunders Company, Philadelphia, 1994, p.1821.

XVIII. SUMMARY OF THE PROTOCOL

	Total count	NSB	Calibrator (0-7)	Controls	Samples (extracted)
Calibrator	-	-	100 μl	-	-
Controls	-	-	-	100 μl	-
Sample (extracted)	-	-	-	-	200 μl
Buffer	-	300 μl	100 μl	100 μl	-
Antiserum	-	-	100 μl		
125I Tracer	100 μl				
30 min at 37°C (water bath)					
PEG	-	1 ml			
Mix and incubate 20 min at room temperature					
Centrifuge at least 20 min (1800 g; 4 - 12 °C)					
Aspirate or decant; count for 1 minute in a gamma counter					

DIAsource Catalogue Nr : KIP19900	P.I. Number : 1701100/en	Revision nr : 120418/1
--------------------------------------	-----------------------------	---------------------------

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	Used symbols
	Consult instructions for use
	Storage temperature
	Use by
LOT	Batch code
REF	Catalogue number
CONTROL	Control
I V D	In vitro diagnostic medical device
	Manufacturer
	Contains sufficient for <n> tests
WASH SOLN CONC	Wash solution concentrated
CAL 0	Zero calibrator
CAL N	Calibrator #
CONTROL N	Control #
Ag 125I	Tracer
Ab 125I	Tracer
Ag 125I CONC	Tracer concentrated
Ab 125I CONC	Tracer concentrated
	Tubes
INC BUF	Incubation buffer
ACETONITRILE	Acetonitrile
SERUM	Serum
DIL SPE	Specimen diluent
DIL BUF	Dilution buffer
ANTISERUM	Antiserum
IMMUNOADSORBENT	Immunoabsorbent
DIL CAL	Calibrator diluent
REC SOLN	Reconstitution solution
PEG	Polyethylene glycol
EXTR SOLN	Extraction solution
ELU SOLN	Elution solution
GEL	Bond Elut Silica cartridges
PRE SOLN	Pre-treatment solution
NEUTR SOLN	Neutralization solution
TRACEUR BUF	Tracer buffer
UJT	Microtiterplate
Ab HRP	HRP Conjugate
Ag HRP	HRP Conjugate
Ab HRP CONC	HRP Conjugate concentrate
Ag HRP CONC	HRP Conjugate concentrate
CONJ BUF	Conjugate buffer
CHROM TMB CONC	Chromogenic TMB concentrate
CHROM TMB	Chromogenic TMB solution
SUB BUF	Substrate buffer
STOP SOLN	Stop solution
INC SER	Incubation serum
BUF	Buffer
Ab AP	AP Conjugate
SUB PNPP	Substrate PNPP
BIOT CONJ CONC	Biotin conjugate concentrate
AVID HRP CONC	Avidine HRP concentrate
ASS BUF	Assay buffer
Ab BIOT	Biotin conjugate
Ab	Specific Antibody
SAV HRP CONC	Streptavidin HRP concentrate
NSB	Non-specific binding
2nd Ab	2nd Antibody
ACID BUF	Acidification Buffer
DIST	Distributor
TRAY	Incubation trays
PMSF	PMSF solution
	Protect from light
STRIP	Dot Strip
SUB	Substrate
EXTR SOLN CONC	Extraction Buffer Concentrate
CART	Cartridge
SAV HRP	Streptavidin HRP
PIPETTE	Pipette
WASH SOLN	Wash buffer

For Informative/Research Purposes Only

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