



## Evaluation ( Total AGT vs Intact AGT )

Date Prepared: 29th August 2019

Sample (Healthy Subject)	Diluted by 4,000 fold (ng/mL)			Diluted by 10,000 fold (ng/mL)		
	27412 (Total AGT)	27742 (Intact AGT)	Intact/Total (ratio)	27412 (Total AGT)	27742 (Intact AGT)	Intact/Total (ratio)
Serum 42	36,356	16,176	44.5%	39,280	15,760	40.1%
Serum 43	25,632	14,544	56.7%	26,230	14,080	53.7%
Serum 44	25,704	13,820	53.8%	28,130	14,430	51.3%
EDTA Plasma 42	28,612	14,168	49.5%	30,050	14,890	49.6%
EDTA Plasma 43	25,284	13,780	54.5%	26,790	14,830	55.4%
EDTA Plasma 44	21,240	11,892	56.0%	22,740	11,100	48.8%



Distributed By:  
**IBL-America, Inc.**  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432, USA  
[info@ibl-america.com](mailto:info@ibl-america.com)  
(888) 523 1246