


Antibody Specificity Evaluation

Rat new born brain extract

(1 mg/ml each)



Comment:

3 Rat new born brain extracts were used for W.B. Evaluation and targeted bands were observed at 62 and 66kDa using #11098 Ab and the bands were observed at 66 kDa using #29060.

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246

Email: info@IBL-America.com

Web: www.IBL-America.com