


## The concentrations in Serum and EDTA plasma of healthy subjects

Sample	ApoA5 ng/mL
5	192.2
6	189.2
8	396.8
11	107.5
12	132.7
28	213.3
29	80.2
30	278.6
32	263.6
33	334.5
34	327.8
35	217.6
36	162.9
37	295.8
38	177.0
39	202.0
40	153.3

Sample	ApoA5 ng/mL
5	217.7
6	264.3
8	396.8
11	149.7
12	145.0
28	205.0
29	110.4
30	187.8
32	272.6
33	335.8
34	345.4
35	207.9
36	214.3
37	277.9
38	221.6
39	218.2
40	155.9



Distributed by:



Immuno-Biological Laboratories, Inc.  
 8201 Central Ave NE, Suite P  
 Minneapolis, MN 55432

Toll-Free: 888-523-1246  
 Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
 Web: [www.IBL-America.com](http://www.IBL-America.com)