

Code No.27262 Human Periostin Assay Kit – IBL

Correlation data of Healthy Volunteer samples (Serum vs Heparin Plasma)

Serum

Sample	Dilute	Value	Result	Periostin ng/mL
1	10	0.449	1.74	17.35
2	10	0.504	1.99	19.88
4	10	1.311	6.13	61.30
5	10	1.170	5.36	53.61
7	10	1.824	9.04	90.44
6	10	0.577	2.33	23.31
8	10	1.059	4.77	47.67
10	10	0.415	1.58	15.81

Heparin Plasma

Sample	Dilute	Value	Result	Periostin ng/mL
1	10	0.439	1.69	16.89
2	10	0.452	1.75	17.48
4	10	1.155	5.28	52.80
5	10	1.171	5.37	53.66
7	10	1.702	8.34	83.36
6	10	0.566	2.28	22.79
8	10	1.000	4.46	44.56
10	10	0.401	1.52	15.18

