


Panel plasma and serum

	Plasma	Serum
11	105.8	83.6
12	67.0	57.5
13	135.9	120.2
14	151.3	139.7
15	359.9	277.8
16	254.7	224.2
17	344.5	267.9
18	365.5	212.4
19	92.6	81.7
20	151.8	137.4
21	487.6	433.2
22	416.7	379.1
23	160.2	165.7
24	462.2	347.1
25	148.3	138.4
26	189.8	167.9
27	120.7	131.8
28	112.4	97.7
29	61.5	56.3
30	170.4	157.7
Ave	217.9	183.9
SD	136.51	107.51
R	0.968	

ANGPTL8



Distributed by:



Immuno-Biological Laboratories, Inc.
 8201 Central Ave NE, Suite P
 Minneapolis, MN 55432

Toll-Free: 888-523-1246
 Email: info@IBL-America.com
 Web: www.IBL-America.com