



Antibody Specificity Evaluation

#28065 Anti-Human GBF1
(T1337 Phosphorylated) Rabbit
IgG Affinity Purify

#28063 Anti-Human
GBF1 (1844) Rabbit
IgG Affinity Purify



WB #28063 (Lot.0B-828)



W.B. #28065 (Lot: 0B-828)

Date Reported:	2 nd July 2020
Conducted by	K. Ishikawa
Result:	Pass

Test Details:

Cultured Huh7 cell line with stimulation 2-DG (-)/(+) and conducted W.B. using #28063 and #28065.

Comment:

The targeted bands were confirmed around expected size. As same as the reference (PMID: 18063581), 2-DG (+) lanes have stronger bands.

Lanes:

- 1: HuH-7 2DG(-)
- 2: HuH-7 2DG(+)
- 3: HuH-7 2DG(-) IP/#28063 (Lot: 0B-828)
- 4: HuH-7 2DG(+) IP/#28063(Lot: 0B-828)
- 5: HuH-7 2DG(-) IP/#28063 (Lot: 0B-806)
- 6: HuH-7 2DG(+) IP/#28063 (Lot: 0B-806)

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246

Email: info@IBL-America.com

Web: www.IBL-America.com