



# COX Antibodies

COX Antibodies							
Code No.		Name		Volume	IHC (※ F/P)	WB	
18511	①	Anti-Human	COX-1 (C24)	Rabbit IgG Affinity Purify	100 µg	—	○ 2-5 µg/mL
18521	②	Anti-Rat	COX-1 (C448)	Rabbit IgG Affinity Purify	100 µg	○ (※ MW) 5 µg/mL	○ 2-5 µg/mL
10211	③	Anti-Human	COX-2 (13H14)	Mouse IgG MoAb	100 µg	○ (※ MW) 2-4 µg/mL	○ 2-5 µg/mL
18515	④	Anti-Human	COX-2 (C295)	Rabbit IgG Affinity Purify	100 µg	○ (※ MW) 2-10 µg/mL	○ 2-5 µg/mL
18516	⑤	Anti-Human	COX-2 (C631)	Rabbit IgG Affinity Purify	100 µg	○ (※ MW) 2-5 µg/mL	○ 2-5 µg/mL
18525	⑥	Anti-Rat	COX-2 (C580)	Rabbit IgG Affinity Purify	100 µg	—	○ 2-5 µg/mL
18955	⑦	Anti-Rat	COX-2	Rabbit IgG Affinity Purify	100 µg	○ (※ MW) 5 µg/mL	○ 5 µg/mL

※ F/P : Can be Applied to Formalin Fixed Paraffin Embedded Tissue ※MW:MW Pre-treatment is necessary.



③、④、⑤-human lung adenocarcinoma,  
②、⑦- erbB2 transgenic mouse gallbladder carcinoma provided from Dr. Kiguchi (Anderson Cancer Center, Texas, USA) and Dr. Kawamoto (University of Tsukuba)

**Manufacturer :**

**Immuno-Biological Laboratories Co., Ltd.**  
 5-1 Aramachi, Takasaki-shi, Gunma 370-0831, JAPAN  
 TEL: 027-310-8040  
 FAX: 027-310-8045  
 E-mail: do-ibl@ibl-japan.co.jp  
 URL: http://www.ibl-japan.co.jp



*Distributed By:*

**Immuno-Biological Laboratories, Inc. (IBL-America)**  
 8201 Central Ave NE, Suite P, Minneapolis, MN 55432  
 Phone: (888) 523-1246  
 Fax: (763) 780-2988  
 Email: info@ibl-america.com  
 Website: www.ibl-america.com