

Product information

Information about other products is available at: www.ibl-america.com



Mouse Control

Mouse control suitable for internal quality control of the following IBL-America ELISA:

TSH mouse ELISA

(IB79180)

RUO

REF

IBDEV88MC

1. MATERIALS SUPPLIED, STORAGE AND STABILITY

Components	Content	Preparation	Storage at	Shelf life
Mouse Control 1/2: Two levels of TSH in serum	1 Set 2 vials à 1 ml	lyophilized; reconstitute each vial with 1 ml dist. H ₂ O 30 minutes before use	2 - 8 °C; after reconstitution store in aliquots at <-20 °C	6 months after reconstitution or until the expiration date

2. APPLICATION

According to the respective package insert of the IBL-America ELISA kits

IBDEV8822 TSH mouse

the control can be used for internal quality control.

The procedure calls for **equivalent sample volume** per test.

3. TARGET VALUES

The values provided here are intended as guidelines. Because of many sources of inter-laboratory variations each laboratory should develop its own expected values.

Parameter		Mean (ng/ml)	Target Range (ng/ml)
TSH	Control 1	1.2	0.7 - 1.7
	Control 2	2.6	1.6 - 3.7

SYMBOLS USED WITH IBL-AMERICA ASSAYS

Symbol	English	Deutsch	Française	Espanol	Italiano
	European Conformity	CE-Konformitäts-kennzeichnung	Conforme aux normes européennes	Conformidad europea	Conformità europea
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Consulte las Instrucciones	Consultare le istruzioni per l'uso
	In vitro diagnostic device	In-vitro-Diagnostikum	utilisation Diagnostic in vitro	Diagnóstico in vitro	Per uso Diagnostica in vitro
	For research use only	Nur für Forschungszwecke	Seulement dans le cadre de recherches	Sólo para uso en investigación	Solo a scopo di ricerca
	Catalogue number	Katalog-Nr.	Référence	Número de catálogo	No. di catalogo
	Lot. No. / Batch code	Chargen-Nr.	No. de lot	Número de lote	Lotto no
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Note warnings and precautions	Warnhinweise und Vorsichtsmaßnahmen beachten	Avertissements et mesures de précaution font attention	Tiene en cuenta advertencias y precauciones	Annoti avvisi e le precauzioni
	Storage Temperature	Lagerungstemperatur	Température de conservation	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Fecha de caducidad	Data di scadenza
	Legal Manufacturer	Hersteller	Fabricant	Fabricante	Fabbricante
<i>Distributed by</i>	Distributed by	Vertrieb durch	Distribution par	Distribución por	Distribuzione da parte di
	Version	Version	Version	Versión	Versione
	Single-use	Einmalverwendung	À usage unique	Uso único	Uso una volta