

## Total Vitamin D derivative - Amino



### PRODUCT SPECIFICATION

Cat # **5019703**

Batch\*: *[Batch dependent]*

Size 1 mg/ vial

Format: Dry

This product has been dried down in its pure form, and not lyophilized from a buffer solution. Therefore, small quantities might not be easy to visualize, being stuck to the bottom of the vial. Please follow carefully the instruction for dissolution, for a complete recovery of the product.

Purity:  $\geq 95\%$  (at 215 nm) and  $\geq 95\%$  (at 266 nm)

Reconstitution: Ethanol

Storage: -20°C before reconstitution

-20°C after reconstitution in inert gas environment

-80°C for long storage in inert gas environment

Expiry date: *[Batch dependent]*

\* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature

Distributed by:



Immuno-Biological Laboratories, Inc.  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432

Toll-Free: 888-523-1246  
Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Web: [www.IBL-America.com](http://www.IBL-America.com)