

Total Vitamin D derivative - Biotin conjugate



PRODUCT SPECIFICATION

Cat # **5019708**

Batch*: *[Batch dependent]*
Size: 1 mg/ vial
Format: Dry

This product has been dried down in its pure form, and not lyophilized from a buffer solution. Therefore, small quantities might not be easy to visualize, being stuck to the bottom of the vial. Please follow carefully the instruction for dissolution, for a complete recovery of the product.

Purity: $\geq 95\%$ (at 215 nm) and $\geq 95\%$ (at 266 nm)
Reconstitution: Ethanol
Storage: -80°C
Expiry date: *[Batch dependent]*

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com