

Code No. 10096

Anti-Neurofilament, 160K, 210K, (N-21-8-4) Mouse Ascites MoAbVolume : 0.5 mL

Introduction : Neurofilament is intermediate filament which expressed in nerve cells. And the protein consists of three subunits, 200 kDa, 160 kDa and 6.8 kDa. These subunits exists in phosphorylated form in axonal transport. This antibody reacts with the 200 kDa and 160 kDa subunits.

Antigen : Bovine Brain Neurofilament

Source : Mouse-Mouse hybridoma (abdominal)
(X63-Ag8.653 × BALB/c spleen cells)

Clone : N-21-8-4

Subclass : IgG₁

Purification : Purified with Protein A

Form : Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

How to use : 0.5 mL deionized water will be added to the product

Dilution : PBS (pH7.4) containing 1% BSA

Stability : Lyophilized product, 5 years at 2 - 8 °C
: Solution, 2 years at -20 °C

Application : This antibody can be used for immunohistochemistry with formalin fixed paraffin embedded tissues by several techniques such as Avidin Biotin Complex (ABC) method. The optimal dilution is about X200, however, the dilution rate should be optimized by each laboratory.

Specificity : Reacts with Human Neurofilament 160 kDa and 210 kDa

For research use only, not for use in diagnostic procedures.

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com