

Code No. 10098

**Anti-Human  
Platelets (KOR-P77) Mouse IgG MoAb**

Volume : 100  $\mu$ g

---

**Antigen** : Human Platelets

**Source** : Mouse-Mouse hybridoma (Supernatant)  
(NS-1 8Ag  $\times$  BALB/c spleen cells)

**Clone** : KOR-P77

**Subclass** : IgG<sub>1</sub>

**Purification** : Purified with Protein A

**Form** : Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use** : 1.0 ml distilled water will be added to the product

**Dilution** : PBS (pH7.4) containing 1% BSA

**Stability** : Lyophilized product, 5 years at 2 – 8 °C  
: Solution, 2 years at –20 °C

**Application** : Search for megakaryocytic leukemia

**Specificity** : Non-crossreacts with gp1b and gp2b-3a

---

*For research use only, not for use in diagnostic procedures.*

Distributed by:



Immuno-Biological Laboratories, Inc.  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432

Toll-Free: 888-523-1246  
Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Web: [www.IBL-America.com](http://www.IBL-America.com)