

Code No. 10133

**Anti-Human
KGF (8A4) Mouse IgG MoAb**

Volume : 200 µg

Lot No : 9K-918

Antigen : Recombinant Human KGF (*E. coli*)

Source : Mouse-Mouse hybridoma, supernatant

Clone : 8A4

Subclass : IgG₁

Purification : Affinity Purified with anti-mouse IgG

Form : Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use : 1 ml distilled water will be added to the product

Dilution : PBS (pH7.4) containing 1% BSA

Stability : Lyophilized product, 5 years at 2 – 8
: Solution, 2 years at –20

Application : This antibody can be used for western blotting in concentration of 1 ~ 5 µg/ml.

For Non-Clinical Research Use Only

IBL Co., LTD. 5-1 Aramachi, Takasaki-Shi, Gunma, 370-0831, JAPAN
TEL: +81 (0)27-310-8040 FAX: +81 (0)27-310-8045 URL: <http://www.ibl-japan.co.jp> E-Mail: do-ibl@ibl-japan.co.jp