

Code No. 11041

**Anti-Human
G-CSF (4-12-2) Mouse IgG MoAb**Volume : 1 mL

- Introduction** : G-CSF (Granulocyte Colony-Stimulating Factor) is cytokine which accelerates differentiation and multiplication of neutrophilic precursor cell. It is reported that G-CSF is produced in fibroblast and endothelial cell, and that the multiplication of G-CSF is accelerated in leukemia and some solid tumor which has G-CSF receptor.
- Antigen** : Recombinant Human G-CSF
- Source** : Mouse-Mouse hybridoma
(X63 - Ag 8.653 × BALB/c mouse spleen cells)
- Clone** : 4-12-2 **Subclass** : IgG₁
- Purification** : Affinity purified with Protein A
- Form** : Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃
- How to use** : 1.0 mL deionized water will be added to the product.
- Stability** : Lyophilized product, 5 years at 2 - 8 °C
: Solution, 2 years at -20 °C
- Application** : This antibody can be used for immunohistochemistry with formalin fixed paraffin embedded tissues after microwave pretreatment. The recommended dilution is 100-fold, however, the dilution rate should be optimized by each laboratory.

For research use only, not for use in diagnostic procedures.

Distributed by:

Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com