


**25(OH) Vitamin D2/D3
(Mouse monoclonal Antibody anti Vitamin D2 / D3)**



**For Research Use Only
Not for Use in Diagnostic Procedures**

PRODUCT ANALYSIS SHEET

Cat # 53.197.06

Clone:	LMBP7013CB
Batch:	
Size	mg / ml
Format:	Purified , Liquid
Concentration:	mg/ml in PBS 0.1M pH 7.4
Preservative:	NaN3 0.5 g/L
Purification:	Purified by Protein A affinity chromatography
Isotype:	Ig G1 kappa
Isoelectric point:	N/A
Immunogen:	Vitamin D derivate – BSA
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specificity :	Recognizes natural 25OH Vitamin D2 and Vitamin D3
Application:	Elisa
Storage & Stability:	4°C/
Update:	June 28, 2017

Date & Signature