

## mAb Hu 25(OH)-Vit D Total Capture, Purified Unconjugated



### PRODUCT SPECIFICATION

Cat # **5319716**

|                          |   |
|--------------------------|---|
| Clone:                   | LMBP 7012CB                                       |
| Batch*:                  | <i>[Batch dependent]</i>                          |
| Size                     | <i>[≥ 1 mg]</i>                                   |
| Format:                  | Purified Unconjugated, Liquid                     |
| Concentration:           | <i>[≥ 1 mg/mL, Batch dependent]</i>               |
| Buffer:                  | PBS 0,1M pH 7,4                                   |
| Preservative:            | NaN3 0,5 g/L                                      |
| Purification:            | Protein A affinity chromatography                 |
| Isotype:                 | IgG1 K  |
| Isoelectric point:       | 7.35  |
| Immunogen:               | Vitamin D derivate – BSA                          |
| Mouse / fusion partners: | Balb/c (Charles River)/ NSO                       |
| Specifications :         | Recognizes natural 25OH Vitamin D2 and Vitamin D3 |
| Applications:            | Immunoassays                                      |
| Storage:                 | 2-8°C   |
| Expiry date:             | <i>[Batch dependent]</i>                          |

\* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature

Distributed by:



Immuno-Biological Laboratories, Inc.  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432

Toll-Free: 888-523-1246  
Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Web: [www.IBL-America.com](http://www.IBL-America.com)