

mAb Hu 25(OH)-Vit D Total Capture, Purified Unconjugated



PRODUCT SPECIFICATION

Cat # **5319726**

Clone:	LMBP 7011CB
Batch*:	<i>[Batch dependent]</i>
Size	<i>[≥ 1 mg]</i>
Format:	Purified Unconjugated, Liquid
Concentration:	<i>[≥ 1 mg/mL, Batch dependent]</i>
Buffer:	PBS 0,1M pH 7,4
Preservative:	NaN3 0,5 g/L
Purification:	Protein A affinity chromatography
Isotype:	IgG1 K
Immunogen:	Vitamin D derivate – BSA
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specifications :	Recognizes natural 25OH Vitamin D2 and Vitamin D3
Applications:	Immunoassays
Storage:	2-8°C
Expiry date:	<i>[Batch dependent]</i>

* Only the figures before the last slash should be taken into account for the batch identification.

Date

Signature

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com