

Product Information

Human DAN/NBL1 antibody

Catalog Number: IBATGA0267



PRODUCT INFORMATION

Catalog number

ATGA0267

Clone No.

AT38G8

Product type

Monoclonal Antibody

UnitProt No.

P41271

NCBI Accession No.

NP_005371

Alternative Names

Neuroblastoma suppressor of tumorigenicity 1, D1S1733E, DAN, DAND1, NB, NO3

PRODUCT SPECIFICATION

Antibody Host

Mouse

Reacts With

Human

Concentration

1mg/ml (determined by BCA assay)

Formulation

Liquid in. Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% glycerol

Immunogen

Recombinant human NBL1 (18-181aa) purified from E. coli

Isotype

IgG1 kappa

Purification Note

By protein-A affinity chromatography

Application

ELISA, WB, ICC/IF, FACS

Usage

The antibody has been tested by ELISA, Western blot, ICC/IF and FACS analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:

Immuno-Biological Laboratories, Inc. (IBL-America)

8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA

Phone: (888) 523-1246

Fax.: (763) 780-2988

Email: info@ibl-america.com

Web: www.ibl-america.com

Product Information

Human DAN/NBL1 antibody

Catalog Number: IBATGA0267



Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

BACKGROUND

Description


NBL1 is the founding member of the evolutionarily conserved CAN (Cerberus and DAN) family of proteins, which contain a domain resembling the CTCK (C-terminal cystine knot-like) motif found in a number of signaling molecules. These proteins are secreted, and act as BMP (bone morphogenetic protein) antagonists by binding to BMPs and preventing them from interacting with their receptors. They may thus play an important role during growth and development and alternatively spliced transcript variants have been identified for this gene. Also, as a tumor suppressor gene of neuroblastoma may play an important role in preventing cells from entering the final stage (G1/S) of the transformation process. Read-through transcripts between this locus and the upstream mitochondrial inner membrane organizing system 1 have been observed.

General References

Olakowski, M., et al. (2009) Folia Histochem Cytobiol 47(2): 249-55.
Shaikh Ibrahim, Z., et al. (2011) Int J Mol Med 28(4): 605-11.


DATA

Western blot analysis (WB)



The cell lysates(30ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human NBL1 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.
Lane 1.: Jurkat cell lysate

Immunofluorescence (ICC/IF)



ICC/IF analysis of NBL1 in Jurkat cells. The cell was stained with ATGA0267 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.




Manufactured for:
Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax.: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com

Product Information


Human DAN/NBL1 antibody

Catalog Number: IBATGA0267



ICC/IF analysis of NBL1 in A549 cells. The cell was stained with ATGA0267 (1:100). The secondary antibody (green) was used Alexa Fluor 488. DAPI was stained the cell nucleus (blue).

Flow cytometry (FACS)



Flow cytometry analysis of NBL1 in Jurkat cell line, staining at 2-5ug for 1×10^6 cells (red line). The secondary antibody used goat anti-mouse IgG Alexa fluor 488 conjugate. Isotype control antibody was mouse IgG (black line).

For research use only. This product is not intended or approved for human, diagnostics or veterinary use.



Manufactured for:
Immuno-Biological Laboratories, Inc. (IBL-America)
8201 Central Ave. NE, Suite P, Minneapolis, Minnesota 55432, USA
Phone: (888) 523-1246 Fax: (763) 780-2988
Email: info@ibl-america.com Web: www.ibl-america.com