

Code No. 16201

**Anti-Human
Endothelin-A Receptor (A405) Rabbit IgG Affinity Purify**Volume : 100 µg

Introduction : Two kinds of receptors for Endothelin have been found, namely Endothelin-A (ET-A) receptor and Endothelin-B (ET-B) receptor. ET-A receptor is ET-1 specific and expressed in lung, aorta, atrium, liver, muscle, cerebral cortex, cerebellum, kidney, adrenal, duodenum, colon and placenta at the mRNA or cDNA level, while ET-B receptor lacks isopeptide specificity and is found in liver, lung, placenta, kidney, heart and brain. Cellular localization of ET-receptors in these organs or tissues, however, has not been elucidated, although stromal cells, but not epithelial cells of the breast, express ET-1 receptor *in vitro*.

Antigen : Synthetic peptide of the C terminal part of human Endothelin-A receptor (IQWKNHEQNNHNT E)

Purification : Purified with antigen peptide

Form : Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

How to use : 1.0 mL deionized water will be added to the product (the conc. comes up 100 µg /mL)

Stability : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C

Application : This antibody can be used for immunohistochemistry with formalin fixed paraffin embedded tissues by several method. The optimal concentration is about 5 µg/mL, however, the concentration should be optimized by each laboratory.
: This antibody can be used for western blotting in concentration of about 5 µg /mL.

Specificity : Cross-reacts with bovine and rat.

Reference : 1. Hiraki H., Hoshi N., Hasegawa H., Tanigawa T., Emura I., Seito T., Yamaki T., Fukuda T., Watanabe K., and Suzuki T. Regular immunohistochemical localization of endothelin-1 and endothelin-B receptor in normal, hyperplastic and neoplastic

For research use only, not for use in diagnostic procedures.

Distributed by:



Immuno-Biological Laboratories, Inc. Toll-Free: 888-523-1246

8201 Central Ave NE, Suite P

Minneapolis, MN 55432

Email: info@IBL-America.comWeb: www.IBL-America.com