

Code No. 17312

**Rabbit IgG**

Volume : 1 mg

---

**Introduction** : The product is normal rabbit IgG purified from normal rabbit serum.

**Purification** : Purified with protein A column

**Form** : Lyophilized product in PBS (without additives)

**How to use** : 1.0 mL deionized water will be added to the product

**Dilution** : PBS (pH7.4)

**Stability** : Lyophilized product, 5 years at 2 - 8 °C  
: Solution, 2 years at -20 °C

**Application** : Negative control in immunohistochemistry, western blotting or some tests.

---

*For Non-Clinical Research Use Only*