

Code No. 17507

**Anti-Human
IgG (Fc) Goat IgG Fab'- HRP**

Volume : 0.2 mg

Antigen : Human IgG (Fc)

Label : HRP

Purification : Purified with antigen

Form : Lyophilized product from 1 % BSA in PB (pH6.5)
(containing ProClin[®]300 as preservative)

How to use : 1.0 mL deionized water will be added to the product

Stability : Lyophilized product, 5 years at -20 °C
: Solution, 2 years at -20 °C

Application : This antibody can be used for western blotting in concentration of 10 - 100 ng/mL. (ECL)
: This antibody can be used for ELISA or cell ELISA in concentration of 250 - 500 ng/mL.

Specificity : The antibody is absorbed by goat IgG, mouse IgG, rabbit IgG and rat IgG.

For Non-Clinical Research Use Only