

Code No. 17508

**Anti-Mouse
IgG (H+L) Rabbit IgG Fab'- HRP**Volume : 0.2 mG
Lot No : 0H-319

Antigen : Mouse IgG (H+L)

Label : HRP

Purification : Purified with antigen peptide

Form : Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use : 1.0 ml distilled water will be added to the product

Dilution : PBS (pH7.4) containing 1% BSA

Stability : Lyophilized product, 5 years at 2 – 8
: Solution, 2 years at –20

Application : This antibody can be used for immunohistochemistry in concentration of 1 ~ 10 µg/ml.
: This antibody can be used for western blotting in concentration of 1 ~ 10 µg/ml.

Specificity : The antibody is absorbed by goat IgG, rabbit IgG and human IgG, but is not absorbed by rat IgG.

For Non-Clinical Research Use Only

IBL Co., LTD. 5-1 Aramachi, Takasaki-Shi, Gunma, 370-0831, JAPAN
TEL: +81 (0)27-310-8040 FAX: +81 (0)27-310-8045 URL: <http://www.ibl-japan.co.jp> E-Mail: do-ibl@ibl-japan.co.jp