

Code No. 17567

**Anti-Human  
IgG (Fc) Goat IgG Affinity Purify**

Volume : 1 mg  
Lot No :

---

Antigen : Human IgG (Fc)

Immunized animal :Goat

Purification : Affinity Purified with antigen

Process of absorption : This antibody was absorbed with Goat IgG, Rabbit IgG, Mouse IgG, Rat IgG.

Form : Lyophilized product from PBS containing 1% BSA and 0.05% NaN<sub>3</sub>.

How to use : 1mL deionized water will be added to the product.

Application : This antibody can be used for IHC (1 - 10µg/mL) and WB (1 - 10µg/mL).

Method of Dilution :Dilute with PBS containing 1% BSA (pH 7.4)

Stability : Lyophilized product, 5 years at 2 - 8°C  
: Solution, 2years at -20°C

---

*For research use only, not for use in diagnostic procedures.*

Distributed by:



Immuno-Biological Laboratories, Inc.  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432

Toll-Free: 888-523-1246  
Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Web: [www.IBL-America.com](http://www.IBL-America.com)