

Code No. [17568](#)

**Anti-Mouse
IgG (H+L) Rabbit IgG Affinity Purify**

Volume : 1 mg

- Antigen** : Mouse IgG (H+L)
- Purification** : Purified with antigen
- Form** : Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃
- How to use** : 1.0 mL deionized water will be added to the product
- Dilution** : PBS (pH7.4) containing 1% BSA
- Stability** : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C
- Application** : This antibody can be used for immunohistochemistry in concentration of 1 - 10 µg/mL.
: This antibody can be used for western blotting in concentration of 1 - 10 µg/mL
- Specificity** : The antibody is absorbed by goat IgG, rabbit IgG and human IgG, but is not absorbed by rat IgG.

For research use only, not for use in diagnostic procedures.

Distributed by:



Immuno-Biological Laboratories, Inc.
8201 Central Ave NE, Suite P
Minneapolis, MN 55432

Toll-Free: 888-523-1246
Email: info@IBL-America.com
Web: www.IBL-America.com