

Code No. 17662

**Anti-Mouse  
IgG (H+L) Goat IgG Affinity Purify**Volume : 1mg  
Lot No. :

- 
- Antigen** : Mouse IgG (H+L)
- Purification** : Purified with antigen
- Form** : Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>
- How to use** : 1.0 mL deionized water will be added to the product
- Stability** : Lyophilized product, 5 years at 2 - 8 °C  
: Solution, 2 years at -20 °C
- Application** : This antibody can be used for immunohistochemistry in concentration of 1 - 10 µg/mL.  
: This antibody can be used for western blotting in concentration of 1 - 10 µg/mL.
- Specificity** :The antibody is absorbed by rat IgG, goat IgG, rabbit IgG and human IgG.  
Cross-reactivity to goat IgG, rabbit IgG and human IgG is less than 0.1 %.  
Cross-reactivity to rat IgG is about 0.5 %.

---

*For research use only, not for use in diagnostic procedures.*

Distributed by:

Immuno-Biological Laboratories, Inc.  
8201 Central Ave NE, Suite P  
Minneapolis, MN 55432Toll-Free: 888-523-1246  
Email: [info@IBL-America.com](mailto:info@IBL-America.com)  
Web: [www.IBL-America.com](http://www.IBL-America.com)