

Code No. 18203

**Anti-Human
Big Endothelin-1 (22-38) Rabbit IgG Affinity Purify**Volume : 50 µg

Introduction : Endothelins are 21 amino acid peptides that are produced ubiquitously by vascular endothelial cells, smooth muscle cells, and other cells in different organs. Endothelins are secreted as big-endothelins that are converted to active proteins by the endothelin-converting enzyme. These peptides possess many biological activities, such as vasoconstriction and mitogenesis, and are involved in numerous physiological and pathophysiological processes in humans. Elevated plasma levels of endothelin have been associated with heart failure, and increased immunoreactivity for endothelin is observed in transplant coronary artery disease. Endothelin-1 (ET-1) is known as the most potent vasoconstrictor substance.

Antigen : Synthetic peptides of human big endothelin-1 (22-38)

Purification : Purified with antigen peptide

Form : Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

How to use : 1.0 mL deionized water will be added to the product (the conc. comes up 50 µg /mL)

Stability : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C

Application : This antibody can be used for immunohistochemistry with formalin fixed paraffin embedded tissues by several techniques such as SAB Method. The optimal concentration is 3 - 5 µg/mL, however, the concentration should be optimized by each laboratory.

Specificity :	Endothelin-1	0 %
	Endothelin-2	0 %
	Endothelin-3	0 %
	Big Endothelin-1	100.0 %
	Big Endothelin-2	<0.1 %
	Big Endothelin-3	<0.1 %
	VIC	<0.1 %

For research use only, not for use in diagnostic procedures.

Distributed by:



Immuno-Biological Laboratories, Inc. Toll-Free: 888-523-1246
8201 Central Ave NE, Suite P Email: info@IBL-America.com
Minneapolis, MN 55432 Web: www.IBL-America.com